# Veritor™ System

For Rapid Detection of Flu A+B

CLIA-waived kit configured for testing nasal and nasopharyngeal swab samples freshly collected, processed and dispensed directly onto assay test device.

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For most current information regarding the Veritor System product, refer to [BD.com/ds/VeritorSystem](http://BD.com/ds/VeritorSystem)

## For Rapid Detection of Flu A+B

For use with nasal and nasopharyngeal swab specimens.

### CLIA Complexity-WAIVED

For *in vitro* diagnostic use only.

### R<sub>x</sub> Only

A Certificate of Waiver is required to perform this test in a CLIA waived setting. To obtain a Certificate of Waiver, please contact your state health department. Additional CLIA waiver information is available at the Centers for Medicare and Medicaid website at [www.cms.hhs.gov/CLIA](http://www.cms.hhs.gov/CLIA) or from your state health department.

Failure to follow the instructions or modification to the test system instructions will result in the test no longer meeting the requirements for waived category.

### INTENDED USE

The **BD Veritor** System for Rapid Detection of Flu A+B is a rapid chromatographic immunoassay for the direct and qualitative detection of influenza A and B viral nucleoprotein antigens from nasal and nasopharyngeal swabs of symptomatic patients. The **BD Veritor** System for Rapid Detection of Flu A+B (also referred to as the **BD Veritor** System and **BD Veritor** System Flu A+B) is a differentiated test, such that influenza A viral antigens can be distinguished from influenza B viral antigens from a single processed sample using a single device. The test is to be used as an aid in the diagnosis of influenza A and B viral infections. A negative test is presumptive and it is recommended that these results be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay. Outside the U.S., a negative test is presumptive and it is recommended that these results be confirmed by viral culture or a molecular assay cleared for diagnostic use in the country of use. FDA has not cleared this device for use outside of the U.S. Negative test results do not preclude influenza viral infection and should not be used as the sole basis for treatment or other patient management decisions. The test is not intended to detect influenza C antigens.

Performance characteristics for influenza A and B were established during January through March of 2011 when influenza viruses A/2009 H1N1, A/H3N2, B/Victoria lineage, and B/Yamagata lineage were the predominant influenza viruses in circulation according to the *Morbidity and Mortality Weekly Report* from the CDC entitled "Update: Influenza Activity—United States, 2010–2011 Season, and Composition of the 2011–2012 Influenza Vaccine." Performance characteristics may vary against other emerging influenza viruses.

If infection with a novel influenza virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to the state or local health department for testing. Virus culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.

### SUMMARY AND EXPLANATION

Influenza illness classically presents with sudden onset of fever, chills, headache, myalgias, and a non-productive cough. Epidemics of influenza typically occur during winter months with estimated 114,000 hospitalizations<sup>1</sup> and 36,000 deaths<sup>2</sup> per year in the U.S. Influenza viruses can also cause pandemics, during which rates of illness and death from influenza-related complications can increase dramatically.

Patients who present with suspected influenza may benefit from treatment with an antiviral agent especially if given within the first 48 hours of onset of illness. It is important to rapidly distinguish influenza A from influenza B in order to allow physicians a choice in selective antiviral intervention. Moreover, it is important to determine if influenza A or B is causing symptomatic disease in a particular institution (e.g., nursing home) or community, so that appropriate preventative intervention can be taken for susceptible individuals. It is therefore important to not only rapidly determine whether influenza is present, but also which type of influenza virus is present.<sup>3</sup>

Diagnostic tests available for influenza include rapid immunoassay, immunofluorescence assay, polymerase chain reaction (PCR), serology, and viral culture.<sup>4–11</sup> Immunofluorescence assays entail staining of specimens immobilized on microscope slides using fluorescent-labeled antibodies for observation by fluorescence microscopy.<sup>6,12,13</sup> Culture methods employ initial viral isolation in cell culture, followed by hemadsorption inhibition, immunofluorescence, or neutralization assays to confirm the presence of the influenza virus.<sup>13–15</sup>

The **BD Veritor** System for Rapid Detection of Flu A+B (also referred to as the **BD Veritor** System and **BD Veritor** System Flu A+B) is a chromatographic immunoassay to detect influenza A or B nucleoprotein antigens from respiratory specimens of symptomatic patients with a time to result of 10 minutes. The speed and simplified workflow of the **BD Veritor** System for Rapid Detection of Flu A+B makes it applicable as a "STAT" influenza A and B antigen detection test providing relevant information to assist with the diagnosis of influenza.

All **BD Veritor** System Flu A+B test devices are interpreted by a **BD Veritor** System Instrument, either a **BD Veritor** Reader or **BD Veritor** Plus Analyzer (the "Analyzer"). When using the **BD Veritor** Plus Analyzer, workflow steps depend on the selected operational mode and the Analyzer configuration settings. In **Analyze Now mode**, the instrument evaluates assay

devices after manual timing of their development. In **Walk Away mode**, devices are inserted immediately after application of the specimen, and timing of assay development and analysis is automated. Additionally, connection of the Analyzer to a printer is possible if desired. Additional result documentation capabilities are possible with the integration of the **BD Veritor InfoScan** ("InfoScan") or **BD Veritor InfoSync** ("InfoSync") module. Please refer to the Analyzer *Instructions for Use* for details on how to implement these features. The **BD Veritor InfoSync** is not available in all regions.

**PRINCIPLES OF THE PROCEDURE**

The **BD Veritor** System for Rapid Detection of Flu A+B is a qualitative, digital immunoassay for the detection of influenza A and B viral antigens in samples processed from respiratory specimens. When specimens are processed and added to the test device, influenza A or B viral antigens bind to anti-influenza antibodies conjugated to detector particles in the A+B test strip. The antigen-conjugate complex migrates across the test strip to the reaction area and is captured by the line of antibody on the membrane. A positive result for influenza A is determined by the **BD Veritor** System Instrument when antigen-conjugate is deposited at the Test "A" position and the Control "C" position on the **BD Veritor** System Flu A+B assay device. A positive result for influenza B is determined by the **BD Veritor** System Instrument when antigen-conjugate is deposited at the Test "B" position and the Control "C" position in the **BD Veritor** System Flu A+B assay device. The instrument analyzes and corrects for non-specific binding and detects positives not recognized by the unaided eye to provide an objective digital result.

**REAGENTS**

The following components are included in the **BD Veritor** System for Rapid Detection of Flu A+B kit:

|                                          |                              |                                                                                                                                                                                                          |
|------------------------------------------|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BD Veritor</b> System Flu A+B Devices | 30 devices                   | Foil pouched device containing one reactive strip. Each strip has two test lines of monoclonal antibody specific to either influenza A or B viral antigen and murine monoclonal control line antibodies. |
| <b>RV Reagent D</b>                      | 30 tubes with 400 µL reagent | Detergent with < 0.1% sodium azide                                                                                                                                                                       |
| Flexible minitip flocked swab            | 30 each                      | Swab for nasopharyngeal or nasal collection                                                                                                                                                              |
| Control A+/B- Swab                       | 1 each                       | Flu A Positive and Flu B Negative Control Swab, influenza A antigen (inactive recombinant nucleoprotein) with < 0.1% sodium azide                                                                        |
| Control B+/A- Swab                       | 1 each                       | Flu A Negative and Flu B Positive Control Swab, influenza B antigen (inactive recombinant nucleoprotein) with < 0.1% sodium azide                                                                        |

**Materials Required But Not Provided:** **BD Veritor™** System Reader (Cat. No. 256055) or **BD Veritor™** Plus Analyzer (Cat. No. 256066), Timer, Tube Rack for specimen testing.

**Optional Equipment:** **BD Veritor™** InfoScan Module (Cat. No. 256068), **BD Veritor™** InfoSync Module (Cat. No. 256067), USB Printer cable for **BD Veritor™** Analyzer (Cat. No. 443907), Epson Printer model TM-T20 II.

**Warnings and Precautions:**

**Warning**

**H302** Harmful if swallowed. **H402** Harmful to aquatic life. **H412** Harmful to aquatic life with long lasting effects.  
**P273** Avoid release to the environment. **P264** Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

- For *in vitro* Diagnostic Use.
- Test results are not meant to be visually determined. **All test results must be determined using the BD Veritor System Instrument.**
- If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to the state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.
- Pathogenic microorganisms, including hepatitis viruses, Human Immunodeficiency Virus and novel influenza viruses, may be present in clinical specimens. "Standard Precautions"<sup>16-19</sup> and institutional guidelines should be followed in handling, storing and disposing of all specimens and all items contaminated with blood and other body fluids.
- Dispose of used **BD Veritor** System test devices as biohazardous waste in accordance with federal, state and local requirements.
- Reagents contain sodium azide, which is harmful if inhaled, swallowed or exposed to skin. Contact with acids produces very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
- For optimal results, use the flocked swabs provided with the kit for specimen collection.

8. Other than the flocked swabs that are used for specimen collection, kit components should not make contact with the patient.
9. Do not use kit components beyond the expiration date.
10. Do not reuse the device.
11. Do not use the kit if the Control A+/B- swab and Control B+/A- swab do not yield appropriate results.
12. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
13. To avoid erroneous results, swab specimens must be processed as indicated in the assay procedure section.
14. Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.
15. FluMist® is made from attenuated live flu virus and although the concentration tested (1%) was non-interfering, it is possible when tested with higher concentrations that an influenza A and/or influenza B false positive may occur.

**Storage and Handling:** Kits may be stored at 2–30 °C. DO NOT FREEZE. **Reagents and devices must be at room temperature (15–30 °C) when used for testing.**


## **SPECIMEN COLLECTION**

Acceptable specimens for testing with the **BD Veritor** System Flu A+B test include nasal swabs and nasopharyngeal (NP) swabs. Freshly collected specimens should be processed within 1 hour. It is essential that correct specimen collection and preparation methods be followed. Specimens obtained early in the course of the illness will contain the highest viral titers.

Inadequate specimen collection, improper specimen handling and/or transport may yield a false negative result; therefore, specimen collection requires specific training and guidance due to the importance of specimen quality to accurate test results.

### **Proper Nasal Swab Sample Collection**

1. The **BD Veritor** System Kit includes swabs with a flocked tip for nasal specimen collection.
2. Insert the swab into one nostril of the patient.
3. Rotate the swab two complete 360-degree turns; pressing firmly against the nasal mucosa to help ensure the swab obtains both cells and mucus.


4. Withdraw the swab from the nasal cavity. The sample is now ready for processing using the **BD Veritor** System Kit.

#### Proper Nasopharyngeal Swab Sample Collection

1. The **BD Veritor** System Kit includes swabs with a flocked tip for nasopharyngeal specimen collection.
2. Insert the swab into one nostril of the patient, reaching the surface of the posterior nasopharynx.
3. Rotate the swab over the surface of the posterior nasopharynx.
4. Withdraw the swab from the nasal cavity. The sample is now ready for processing using the **BD Veritor** System Kit.

#### DOs and DON'Ts of Sample Collection

- Do collect sample as soon as possible after onset of symptoms
- Do test sample immediately
- BD recommends flocked swabs which are provided in the **BD Veritor** System Flu A+B Kit
- Do not use cotton tips and wood shafts
- Do not use calcium alginate swabs

TEST PROCEDURE

Nasal and Nasopharyngeal Swab Test Procedure

NOTES:


- Reagents, specimens and devices must be at room temperature (15–30 °C) prior to testing.
- The CLIA-waived **BD Veritor** System for Rapid Detection of Flu A+B kit is only intended for nasal and nasopharyngeal swab specimens that are collected and tested directly (i.e., dry swabs that have **NOT** been placed in transport media). The kit includes a pre-diluted process reagent in a ready to use “unitized” tube. This CLIA-waived kit IS NOT INTENDED for testing liquid samples such as wash or aspirate samples or swabs in transport media as results can be compromised by over dilution.

Prepare for Testing

The following steps assume that users of a **BD Veritor** Plus Analyzer have chosen and set all configuration options, and that the Analyzer is ready to use. To choose or change these settings, see the **BD Veritor** Plus Analyzer *Instructions for Use*, section 4.7. A printer is not necessary to display results. However, if your facility has chosen to connect the Analyzer to a printer, check that the printer is plugged into a power source, paper supply is adequate and any necessary network connections are enabled before testing.

Step 1: For each patient specimen:

- Remove one **RV Reagent D** tube/tip and one **BD Veritor** System Flu A+B device from its foil pouch immediately before testing.
- Label with patient's name or ID number.
- Place the labeled **RV Reagent D** tube(s) in the designated area of the tube rack.

Prepare the Sample

Step 2:

- Remove and discard the cap from the **RV Reagent D** tube corresponding to the sample to be tested.


Step 3:

- Insert the patient sample swab all the way into the **RV Reagent D** tube and swirl it against the inside wall three (3) times.


Step 4:

- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Properly discard of the swab.

Step 5:

- Press the attached tip firmly onto the **RV Reagent D** tube containing the processed sample (threading/twisting not required).


- Vortex or mix thoroughly by swirling or flicking the bottom of the tube before adding to assay device.
- Note: Do not use tips from any other product, including other products from BD or other manufacturers.

| After step 5, choose from the model and workflow option below before continuing to step 6: |                                                                                                                                                                                  |                                                                                              |                                                                                                                                          |                                                                                              |
|--------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|
|                                                                                            | BD Veritor Reader or Analyzer in <b>Analyze Now</b> mode                                                                                                                         | BD Veritor Plus Analyzer in <b>Walk Away</b> mode                                            | BD Veritor Plus Analyzer with either the BD Veritor InfoScan or InfoSync module In <b>Analyze now</b> mode---or--- <b>Walk Away</b> mode |                                                                                              |
| Instructions in section:                                                                   |  <b>A</b>  | <b>B</b>  | <b>C</b>                                              | <b>D</b>  |

**Step 6A: Adding the specimen**

- Invert the **RV Reagent D** tube and hold the tube vertically (approximately one inch above the labeled **BD Veritor** System Flu A+B device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled **BD Veritor** System Flu A+B device.


**NOTE:** Squeezing the tube too close to the tip may cause leakage.

**Step 7A: Timing development**

- After adding the sample, allow the test to run for 10 minutes before inserting into the **BD Veritor** Instrument.
- **CAUTION: incorrect results may occur if development time is less than 10 minutes.** Some lines may appear on the device sooner. Do not read device visually.
- **NOTE:** If running test under laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.

**Step 8A: Using the BD Veritor Instrument**

- During incubation time, turn the **BD Veritor** Instrument on by pressing the power button once.
- Insert assay device when the 10-minute assay development time is complete. Follow the on-screen prompt to complete the procedure.
- The status of the assay analysis process appears in the display window, followed by the result display.

**Step 9A: Record the Result**


- When analysis is complete, the test result appears in the display window. Record the result and discard the test device appropriately.

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).**


To use Walk Away mode - connect the AC power adapter to the Analyzer and a power source

**Step 6B: Starting Walk Away Mode**

- Turn the Analyzer on by pressing the blue power button once
- When the display window reads: "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE"
  - **Double-click** the blue power button.

**Step 7B: Adding the specimen**

- When the display window reads "ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY"
  - Invert the tube, holding it vertically (approximately one inch above the **BD Veritor System Flu A+B** device sample well).
  - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled **BD Veritor System Flu A+B** device.

**NOTE:** Squeezing the tube too close to the tip may cause leakage.

**CAUTION:** A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.


**Step 8B: Starting the development and reading sequence**

- Insert the test device into the slot on the right side of the Analyzer.

**The test device must remain horizontal to prevent spilling the specimen out of the sample well**

- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- The display window shows the remaining analysis time.

**Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**

**Step 9B: Record the Result**


- When analysis is complete, the test result appears in the display window. Record the result and discard the test device appropriately.

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).**

**Step 6C: Adding the specimen**

- Invert the tube, holding it vertically (approximately one inch above the **BD Veritor** System Flu A+B device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled **BD Veritor** System Flu A+B device.

**NOTE:** Squeezing the tube close to the tip may cause leakage.

**Step 7C: Timing development**


- Allow the test to develop for 10 minutes.
- **CAUTION: incorrect results may occur if development time is less than 10 minutes.** Some lines may appear on the device sooner. Do not read device visually.
- If running the test in a laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.

**Step 8C: Using the Analyzer**


During the incubation time, turn the Analyzer on by pressing the blue button once.

The display window briefly shows “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration steps. Ignore this message and postpone this process when an assay is awaiting analysis.

- When assay development time is complete and the Analyzer display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”:
  - Insert the **BD Veritor** System Flu A+B device into the slot on the right side of the Analyzer.

**Step 9C: Using the Bar Code scanner**

- Follow the prompts on the display window to complete any required barcode scans of:
  - OPERATOR ID
  - SPECIMEN ID and/or
  - KIT LOT NUMBER

- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 8C. To restart this step, remove and reinsert the test device to initiate a new reading sequence.
- Move barcodes slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the kit Lot Number and expiration date in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

After required scans are completed, the Analyzer displays a countdown timer and test analysis begins.

- **Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**
- When analysis is complete a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.

**If a printer is not connected, record the result before removing the assay device.**

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).**

**Step 10C: Remove the test device**

- Remove and then discard the test device appropriately. The display will show “INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE” to indicate the Analyzer is ready to perform another test.

If the **BD Veritor** InfoSync module is installed the ENVELOPE symbol will appear to indicate that results are transmitting.

- In the event that the Analyzer does not detect adequate cellular network strength while the ENVELOPE symbol is still displayed, it will queue all results to be transmitted and continuously attempt to transmit them. If it is powered off during this time, it will attempt to transmit as soon as power is restored.

To use Walk Away mode - connect the AC power adapter to the Analyzer and a power source

**Step 6D: Starting Walk Away mode**


- Turn the Analyzer on by pressing the blue power button once.

The display window will briefly show "SCAN CONFIG BARCODE." This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration steps. Ignore this message and postpone this process when an assay is awaiting analysis.

- When the display window reads:  
"INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE"  
– **Double-click** the blue power button.

**Step 7D: Using the Bar Code scanner**


- Follow the prompts on the display window to complete any required barcode scans of:
  - OPERATOR ID
  - SPECIMEN ID and/or
  - KIT LOT NUMBER


- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 6D. To restart this step, double-click the power button.
- Move barcodes slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the kit Lot Number and expiration date in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

**Step 8D: Add the specimen to the test device**

- When the display window reads: ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY:
  - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
  - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device. **NOTE: Squeezing the tube close to the tip may cause leakage.**
  - **CAUTION: A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.**

**Step 9D: Starting the development and reading sequence**

- Insert the test device into the slot on the right side of the Analyzer. **The test device must remain horizontal to prevent spilling the specimen out of the sample well.**
- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- The display window shows the remaining analysis time.

**Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**

- When analysis is complete, a result appears in the display window. If configured to display, the Specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed. **If a printer is not connected, record the result before removing the assay device.**

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).**

**Step 10D: Removing the test device**

- Remove and then discard the test device appropriately. The display will show INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE to indicate the Analyzer is ready to perform another test. Note that the Analyzer returns to Analyze Now mode at the conclusion of each read sequence.

If the BD Veritor InfoSync module is installed the ENVELOPE symbol will appear to indicate that results are transmitting.

- In the event that the Analyzer does not detect adequate cellular network strength while the ENVELOPE symbol is still displayed, it will queue all results to be transmitted and continuously attempt to transmit them. If it is powered off during this time, it will attempt to transmit as soon as power is restored.

#### OPTIONAL TEST PROCEDURE: Testing for INFLUENZA A+B and RSV using a single NP swab

**Note:** The BD Veritor™ System for Rapid Detection of RSV (Cat. # 256038) is required for this procedure in addition to the BD Veritor™ System for Rapid Detection of Flu A+B (Cat. # 256045).

**IMPORTANT NOTE: THE SAMPLE TO BE TESTED IN THE RSV KIT MUST BE FROM A PATIENT LESS THAN 6 YEARS OF AGE AS INDICATED IN THE BD VERITOR RSV POC KIT PACKAGE INSERT. THE PROCESSED SAMPLE SHOULD BE TESTED WITHIN 15 MINUTES.**

This procedure allows for use of the remaining processed sample from Step 5 above to test additionally for RSV. When using this optional test procedure, the sample may be used up to 15 minutes after initial processing.

1. Collect NP swab from the patient and follow Steps 1-5 of the test procedure above as instructed for the **BD Veritor** System for Rapid Detection of Flu A+B.
2. Using the sample from Step 5, continue the test procedure using the test device for RSV.
3. Refer to the product insert for **BD Veritor**™ System for Rapid Detection of RSV, (Cat. # 256038) for the test procedure and full description of the **BD Veritor** RSV test. Follow the Instructions in the insert and the Instrument on-screen prompts to complete the test procedure and obtain results. Refer to the product insert for the **BD Veritor**™ System RSV Kit (Cat. # 256038) for result interpretation.

#### INTERPRETATION OF RESULTS

The **BD Veritor** System Instrument (purchased separately) must be used for interpretation of all test results. Operators should not attempt to interpret assay results directly from the test strip contained within the **BD Veritor** System Flu A+B assay device. With some specimens, up to four lines may be visible on the test device. The Instrument will appropriately interpret the result.

| Display                  | Interpretation                                          |
|--------------------------|---------------------------------------------------------|
| FLU A: +<br>FLU B: -     | Positive Test for Flu A (influenza A antigen present)   |
| FLU A: -<br>FLU B: +     | Positive Test for Flu B (influenza B antigen present)   |
| FLU A: -<br>FLU B: -     | Negative Test for Flu A and Flu B (no antigen detected) |
| RESULT INVALID           | Result Invalid. Repeat the test.                        |
| POSITIVE CONTROL INVALID | Test Invalid. Repeat the test.                          |
| NEGATIVE CONTROL INVALID | Test Invalid. Repeat the test.                          |

**Invalid Test** – If the test is invalid, the **BD Veritor** System Instrument will display “RESULT INVALID” or “POSITIVE CONTROL INVALID” or “NEGATIVE CONTROL INVALID” and the test or control must then be repeated. Because true dual positives are exceptionally rare, the **BD Veritor** System Instrument reports dual positive influenza A and influenza B results as “Result Invalid.” Specimens generating a “Result Invalid” should be retested. Upon retesting, if the specimen produces a “Result Invalid” the user may want to consider other methods to determine whether the sample is positive or negative for influenza virus. If either the POSITIVE or NEGATIVE “CONTROL INVALID” reading recurs, contact BD Technical Support.

#### REPORTING OF RESULTS

|                      |                                                                                                                                                                                                                                                                                                                                 |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Positive Test</b> | Positive for the presence of influenza A or influenza B antigen. A positive result may occur in the absence of viable virus.                                                                                                                                                                                                    |
| <b>Negative Test</b> | Negative for the presence of influenza A or influenza B antigen. Infection due to influenza cannot be ruled out because the antigen present in the sample may be below the detection limit of the test. It is recommended that these results be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay. |
| <b>Invalid Test</b>  | Test result is inconclusive. Do not report results. Repeat the test.                                                                                                                                                                                                                                                            |

## QUALITY CONTROL

**ATTENTION:** To utilize the BD Veritor Plus Analyzer's QC documentation capability, the Analyzer must be equipped with either the BD Veritor InfoScan or InfoSync module and specimen barcode scanning must be enabled. Please refer to the Analyzer *Instructions for Use*, section 4, to choose or change this configuration, and for specific QC procedure steps.

Each **BD Veritor** System Flu A+B device contains both positive and negative internal/procedural controls:

1. The internal positive control validates the immunological integrity of the device, proper reagent function, and assures correct test procedure.
2. The membrane area surrounding test lines functions as a background check on the assay device.

**The BD Veritor System Instrument evaluates the positive and negative internal/procedural controls after insertion of each BD Veritor System test device. The BD Veritor Instrument also prompts the operator if a quality issue occurs during assay analysis. Failure of the internal/procedural controls will generate an invalid test result. NOTE: The internal controls do not assess proper sample collection technique.**

### External Positive and Negative Controls:

Flu A positive/B negative and Flu B positive/A negative control swabs are supplied with each kit. These controls provide additional quality control material to assess that the test reagents and the **BD Veritor** System Instrument perform as expected. Prepare kit control swabs and test using the same procedure (either **Analyze Now** or **Walk Away** mode) as used for patient specimen swabs. When using the barcode scanning feature to document QC procedures, scan the barcode on the control swab packaging when prompted for a Specimen ID.

**Your laboratory's standard Quality Control procedures and applicable local, state and/or federal regulations or accreditation requirements dictate the performance of external quality control procedures.**

BD recommends controls be run once for:

- each new kit lot,
- each new operator,
- each new shipment of test kits,
- as required by site quality control procedures and in accordance with local, state and federal regulations or accreditation requirements.

### Control Swab Test Procedure – Preparing the sample

1. Insert the swab all the way into the appropriately labeled **RV Reagent D** tube and vigorously plunge the swab up and down in the fluid for a minimum of 15 seconds.
2. Continue processing the swab according to the Test Procedure for Nasal and Nasopharyngeal swabs above beginning at Step 4 **"Remove the swab."**

**If the kit controls do not perform as expected, do not report patient results. Contact BD Technical Services at 1.800.638.8663.**

## LIMITATIONS OF THE PROCEDURE

- Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- The contents of this kit are to be used for the qualitative detection of influenza type A and B antigens from nasal swab and nasopharyngeal swab specimens.
- The **BD Veritor** System for Rapid Detection of Flu A+B is capable of detecting both viable and non-viable influenza particles. The **BD Veritor** System for Rapid Detection of Flu A+B performance depends on antigen load and may not correlate with other diagnostic methods performed on the same specimen.
- Results from the **BD Veritor** System for Rapid Detection of Flu A+B test should be correlated with the clinical history, epidemiological data and other data available to the clinician evaluating the patient.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly; therefore, a negative test result does not eliminate the possibility of influenza A or influenza B infection, and should be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay.
- Positive test results do not rule out co-infections with other pathogens.
- Positive test results do not identify specific influenza A virus subtypes.
- Negative test results are not intended to rule in other non-influenza viral or bacterial infections.
- Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false positive results during periods of little/no influenza activity when disease prevalence is low. False negative test results are more likely during peak influenza activity when prevalence of disease is high.
- This device has been evaluated for use with human specimen material only.

- Monoclonal antibodies may fail to detect, or detect with less sensitivity, influenza A viruses that have undergone minor amino acid changes in the target epitope region.
- The analytical reactivity of this device has not been established for avian or swine origin influenza strains other than those included in the "strain reactivity" tables.
- The performance of this test has not been evaluated for use in patients without signs and symptoms of respiratory infection.
- The **BD Veritor** System Instrument reports dual positive influenza A and influenza B results as "Result Invalid." True dual positives are exceptionally rare. Specimens generating a "Result Invalid" should be retested. Upon retesting, if the specimen produces a "Result Invalid" the user may want to consider other methods to determine whether the sample is positive or negative for influenza virus.

## EXPECTED VALUES

The rate of positivity observed in respiratory testing will vary depending on the method of specimen collection, handling/transport system employed, detection method utilized, the time of year, age of the patient, geographic location and most importantly, local disease prevalence. The overall prevalence observed with an FDA-cleared Influenza A and B molecular assay in the U.S. during the 2010-2011 clinical study was 29.9% for Influenza A and 19.7% for influenza B. The overall prevalence observed with the same FDA-cleared Influenza A and B molecular assay in Japan during the 2010-2011 clinical study was 32.2% for Influenza A and 31.7% for influenza B.

## PERFORMANCE CHARACTERISTICS

### Clinical Performance:

Performance characteristics for the **BD Veritor** System for Rapid Detection of Flu A+B test were established in multi-center Point-of-Care (POC) studies conducted at five U.S. trial sites and eight Japan trial sites during the 2010-2011 respiratory season. A total of 736 prospective specimens (515 in the U.S. and 221 in Japan) were tested using the **BD Veritor** System for Rapid Detection of Flu A+B test. These specimens consisted of nasal and nasopharyngeal swabs from symptomatic patients. In the U.S., 54% of the samples were from females and 46% from males. 20.3% of the samples were from patients less than or equal to 5 years of age, 40.8% were from patients in the 6-21 year age group, 35.6% were from 22-59 years of age, and the remaining 3.3% were obtained from people greater than or equal to age 60. In Japan, 43.3% of the samples were from females and 56.7% from males. 27.3% of the samples were from patients less than or equal to 5 years of age, 58.4% were from patients in the 16-21 year age group, 13.1% from 22-59 years of age, and 1.3% were obtained from people greater than or equal to age 60.

The performance of the **BD Veritor** System for Rapid Detection of Flu A+B test at the U.S. sites were compared to an FDA-cleared Influenza A and B molecular assay (PCR).

### Explanation of Terms:

PPA: Positive Percent Agreement =  $a / (a+c) \times 100\%$

NPA: Negative Percent Agreement =  $d / (b+d) \times 100\%$

P: Positive

N: Negative

C.I: Confidence Interval

| New Test Method | Comparator Method |       |
|-----------------|-------------------|-------|
|                 | P                 | N     |
| P               | a                 | b     |
| N               | c                 | d     |
| Total           | (a+c)             | (b+d) |

The performance is presented in Table 1 through Table 4 below.

**Table 1: Summary of the Performance Data for the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for All Swabs - All Sites**

| POC: <b>BD</b> Flu A      | Reference PCR |     |       |
|---------------------------|---------------|-----|-------|
|                           | P             | N   | Total |
| P                         | 189           | 13  | 202   |
| N                         | 37            | 497 | 534   |
| Total                     | 226           | 510 | 736   |
| Reference Method: PCR     |               |     |       |
| PPA: 83.6% (76.1%, 89.1%) |               |     |       |
| NPA: 97.5% (95.7%, 98.5%) |               |     |       |

| POC: <b>BD</b> Flu B      | Reference PCR |     |       |
|---------------------------|---------------|-----|-------|
|                           | P             | N   | Total |
| P                         | 139           | 10  | 149   |
| N                         | 32            | 555 | 587   |
| Total                     | 171           | 565 | 736   |
| Reference Method: PCR     |               |     |       |
| PPA: 81.3% (71.1%, 88.5%) |               |     |       |
| NPA: 98.2% (95.7%, 99.3%) |               |     |       |

Wald 95% Confidence intervals corrected for over-dispersion, where needed, due to potential variability between sites.

**Table 2: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for All Swabs - U.S. Sites**

|                                                                                                 | Reference PCR |     |       |                                                                                               | Reference PCR |     |       |
|-------------------------------------------------------------------------------------------------|---------------|-----|-------|-----------------------------------------------------------------------------------------------|---------------|-----|-------|
| POC: <b>BD</b> Flu A                                                                            | P             | N   | Total | POC: <b>BD</b> Flu B                                                                          | P             | N   | Total |
| P                                                                                               | 122           | 8   | 130   | P                                                                                             | 75            | 2   | 77    |
| N                                                                                               | 33*           | 352 | 385   | N                                                                                             | 26**          | 412 | 438   |
| Total                                                                                           | 155           | 360 | 515   | Total                                                                                         | 101           | 414 | 515   |
| Reference Method: PCR<br>PPA: 78.7% (95% C.I.: 71.6–84.4%)<br>NPA: 97.8% (95% C.I.: 95.7–98.9%) |               |     |       | Reference Method: PCR<br>PPA: 74.3% (95% C.I.: 65–81.8%)<br>NPA: 99.5% (95% C.I.: 98.3–99.9%) |               |     |       |

\*Of the 33 PCR positive, **BD Veritor** negative Influenza A specimens, eight were positive in the **BD Veritor** assay using a second swab specimen (reference method specimen) collected from the same patient.

\*\*Of the 26 PCR positive, **BD Veritor** negative Influenza B specimens, six were positive in the **BD Veritor** assay using a second swab specimen (reference method specimen) collected from the same patient.

**Table 3: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for Nasopharyngeal Swabs - U.S. Sites**

|                                                                                                 | Reference PCR |     |       |                                                                                                 | Reference PCR |     |       |
|-------------------------------------------------------------------------------------------------|---------------|-----|-------|-------------------------------------------------------------------------------------------------|---------------|-----|-------|
| POC: <b>BD</b> Flu A                                                                            | P             | N   | Total | POC: <b>BD</b> Flu B                                                                            | P             | N   | Total |
| P                                                                                               | 53            | 5   | 58    | P                                                                                               | 22            | 1   | 23    |
| N                                                                                               | 18            | 135 | 153   | N                                                                                               | 8             | 180 | 188   |
| Total                                                                                           | 71            | 140 | 211   | Total                                                                                           | 30            | 181 | 211   |
| Reference Method: PCR<br>PPA: 74.6% (95% C.I.: 63.4–83.3%)<br>NPA: 96.4% (95% C.I.: 91.9–98.5%) |               |     |       | Reference Method: PCR<br>PPA: 73.3% (95% C.I.: 55.6–85.8%)<br>NPA: 99.4% (95% C.I.: 96.9–99.9%) |               |     |       |

**Table 4: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for Nasal Swabs – U.S. Sites**

|                                                                                                 | Reference PCR |     |       |                                                                                                 | Reference PCR |     |       |
|-------------------------------------------------------------------------------------------------|---------------|-----|-------|-------------------------------------------------------------------------------------------------|---------------|-----|-------|
| POC: <b>BD</b> Flu A                                                                            | P             | N   | Total | POC: <b>BD</b> Flu B                                                                            | P             | N   | Total |
| P                                                                                               | 69            | 3   | 72    | P                                                                                               | 53            | 1   | 54    |
| N                                                                                               | 15            | 217 | 232   | N                                                                                               | 18            | 232 | 250   |
| Total                                                                                           | 84            | 220 | 304   | Total                                                                                           | 71            | 233 | 304   |
| Reference Method: PCR<br>PPA: 82.1% (95% C.I.: 72.6–88.9%)<br>NPA: 98.6% (95% C.I.: 96.1–99.5%) |               |     |       | Reference Method: PCR<br>PPA: 74.6% (95% C.I.: 63.4–83.3%)<br>NPA: 99.6% (95% C.I.: 97.6–99.9%) |               |     |       |

The performance of the **BD Veritor** System for Rapid Detection of Flu A+B test at the Japan sites were also compared to the same FDA-cleared Influenza A and B molecular assay (PCR) and are presented in Table 5 through Table 7.

**Table 5: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for all Swabs - Japan Sites**

|                                                                                                 | Reference PCR |     |       |                                                                                               | Reference PCR |     |       |
|-------------------------------------------------------------------------------------------------|---------------|-----|-------|-----------------------------------------------------------------------------------------------|---------------|-----|-------|
| POC: <b>BD</b> Flu A                                                                            | P             | N   | Total | POC: <b>BD</b> Flu B                                                                          | P             | N   | Total |
| P                                                                                               | 67            | 5   | 72    | P                                                                                             | 64            | 8   | 72    |
| N                                                                                               | 4             | 145 | 149   | N                                                                                             | 6             | 143 | 149   |
| Total                                                                                           | 71            | 150 | 221   | Total                                                                                         | 70            | 151 | 221   |
| Reference Method: PCR<br>PPA: 94.4% (95% C.I.: 86.4–97.8%)<br>NPA: 96.7% (95% C.I.: 92.4–98.6%) |               |     |       | Reference Method: PCR<br>PPA: 91.4% (95% C.I.: 82.5–96%)<br>NPA: 94.7% (95% C.I.: 89.9–97.3%) |               |     |       |

**Table 6: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for all Nasopharyngeal Swabs - Japan Sites**

|                                                                                                 | Reference PCR |    |       |                                                                                               | Reference PCR |    |       |
|-------------------------------------------------------------------------------------------------|---------------|----|-------|-----------------------------------------------------------------------------------------------|---------------|----|-------|
| POC: <b>BD</b> Flu A                                                                            | P             | N  | Total | POC: <b>BD</b> Flu B                                                                          | P             | N  | Total |
| P                                                                                               | 30            | 1  | 31    | P                                                                                             | 38            | 2  | 40    |
| N                                                                                               | 2             | 83 | 85    | N                                                                                             | 1             | 75 | 76    |
| Total                                                                                           | 32            | 84 | 116   | Total                                                                                         | 39            | 77 | 116   |
| Reference Method: PCR<br>PPA: 93.8% (95% C.I.: 79.9–98.3%)<br>NPA: 98.8% (95% C.I.: 93.6–99.8%) |               |    |       | Reference Method: PCR<br>PPA: 97.4% (95% C.I.: 86.8–99.5%)<br>NPA: 97.4% (95% C.I.: 91–99.3%) |               |    |       |

**Table 7: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for Nasal Swabs – Japan Sites**

|                                                                                                 | Reference PCR |    |       |                                                                                                 | Reference PCR |    |       |
|-------------------------------------------------------------------------------------------------|---------------|----|-------|-------------------------------------------------------------------------------------------------|---------------|----|-------|
| POC: <b>BD</b> Flu A                                                                            | P             | N  | Total | POC: <b>BD</b> Flu B                                                                            | P             | N  | Total |
| P                                                                                               | 37            | 4  | 41    | P                                                                                               | 26            | 6  | 32    |
| N                                                                                               | 2             | 62 | 64    | N                                                                                               | 5             | 68 | 73    |
| Total                                                                                           | 39            | 66 | 105   | Total                                                                                           | 31            | 74 | 105   |
| Reference Method: PCR<br>PPA: 94.9% (95% C.I.: 83.1–98.6%)<br>NPA: 93.9% (95% C.I.: 85.4–97.6%) |               |    |       | Reference Method: PCR<br>PPA: 83.9% (95% C.I.: 67.4–92.9%)<br>NPA: 91.9% (95% C.I.: 83.4–96.2%) |               |    |       |

### Reproducibility

The reproducibility of the **BD Veritor** System for Rapid Detection of Flu A+B test was evaluated at three POC sites. The reproducibility panel was composed of 30 simulated influenza A or B samples. These included moderate positive samples, low positive samples (near the assay limit of detection), high negative samples (i.e., containing very low concentrations of virus such that positive results occur ~5% of the time) and negative samples. The panel was tested by two operators at each site for five consecutive days. The results are summarized below.

| Reproducibility Results – Percent of Flu A Positives |                                         |                                         |                                         |                                         |
|------------------------------------------------------|-----------------------------------------|-----------------------------------------|-----------------------------------------|-----------------------------------------|
| Sample                                               | Site 1                                  | Site 2                                  | Site 3                                  | Total                                   |
| High negative H1N1 A                                 | 0% (0/30)<br>(95% C.I.: 0–11.3%)        | 10% (3/30)<br>(95% C.I.: 3.5–25.6%)     | 26.7% (8/30)<br>(95% C.I.: 14.2–44.4%)  | 12.2% (11/90)<br>(95% C.I.: 7–20.6%)    |
| Low positive H1N1 A                                  | 86.7% (26/30)<br>(95% C.I.: 70.3–94.7%) | 96.7% (29/30)<br>(95% C.I.: 83.3–99.4%) | 100% (30/30)<br>(95% C.I.: 88.6–100%)   | 94.4% (85/90)<br>(95% C.I.: 87.6–97.6%) |
| Moderate positive H1N1 A                             | 100% (30/30)<br>(95% C.I.: 88.6–100%)   | 100% (30/30)<br>(95% C.I.: 88.6–100%)   | 100% (30/30)<br>(95% C.I.: 88.6–100%)   | 100% (90/90)<br>(95% C.I.: 95.9–100%)   |
| High negative H3N2 A                                 | 0% (0/30)<br>(95% C.I.: 0–11.3%)        | 10% (3/30)<br>(95% C.I.: 3.5–25.6%)     | 16.7% (5/30)<br>(95% C.I.: 7.3–33.6%)   | 8.9% (8/90)<br>(95% C.I.: 4.6–16.6%)    |
| Low positive H3N2 A                                  | 100% (30/30)<br>(95% C.I.: 88.6–100%)   | 93.3% (28/30)<br>(95% C.I.: 78.7–98.2%) | 96.7% (29/30)<br>(95% C.I.: 83.3–99.4%) | 96.7% (87/90)<br>(95% C.I.: 90.7–98.9%) |
| Moderate positive H3N2 A                             | 100% (30/30)<br>(95% C.I.: 88.6–100%)   | 100% (30/30)<br>(95% C.I.: 88.6–100%)   | 100% (30/30)<br>(95% C.I.: 88.6–100%)   | 100% (90/90)<br>(95% C.I.: 95.9–100%)   |
| Negatives                                            | 0% (0/119)<br>(95% C.I.: 0–3.1%)        | 0.8% (1/119)<br>(95% C.I.: 0.1–4.6%)    | 0% (0/119)<br>(95% C.I.: 0–3.1%)        | 0.3% (1/357)<br>(95% C.I.: 0–1.6%)      |

| Reproducibility Results – Percent of Flu B Positives |                                         |                                         |                                        |                                         |
|------------------------------------------------------|-----------------------------------------|-----------------------------------------|----------------------------------------|-----------------------------------------|
| Sample                                               | Site 1                                  | Site 2                                  | Site 3                                 | Total                                   |
| High negative B                                      | 0% (0/30)<br>(95% C.I.: 0–11.3%)        | 3.3% (1/30)<br>(95% C.I.: 0.6–16.7%)    | 26.7% (8/30)<br>(95% C.I.: 14.2–44.4%) | 10% (9/90)<br>(95% C.I.: 5.4–17.9%)     |
| Low positive B                                       | 73.3% (22/30)<br>(95% C.I.: 55.6–85.8%) | 90% (27/30)<br>(95% C.I.: 74.4–96.5%)   | 90% (27/30)<br>(95% C.I.: 74.4–96.5%)  | 84.4% (76/90)<br>(95% C.I.: 75.6–90.5%) |
| Moderate positive B                                  | 100% (29/29)<br>(95% C.I.: 88.3–100%)   | 96.6% (28/29)<br>(95% C.I.: 82.8–99.4%) | 100% (29/29)<br>(95% C.I.: 88.3–100%)  | 98.9% (86/87)<br>(95% C.I.: 93.8–99.8%) |
| Negatives                                            | 0% (0/210)<br>(95% C.I.: 0–1.8%)        | 1.0% (2/210)<br>(95% C.I.: 0.3–3.4%)    | 0% (0/210)<br>(95% C.I.: 0–1.8%)       | 0.3% (2/630)<br>(95% C.I.: 0.1–1.2%)    |

### Analytical Studies

#### Analytical Sensitivity (Limit of Detection)

The limit of detection (LOD) for the **BD Veritor** System for Rapid Detection of Flu A+B test was established for a total of 8 influenza strains: 5 influenza A and 3 influenza B. The LOD for each strain represents the lowest concentration producing a positivity rate of ≥95% based on testing 20 to 60 replicates.

| Type | Influenza Viral Strain   | Calculated LOD (TCID <sub>50</sub> /mL) | Calculated LOD (EID <sub>50</sub> /mL) | No. Positive / Total | % Positive |
|------|--------------------------|-----------------------------------------|----------------------------------------|----------------------|------------|
| A    | A/Brisbane/10/2007 H3N2  | 7.27 x 10 <sup>2</sup>                  | N/A                                    | 57/60                | 95%        |
| A    | A/Brisbane/59/2007 H1N1  | 3.30 x 10 <sup>2</sup>                  | N/A                                    | 57/60                | 95%        |
| A    | A/California/7/2009 H1N1 | 5.00 x 10 <sup>3</sup>                  | N/A                                    | 57/60                | 95%        |
| A    | A/Victoria/3/75 H3N2     | 3.11 x 10 <sup>3</sup>                  | N/A                                    | 59/60                | 98.3%      |
| A    | A/Anhui/1/2013 H7N9      | N/A                                     | 5.42 x 10 <sup>6</sup>                 | 59/60                | 98.3%      |
| B    | B/Brisbane/60/2008       | 7.42 x 10 <sup>3</sup>                  | N/A                                    | 58/60                | 96.7%      |
| B    | B/Florida/4/2006         | 1.30 x 10 <sup>3</sup>                  | N/A                                    | 58/60                | 96.7%      |
| B    | B/Lee/40                 | 4.44 x 10 <sup>4</sup>                  | N/A                                    | 20/20                | 100%       |

TCID<sub>50</sub>/mL = 50% Tissue Culture Infectious Dose

EID<sub>50</sub>/mL = 50% Egg Infectious Dose

## Strain Reactivity with Influenza A and B Viruses

The **BD Veritor** System for Rapid Detection of Flu A+B test was evaluated using a panel of influenza strains. Each strain was diluted and tested in triplicate until a point where not all of the replicates were positive. The dilution prior to that is provided in the table below as a minimal detected concentration. All influenza A strains showed positive Flu A test results and negative Flu B test results. Conversely, all of the influenza B strains showed positive Flu B test results and negative Flu A test results.

Although this test has been shown to detect novel avian influenza A (H7N9) and H3N2v cultured viruses the performance characteristics of this device with clinical specimens that are positive for novel avian influenza A (H7N9) and H3N2v influenza viruses has not been established. The **BD Veritor** System Flu A+B assay can distinguish between influenza A and B viruses, but it cannot differentiate influenza A subtypes.

| Strain                                   | Subtype | Minimal Detected Concentration                 |
|------------------------------------------|---------|------------------------------------------------|
| A/Brisbane/59/2007                       | H1N1    | 3.3 x 10 <sup>2</sup> TCID <sub>50</sub> /mL*  |
| A/California/7/2009                      | H1N1    | 5.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL*  |
| A/Denver/1/57                            | H1N1    | 4.45 x 10 <sup>4</sup> CEID <sub>50</sub> /mL  |
| A/FM/1/47                                | H1N1    | 7.91 x 10 <sup>4</sup> CEID <sub>50</sub> /mL  |
| A/Fujian-Gulou/1896/2009                 | H1N1    | 4.5 x 10 <sup>5</sup> CEID <sub>50</sub> /mL   |
| A/Mal/302/54                             | H1N1    | 2.22 x 10 <sup>5</sup> CEID <sub>50</sub> /mL  |
| A/New Caledonia/20/1999                  | H1N1    | 2.5 x 10 <sup>3</sup> TCID <sub>50</sub> /mL   |
| A/New Jersey/8/76                        | H1N1    | 1.58 x 10 <sup>3</sup> CEID <sub>50</sub> /mL  |
| A/NWS/33                                 | H1N1    | 1.58 x 10 <sup>4</sup> CEID <sub>50</sub> /mL  |
| A/PR/8/34                                | H1N1    | 6.31 x 10 <sup>2</sup> TCID <sub>50</sub> /mL  |
| A/Solomon Island/03/2006                 | H1N1    | 2.5 x 10 <sup>3</sup> TCID <sub>50</sub> /mL   |
| A/Washington/24/2012                     | H1N1    | 3.16 x 10 <sup>4</sup> EID <sub>50</sub> /mL   |
| A/Weiss/43                               | H1N1    | 7.03 x 10 <sup>6</sup> CEID <sub>50</sub> /mL  |
| A/WS/33                                  | H1N1    | 7.91 x 10 <sup>2</sup> CEID <sub>50</sub> /mL  |
| A/Aichi/2/68                             | H3N2    | 7.91 x 10 <sup>3</sup> CEID <sub>50</sub> /mL  |
| A/Brisbane/10/2007                       | H3N2    | 7.27 x 10 <sup>2</sup> TCID <sub>50</sub> /mL* |
| A/California/02/2014                     | H3N2    | 1.45 x 10 <sup>2</sup> TCID <sub>50</sub> /mL  |
| A/Hong Kong/8/68                         | H3N2    | 8.89 x 10 <sup>4</sup> CEID <sub>50</sub> /mL  |
| A/Moscow/10/99                           | H3N2    | 5.8 x 10 <sup>6</sup> TCID <sub>50</sub> /mL   |
| A/Perth/16/2009                          | H3N2    | 1.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL   |
| A/Port Chalmers/1/73                     | H3N2    | 3.95 x 10 <sup>4</sup> CEID <sub>50</sub> /mL  |
| A/Switzerland/9715293/2013               | H3N2    | 3.25 x 10 <sup>2</sup> TCID <sub>50</sub> /mL  |
| A/Texas/50/2012                          | H3N2    | 1.75 x 10 <sup>3</sup> TCID <sub>50</sub> /mL  |
| A/Wisconsin/67/2005                      | H3N2    | 2.5 x 10 <sup>5</sup> TCID <sub>50</sub> /mL   |
| A/Victoria/3/75                          | H3N2    | 3.11 x 10 <sup>3</sup> TCID <sub>50</sub> /mL* |
| A/Indiana/08/2011                        | H3N2v   | 1 x 10 <sup>4</sup> TCID <sub>50</sub> /mL     |
| A/Indiana/10/2011                        | H3N2v   | 7.9 x 10 <sup>6</sup> CEID <sub>50</sub> /mL   |
| A/Kansas/13/2009                         | H3N2v   | 1.0 x 10 <sup>3</sup> TCID <sub>50</sub> /mL   |
| A/Minnesota/11/2010                      | H3N2v   | 7.9 x 10 <sup>5</sup> CEID <sub>50</sub> /mL   |
| A/Pennsylvania/14/2010                   | H3N2v   | 1.26 x 10 <sup>6</sup> CEID <sub>50</sub> /mL  |
| A/West Virginia/06/2011                  | H3N2v   | 7.9 x 10 <sup>3</sup> TCID <sub>50</sub> /mL   |
| A/Anhui/01/2005                          | H5N1    | 0.512 HA                                       |
| A/Vietnam/1203/2004                      | H5N1    | 0.512 HA                                       |
| A/Northern Pintail/Washington/40964/2014 | H5N2    | 6.28 x 10 <sup>5</sup> EID <sub>50</sub> /mL   |
| A/Pheasant/New Jersey/1355/1998          | H5N2    | 0.256 HA                                       |
| A/Gyrfalcon/Washington/41088-6/2014      | H5N8    | 1.98 x 10 <sup>6</sup> EID <sub>50</sub> /mL   |
| A/Mallard/Netherlands/12/2000            | H7N7    | 0.256 HA                                       |
| A/Anhui/1/2013                           | H7N9    | 5.42 x 10 <sup>6</sup> CEID <sub>50</sub> /mL* |
| A/Chicken/HongKong/G9/1997               | H9N2    | 1.024 HA                                       |

\*Values taken from preceding Analytical Limit of Detection Table

a. EID<sub>50</sub> = 50% Egg Infectious Dose

b. TCID<sub>50</sub> = 50% Tissue Culture Infectious Dose

c. CEID<sub>50</sub> = 50% Chicken Embryo Infectious dose

d. HA = Hemagglutination Assay

| Strain                                         | Minimal Detected Concentration                 |
|------------------------------------------------|------------------------------------------------|
| B/Brazil/178/96                                | 2.32 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Brisbane/33/2008 (Victoria Lineage)          | 2.45 x 10 <sup>5</sup> CEID <sub>50</sub> /mL  |
| B/Brisbane/60/2008                             | 7.42 x 10 <sup>3</sup> TCID <sub>50</sub> /mL* |
| B/Brisbane/72/97                               | 1.00 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Canada/548/99                                | >0.64 HA                                       |
| B/Egypt/393/99                                 | >1.28 HA                                       |
| B/Florida/2/2006                               | 1.08 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Florida/4/2006                               | 1.30 x 10 <sup>3</sup> TCID <sub>50</sub> /mL* |
| B/Fujian/93/97                                 | 3.95 x 10 <sup>5</sup> TCID <sub>50</sub> /mL  |
| B/Fukushima/220/99                             | 9.33 x 10 <sup>2</sup> TCID <sub>50</sub> /mL  |
| B/Guangdong-Liwan/1133/2014 (Yamagata Lineage) | 9.0 x 10 <sup>5</sup> CEID <sub>50</sub> /mL   |
| B/GuangXi/547/98                               | 2.32 x 10 <sup>5</sup> TCID <sub>50</sub> /mL  |
| B/Hawaii/01/97                                 | >6.4 HA                                        |
| B/Hong Kong/5/72                               | 1.11 x 10 <sup>4</sup> CEID <sub>50</sub> /mL  |
| B/Hong Kong/219/98                             | >1 HA                                          |
| B/Hong Kong/259/2010 (Victoria Lineage)        | 1.35 x 10 <sup>6</sup> CEID <sub>50</sub> /mL  |
| B/Jiangsu/10/2003                              | 1.16 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Johannesburg/5/99                            | 3.95 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Lee/40                                       | 4.44 x 10 <sup>4</sup> CEID <sub>50</sub> /mL* |
| B/Lisbon/03/96                                 | >0.08 HA                                       |
| B/Malaysia/2506/2004                           | 5.0 x 10 <sup>4</sup> TCID <sub>50</sub> /mL   |
| B/Maryland/1/59                                | 3.51 x 10 <sup>2</sup> CEID <sub>50</sub> /mL  |
| B/Massachusetts/2/2012 (Yamagata Lineage)      | 1.25 x 10 <sup>6</sup> CEID <sub>50</sub> /mL  |
| B/Mass/3/66                                    | 1.58 x 10 <sup>5</sup> CEID <sub>50</sub> /mL  |
| B/Montana/5/2012                               | 3.14 x 10 <sup>5</sup> EID <sub>50</sub> /mL   |
| B/Ohio/11/96                                   | >0.16 HA                                       |
| B/Ohio/1/05                                    | 1.34 x 10 <sup>5</sup> TCID <sub>50</sub> /mL  |
| B/Phuket/3073/2013                             | 6.08 x 10 <sup>3</sup> TCID <sub>50</sub> /mL  |
| B/Puerto Mont/10427/98                         | 0.02 HA                                        |
| B/Russia/69                                    | 3.9 x 10 <sup>2</sup> TCID <sub>50</sub> /mL   |
| B/Shandong/7/97                                | 1.58 x 10 <sup>6</sup> TCID <sub>50</sub> /mL  |
| B/Shanghai/04/97                               | 1.58 x 10 <sup>5</sup> TCID <sub>50</sub> /mL  |
| B/Shenzhen/135/97                              | 3.16 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Sichuan/116/96                               | 0.016 HA                                       |
| B/Taiwan/2/62                                  | 2.81 x 10 <sup>2</sup> CEID <sub>50</sub> /mL  |
| B/Texas/06/2011 (Yamagata Lineage)             | 6.2 x 10 <sup>5</sup> CEID <sub>50</sub> /mL   |
| B/Texas/02/2013 (Victoria Lineage)             | 2.75 x 10 <sup>4</sup> CEID <sub>50</sub> /mL  |
| B/Utah/09/2014 (Yamagata Lineage)              | 6.3 x 10 <sup>3</sup> CEID <sub>50</sub> /mL   |
| B/Victoria/504/00                              | 4.64 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |
| B/Wisconsin/01/2010 (Yamagata Lineage)         | 7.0 x 10 <sup>2</sup> CEID <sub>50</sub> /mL   |
| B/Yamagata/16/88                               | 9.75 x 10 <sup>3</sup> TCID <sub>50</sub> /mL  |
| B/Yamanashi/166/98                             | 4.88 x 10 <sup>4</sup> TCID <sub>50</sub> /mL  |

\*Values taken from preceding Analytical Limit of Detection Table

a. EID<sub>50</sub> = 50% Egg Infectious Dose

b. TCID<sub>50</sub> = 50% Tissue Culture Infectious Dose

c. CEID<sub>50</sub> = 50% Chicken Embryo Infectious dose

d. HA = Hemagglutination Assay

On January 12, 2017 the US FDA published an order reclassifying antigen based rapid influenza virus detection systems (RIDTs) intended to detect influenza virus antigen directly from clinical specimens from Class I to Class II with special controls. One of these special controls is a requirement for annual analytical reactivity testing of contemporary influenza strains identified by the Centers for Disease control (CDC) using a standardized dilution protocol. Results obtained using the **BD Veritor** System are available for viewing at [www.BD.com/veritorsystem](http://www.BD.com/veritorsystem).

### Analytical Specificity (Cross-reactivity)

The **BD Veritor** System for Rapid Detection of Flu A+B test was evaluated with a total of 51 microorganisms. The 37 bacteria and yeast were tested at a target concentration of approximately  $10^7$  CFU/mL (CFU – Colony Forming Units) with the exception of *Staphylococcus aureus*, which was tested at a final concentration of  $10^6$  CFU/mL. The 14 viruses were evaluated at concentrations of  $10^3$  to  $10^{10}$  TCID<sub>50</sub>/mL. Of the 51 microorganisms tested, none showed cross-reactivity in either the Flu A or Flu B tests.

|                                    |                                           |
|------------------------------------|-------------------------------------------|
| <i>Bacteroides fragilis</i>        | <i>Neisseria sp. (Neisseria perflaus)</i> |
| <i>Bordetella pertussis</i>        | <i>Neisseria subflava</i>                 |
| <i>Candida albicans</i>            | <i>Peptostreptococcus anaerobius</i>      |
| <i>Chlamydia pneumoniae</i>        | <i>Porphyromonas asaccharolyticus</i>     |
| <i>Corynebacterium diphtherium</i> | <i>Prevotella oralis</i>                  |
| <i>Escherichia coli</i>            | <i>Propionibacterium acnes</i>            |
| <i>Fusobacterium nucleatum</i>     | <i>Proteus mirabilis</i>                  |
| <i>Haemophilus influenzae</i>      | <i>Pseudomonas aeruginosa</i>             |
| <i>Haemophilus parainfluenzae</i>  | <i>Serratia marcescens</i>                |
| <i>Kingella kingae</i>             | <i>Staphylococcus aureus</i>              |
| <i>Klebsiella pneumoniae</i>       | <i>Staphylococcus epidermidis</i>         |
| <i>Lactobacillus sp.</i>           | <i>Streptococcus mutans</i>               |
| <i>Legionella sp.</i>              | <i>Streptococcus pneumoniae</i>           |
| <i>Moraxella catarrhalis</i>       | <i>Streptococcus pyogenes</i>             |
| <i>Mycobacterium tuberculosis</i>  | <i>Streptococcus sp. Group C</i>          |
| <i>Mycoplasma pneumoniae</i>       | <i>Streptococcus sp. Group G</i>          |
| <i>Neisseria gonorrhoeae</i>       | <i>Streptococcus salivarius</i>           |
| <i>Neisseria meningitidis</i>      | <i>Veillonella parvula</i>                |
| <i>Neisseria mucosa</i>            |                                           |

|                                    |
|------------------------------------|
| Adenovirus, type 1                 |
| Adenovirus, type 7                 |
| Cytomegalovirus                    |
| Enterovirus                        |
| Epstein Barr Virus                 |
| HSV Type 1                         |
| Human Coronavirus OC43             |
| Human Coronavirus 2229E            |
| Human metapneumovirus (HMPV-27 A2) |
| Human Parainfluenza                |
| Measles virus                      |
| Mumps virus                        |
| Respiratory syncytial virus        |
| Rhinovirus                         |
|                                    |
|                                    |
|                                    |
|                                    |

### Interfering Substances

Various substances were evaluated with the **BD Veritor** System for Rapid Detection of Flu A+B test. These substances included whole blood (2%) and various medications. No interference was noted with this assay for any of the substances tested.

| Substance                | Concentration | Substance                    | Concentration |
|--------------------------|---------------|------------------------------|---------------|
| 4-Acetamidophenol        | 10 mg/mL      | Homeopathic Allergy Medicine | 10 mg/mL      |
| Acetylsalicylic acid     | 20 mg/mL      | Ibuprofen                    | 10 mg/mL      |
| Albuterol                | 0.083 mg/mL   | Loratidine                   | 100 ng/mL     |
| Amantadine Hydrochloride | 500 ng/mL     | Menthol Throat Lozenges      | 10 mg/mL      |
| Ayr Saline Nasal Gel     | 10 mg/mL      | Mometasone                   | 500 ng/mL     |
| Beclomethasone           | 500 ng/mL     | Mupirocin                    | 500 ng/mL     |
| Budesonide               | 500 ng/mL     | Oseltamivir                  | 500 ng/mL     |
| Chlorpheniramine maleate | 5 mg/mL       | Oxymetazoline                | 0.05 mg/mL    |
| Dexamethasone            | 10 mg/mL      | Phenylephrine                | 1 mg/mL       |
| Dextromethorphan         | 10 mg/mL      | Pseudoephedrine HCl          | 20 mg/mL      |
| Diphenhydramine HCl      | 5 mg/mL       | Purified Mucin Protein       | 1 mg/mL       |
| Fexofenadine             | 500 ng/mL     | Ribavirin                    | 500 ng/mL     |
| FluMist                  | 1%            | Rimantadine                  | 500 ng/mL     |
| Flunisolide              | 500 ng/mL     | Three OTC mouthwashes        | 5%            |
| Fluticasone              | 500 ng/mL     | Tobramycin                   | 500 ng/mL     |
| Four OTC nasal sprays    | 10%           | Triamcinolone                | 500 ng/mL     |
| Four OTC throat drops    | 25%           | Whole Blood                  | 2%            |
| Guaiacol Glyceryl Ether  | 20 mg/mL      | Zanamivir                    | 1 mg/mL       |

Of the 44 substances tested in this study, none exhibited interfering reactions when tested with influenza A and influenza B positive samples. Based on the data, the substances tested at the indicated concentration levels did not interfere with the **BD Veritor** System for Rapid Detection of Flu A+B test.

CLIA WAIVER STUDY

As part of a larger prospective study, as described in the Performance Characteristics section above, the accuracy of the **BD Veritor** System for Rapid Detection of Flu A+B test was evaluated at five CLIA waived testing sites. A total of 31 operators representative of CLIA waived sites (intended users) participated in the study. No training on the use of the test was provided. The study included 515 nasal/nasopharyngeal swabs prospectively collected and 80 retrospective archived specimens. The **BD Veritor** System results were compared with results obtained by an FDA cleared molecular influenza A and B assay, the comparator method. Three specimens were excluded due to **BD Veritor** invalid results. The invalid rate was 0.5% (3/598) with 95% C.I: 0.2% to 1.5%.

The positive percent agreement (PPA) and the negative percent agreement (NPA) between the **BD Veritor** results and the comparator method are presented in the tables below (refer to Performance Characteristics section for definition of terms).

| INFLUENZA A                                                                                                     |                 |                         |                 |                         |
|-----------------------------------------------------------------------------------------------------------------|-----------------|-------------------------|-----------------|-------------------------|
| Positive Percent Agreement and Negative Percent Agreement of BD Veritor Flu A+B Test with the Comparator Method |                 |                         |                 |                         |
| Total Number of Samples                                                                                         | PPA             | 95% Confidence Interval | NPA             | 95% Confidence Interval |
| 595                                                                                                             | 82.1% (151/184) | (75.9%, 86.9%)          | 98.1% (403/411) | (96.2%, 99.0%)          |

| INFLUENZA B                                                                                                     |                 |                         |                 |                         |
|-----------------------------------------------------------------------------------------------------------------|-----------------|-------------------------|-----------------|-------------------------|
| Positive Percent Agreement and Negative Percent Agreement of BD Veritor Flu A+B Test with the Comparator Method |                 |                         |                 |                         |
| Total Number of Samples                                                                                         | PPA             | 95% Confidence Interval | NPA             | 95% Confidence Interval |
| 595                                                                                                             | 79.7% (102/128) | (71.9%, 85.7%)          | 99.4% (464/467) | (98.1%, 99.8%)          |

Another study was designed to assess the capability of untrained users to test weakly reactive samples and deliver results with accuracy. The **BD Veritor** System for Rapid Detection of Flu A+B assay was evaluated at three non-laboratory CLIA waived sites using panels of simulated swab samples including two weak positives near the assay cutoff and one negative sample. The positive swab samples were formulated at two levels: a "low positive" sample targeted at the assay limit of detection; and a "high negative" sample targeted just below the assay limit of detection. The panels included two strains of Flu A viruses (A/California/7/2009 and A/Victoria 3/75) and one Flu B virus (B/Lee/40). The swab samples were randomized and masked with respect to their identity. There were two intended users at each of the CLIA waived sites (six operators in total) and each site tested the panel on each of 10 days. The same panels of simulated swab samples were also tested at three clinical laboratory sites as controls. The performance of the **BD Veritor** System with samples near the assay cutoff was acceptable when used by intended users.

The tables below show performance of the test with samples near the cutoff of the assay for influenza A and influenza B in the hands of untrained intended users (across all sites).

| Influenza A Viral Strains              |                          |                         |
|----------------------------------------|--------------------------|-------------------------|
| Sample Type                            | Untrained Intended Users |                         |
|                                        | Percent Detection        | 95% Confidence Interval |
| High Negative A/California/7/2009 H1N1 | 6.7% (4/60)              | (2.6%, 15.9%)           |
| Low Positive A/California/7/2009 H1N1  | 81.7% (49/60)            | (70.1%, 89.4%)          |
| High Negative A/Victoria 3/75 H3N2     | 6.7% (4/60)              | (2.6%, 15.9%)           |
| Low Positive A/Victoria 3/75 H3N2      | 80.0% (48/60)            | (68.2%, 88.2%)          |
| Negative                               | 0% (0/118)*              | (0%, 3.2%)              |

\*Two (2) samples were excluded from the analysis due to errors in data recording.

| Influenza B Viral Strain |                          |                         |
|--------------------------|--------------------------|-------------------------|
|                          | Untrained Intended Users |                         |
| Sample Type              | Percent Detection        | 95% Confidence Interval |
| High Negative B/Lee/40   | 11.7% (7/60)             | (5.8%, 22.2%)           |
| Low Positive B/Lee/40    | 72.4% (42/58)*           | (59.8%, 82.2%)          |
| Negative                 | 0% (0/240)               | (0%, 1.6%)              |

\*Two (2) samples were excluded from the analysis due to errors in data recording.

Using risk analysis as a guide, analytical flex studies were conducted. The studies demonstrated that the test is insensitive to stresses of environmental conditions and potential user errors.

In support of the CLIA waiver, an additional reactivity study was performed at an independent laboratory to demonstrate reactivity of the **BD Veritor** System for the Rapid Detection of Flu A+B with a broad range of contemporary influenza A and influenza B viruses. The **BD Veritor** System yielded positive results with all 18 influenza A viruses and 7 influenza B viruses included in the test panel at acceptable viral load levels.

### Technical Support

For questions, or to report a problem, please call Technical Support at 1.800.638.8663. Test system problems may also be reported to the FDA using the MedWatch reporting system (phone: 1-800 FDA-1088; fax: 1-800 FDA-1078; or <http://www.fda.gov/medwatch>).

### AVAILABILITY

#### Cat. No. Description

- 256045 **BD Veritor™** System for Rapid Detection of Flu A+B, 30 tests
- 256051 **BD Veritor™** System Flu A+B Control Swab Set, 10 pairs of swabs
- 220252 COPAN Flexible Minitip Flocked Swab, 100 swabs
- 256055 **BD Veritor™** System Reader
- 256066 **BD Veritor™** Plus Analyzer
- 256067 **BD Veritor™** InfoSync Module
- 256068 **BD Veritor™** InfoScan Module
- 443907 USB Printer Cable for **BD Veritor™** Analyzer
- 256038 **BD Veritor™** System for Rapid Detection of Respiratory Syncytial Virus (RSV)

### REFERENCES

1. Simonsen L., Fukuda K, Schonberger LB, Cox NJ. Impact of influenza epidemics on hospitalizations. *J. Infect. Dis.* 2000; **181**:831–7
2. Thompson WW, Shay DK, Weintraub E, et al. Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA* 2003; **289**:179–86
3. Treanor, J.J., Hayden, F.G., Vrooman, P.S., et al. 2000. Efficacy and safety of the oral neuraminidase inhibitor oseltamivir in treating acute influenza: a randomized controlled trial. *JAMA.* **283**:1016–1024.
4. Kaiser, L., Couch, R.B., Galasso, G.J., Glezen, W.P., Webster, R.G., Wright, P.F., and Hayden, F.G. 1999. First international symposium on influenza and other respiratory viruses: summary and overview Kapalua, Maui, Hawaii, December 4–6, 1998. *Antiviral Res.*, **42**:149–176
5. Cox, N.J., and Bender, C.A. 1995. The molecular epidemiology of influenza viruses. *Virology*, **6**:359–370.
6. Todd, S.J., Minnich, L., and Waner, J.L. 1995. Comparison of rapid immunofluorescence procedure with TestPack RSV and Directigen Flu A for diagnosis of respiratory syncytial virus and influenza A virus. *J. Clin. Microbiol.***33**:1650–1651.
7. Harris, P.O. 1989. Clinical relevance and efficient detection of seven major respiratory viruses. *ACL.* p. 15–19.
8. McElhaney, J.E., Gravenstein, S., Krause, P., Hooton, J.W., Upshaw, C.M., and Drinka, P. 1998. Assessment of markers of the cell-mediated immune response after influenza virus infection in frail older adults. *Clin. Diag. Lab. Immunol.* **5**:840–844.
9. Fan, J., Henrickson, K.J., and Savatski, L.L. 1998. Rapid simultaneous diagnosis of infections with respiratory syncytial viruses A and B, influenza viruses A and B, and human parainfluenza virus types 1, 2, and 3 by multiplex quantitative reverse transcription-polymerase chain reaction-hybridization assay (hexaplex). *Clin. Infect. Disease* **26**:1397–1402.
10. Wright, K.E., Wilson, G.A.R., Novosad, D., Dimock, C., Tan, D., and Weber, J.M. 1995. Typing and subtyping of influenza viruses in clinical samples by PCR. *J. Clin. Microbiol.***33**:1180–1184.
11. Kendal, A.P. 1985. Influenza Viruses. p. 341–357. *Laboratory Diagnosis of Viral Infections*, *In* H. Lennette, (ed.) Marcel Dekker, Inc., New York.
12. McQuillen, J., Madeley, C.R., and Kendal, A.P. 1985. Monoclonal antibodies for the rapid diagnosis of influenza A and B virus infections by immunofluorescence. *Lancet*.*ii*: 911–914.

13. Guenther, S.H., and Linnemann, C.C., Jr. 1988. Indirect immunofluorescence assay for rapid diagnosis of influenza virus. *Laboratory Medicine*. 19:581–583
14. Minnick, L.L., and Ray, C.G. 1986. Early testing of cell cultures for detection of hemadsorbing viruses. *J. Clin. Microbiol.* 25:421–422.
15. Schmidt, N.J., Ota, M., Gallo, D., and Fox, V.L. 1982. Monoclonal antibodies for rapid, strain specific identification of influenza virus isolates. *J. Clin. Microbiol.* 16:763–765.
16. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
17. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53–80.
18. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S Government Printing Office, Washington, D.C.
19. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Office Journal L262*, 17/10/2000, p.021–0045.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or [www.bd.com](http://www.bd.com).

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# Change History

| Revision | Sections | Change Summary                                                                                                                                                                                                                                                                                                                             |
|----------|----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 14       | All      | Changes required after FDA review of K180438 and CR 180086.<br>Correcting verbiage throughout, adding Caution bullets to Sections A, B, C and D.<br>Add new Table 1 and renumber tables.<br>Rename Analyzer “the BD Veritor Plus Analyzer” throughout.<br>Reorganize Catalog numbers in Availability section.<br>Add Change History Table. |

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For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrense / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiiri / Limite inférieure de température / Najniža dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураның төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrense / Dolna granica temperatury / Limite minimo de temperatura / Limitā minimā de temperaturā / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgrens / Sıcaklık alt sınırı / Мінімальна температура / 温度下限

Control / Контроль / Kontrola / Kontrol / Kontrolle / Márturas / Contrôle / Controllo / Бақылау / 컨트롤 / Kontrolé / Kontrolle / Controle / Контроль / kontroll / Kontrolle / 对照

Positive control / Положительный контроль / Pozitívni kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positiveve kontroll / Contrôle positif / Pozitivna kontrola / Pozitiv kontrol / Controllo positivo / Оң бақылау / 양성 컨트롤 / Teigiama kontrolė / Pozitívna kontrola / Positieve controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивный контроль / 阳性对照试剂

Negative control / Отрицательный контроль / Negatívni kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negativeve kontroll / Contrôle négatif / Negativna kontrola / Negativ kontrol / Controllo negativo / Негативтік бақылау / 음성 컨트롤 / Neigiama kontrolė / Negatívna kontrola / Negative controle / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatív kontrol / Негативный контроль / 阴性对照试剂

Method of sterilization: ethylene oxide / Метод на стерилизация: этиленов оксид / Způsob sterilizace: etylenoxid / Steriliseringmethode: ethylenoxid / Sterilisationsmethode: Ethylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Método de esterilización: óxido de etileno / Steriliseringmetode: etüleenoksiid / Méthode de stérilisation: oxyde d'éthylène / Metoda sterilizacji: etilen oksid / Sterilázás módszere: etilén-oxid / Metodo di sterilizzazione: ossido di etilene / Стерилизація едіці — этилен тотығы / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizēšanas metode: etilēnoksis / Gesteriliseerd met behulp van ethyleenoxide / Steriliseringsmethode: etylenoksid / Metoda sterylizacji: tienek etyli / Método de esterilização: óxido de etileno / Metodă de sterilizare: oxid de etilenă / Метод стерилизації: этиленоксид / Metóda sterilizácie: etylenoxid / Metoda sterilizacije: etilen oksid / Steriliseringsmetod: etenoxid / Sterilizasyon yöntemi: etilen oksit / Метод стерилизації: этиленоксидом / 灭菌方法: 环氧乙烷

Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringsmetode: bestraling / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseerimisselä: kiirgus / Méthode de stérilisation : irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Стерилизация әдісі – сәулелік / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringsmetode: bestraling / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metoda de sterilizare: iradiere / Metod sterilizācijas: apstarošana / Metóda sterilizácie: ožarovanie / Metoda sterilizacije: ozračevanje / Steriliseringsmetod: strålning / Sterilizasyon yöntemi: ırradyasyon / Metod sterilizacji: opromieniam / 灭菌方法: 辐射

Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Biologilised riskid / Risques biologiques / Biološki rizik / Biológiallag veszélyes / Rischio biologico / Биологический риск / 생물학적 위험 / Biologinis pavojus / Biologiskie riski / Biologisch risico / Biologisk risiko / Zagrożenia biologiczne / Perigo biológico / Riscuri biologice / Биологическая опасность / Biologické riziko / Biološki rizici / Biologisk risk / Biyolojik Riskler / Биологічна небезпека / 生物学风险

Caution, consult accompanying documents / Внимание, направте справка в придружаващите документи / Pozor! Prostudujte si příložený dokumentaci! / Forsigtig, se ledsagende dokumenter / Achtung, Begleitdokumente beachten / Προσοχή, συμπληρωθείτε το συνοδευτικό έγγραφο / Precaución, consultar la documentación adjunta / Ettevaatust! Lugeka kaasnevat dokumentatsiooni / Attention, consulter les documents joints / Upozorenje, koristi prateću dokumentaciju / Figyelem! Olvassa el a mellékelt tájékoztató! / Attenzione: consultare la documentazione allegata / Абайлаңыз, тиісті құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Dmesio, zürékite pridédamus dokumentus / Piesardziba, skatit pavaddokumentus / Voorzichtig, raadpleeg bijgevoegde documenten / Forsiktig, se vedlagt dokumentasjon / Należy zapoznać się z dołączonymi dokumentami / Cuidado, consulte a documentação fornecida / Atenție, consultați documentele însoțitoare / Внимание: см. прилагаемую документацию / Vstraha, pozri sprievodné dokumenty / Paźnia! Pogledajcie prilożona dokumenta / Obs! Se medföljande dokumentation / Dikkat, birlikte verilen belgeler başvurun / Увага: див. супутню документацию / 小心, 请参阅附带文档。

Upper limit of temperature / Горен лимит на температурата / Horní hranice teploty / Øvre temperaturgrænse / Temperaturobergrenze / Ανώτερο όριο θερμοκρασίας / Limite superior de temperatura / Ülemine temperatuuripiir / Limite supérieure de température / Gornja dozvoljena temperatura / Felső hőmérsékleti határ / Limite superiore di temperatura / Температураның рұқсат етілген жоғарғы шегі / 상한 온도 / Aukščiausia laikymo temperatūra / Augšējā temperatūras robeža / Hoogste temperatuurilimiet / Øvre temperaturgrense / Górná granica temperatury / Limite máximo de temperatura / Limită maximă de temperatură / Верхний предел температуры / Horná hranica teploty / Gornja granica temperature / Øvre temperaturgräns / Sıcaklık üst sınırı / Максимальна температура / 温度上限

Keep dry / Παзете сухо / Skladujte v suchém prostředí / Opbeveares tørt / Trocklagern / Φυλάξτε το στεγνό / Mantener seco / Hoida kuivas / Conserver au sec / Držati na suhom / Száraz helyen tartandó / Tenere all'asciutto / Құрқак күйінде ұста / A se feri de umezala / Не допускать попадания влаги / Uchovávaťe v suchu / Držite na suvom mestu / Förvaras torr / Kuru bir şekilde muhafaza edin / Берегти від вологи / 请保持干燥

Collection time / Време на събиране / Čas odběru / Opsamlingstidspunkt / Entnahmeuhrzeit / Ώρα συλλογής / Hora de recogida / Kogumisaeg / Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жинау уақыты / 수집 시간 / Paėmimo laikas / Savākšanas laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de coleta / Ora colectării / Время сбора / Doba odboru / Vreme prikupljanja / Uppsamlingstid / Toplama zamanı / Час забору / 采集时间

Peel / Обенете / Otevěte zde / Åbn / Abziehen / Αποκολλήστε / Desprendre / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Устігні кабатын алып таста / 벗기 / Plésti čia / Aftitmet / Schillen / Trek av / Oderwać / Destacar / Se dezlipeste / Отклеить / Odtrhně / Oluštiti / Dra isär / Ayırma / Відкрити / 撕下

Perforation / Перфорация / Perforace / Perforering / Διείρση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Тесик тесу / 찢히선 / Perforacija / Perforácia / Perforate / Perforacja / Perforação / Perforare / Перфорация / Perforácia / Perforasjon / Перфорация / 穿孔

Do not use if package damaged / Не използвайте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inhal beschädigter Packungnicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιώ. / No usar si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használja, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Erep paket бүзүлган болса, пайдаланба / 패키지가 손상된 경우 사용 금지 / Jei pakuotė pažeista, nenaudoti / Nelietot, ja iepakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не использовать при повреждении упаковки / Nepoužívejte, ak je obal poškodený / Ne koristite ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包装破损, 请勿使用

Keep away from heat / Παзете от топлина / Nevystavujte přílišnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valgusest / Protéger de la chaleur / Držati dalje od izvora topline / Ovjia a melegtől / Tenere lontano dal calore / Салқын жерде сақта / 열을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no karstuma / Beschermen tegen warmte / Må ikke utsettes for varme / Przechowywać z dala od źródła ciepła / Manter ao abrigo do calor / A se feri de căldură / Не нагревать / Uchovávať mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Iszidan uzak tutun / Берегти від дії тепла / 请远离热源

Cut / Срежете / Odstřihněte / Klip / Schneiden / Κόψτε / Cortar / Découper / Reži / Vágja ki / Tagliare / Кесіңіз / 잘라내기 / Kirpti / Nogrietz / Knippen / Kutt / Odciać / Cortar / Decupați / Отрезать / Odstrihnite / Iseći / Klipp / Kesme / Rozpaziti / 剪下

Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmedatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumiskuupäev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинаған тізбекүні / 수집 날짜 / Paėmimo data / Savākšanas datums / Verzameldatum / Dato prøvetaking / Data pobrania / Data de coleta / Data colectării / Дата сбора / Dátum odboru / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期

  $\mu\text{L}/\text{test}$  /  $\mu\text{L}/\text{тест}$  /  $\mu\text{L}/\text{Test}$  /  $\mu\text{L}/\text{εξέταση}$  /  $\mu\text{L}/\text{prueba}$  /  $\mu\text{L}/\text{teszt}$  /  $\mu\text{L}/\text{테스트}$  /  $\text{mkl}/\text{тест}$  /  $\mu\text{L}/\text{tyrimas}$  /  $\mu\text{L}/\text{pārbaude}$  /  $\mu\text{L}/\text{teste}$  /  $\text{mkl}/\text{anālis}$  /  $\mu\text{L}/\text{检测}$

 Keep away from light / Пазете от светлина / Nevystavujte světlu / Må ikke udsættes for lys / Vor Licht schützen / Κρατήστε το μακριά από το φως / Mantener alejado de la luz / Hoida eemal valgusest / Conserver à l'abri de la lumière / Držati dalje od svetla / Fény nem érheti / Tenere al riparo dalla luce / Қараңғыланған жерде ұста / 빛을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no gaismas / Niet blootstellen aan zonlicht / Må ikke utsettes for lys / Przechowywać z dala od źródła światła / Manter ao abrigo da luz / Feriți de lumină / Хранить в темноте / Uchovávať mimo dosahu svetla / Držite dalje od svetlosti / Får ej utsättas för ljus / Işıktan uzak tutun / Беретти від дії світла / 请远离光线

$\text{H}_2$   Hydrogen gas generated / Образован водород газ / Možnost úniku plyného vodíku / Frembringer hydrogengas / Wasserstoffgas erzeugt / Δημιουργία αερίου υδρογόνου / Producción de gas de hidrógeno / Vesinikgaasi tekitatud / Produit de l'hydrogène gazeux / Sadržai hydrogen vodik / Hidrogén gázt fejleszt / Produzione di gas idrogeno / Газтекес сутері пайда болды / 수소 가스 생성됨 / Išskiria vandenilio dujas / Rodas ūdegradis / Waterstofgas gegenereerd / Hydrogengass generert / Powoduje powstawanie wodoru / Produção de gás de hidrogénio / Generare gaz de hidrogen / Выделение водорода / Vyrobené použitím vodíka / Osloбаda se vodonik / Genererad vätgas / Açığa çıkan hidrojen gazı / Реакція з виділенням водню / 会产生氢气

 # Patient ID number / ИД номер на пациента / ID pacienta / Patientens ID-nummer / Patienten-ID / Αριθμός αναγνώρισης ασθενούς / Número de ID del paciente / Patsiendi ID / No d'identification du patient / Identifikacijski broj pacijenta / Beteg azonosító száma / Numero ID paziente / Пациенттің идентификациялық нөмірі / 환자 ID 번호 / Paciento identifikavimo numeris / Pacienta ID numurs / Identificatienummer van de patiënt / Patientens ID-nummer / Numer ID pacienta / Número da ID do doente / Număr ID pacient / Идентификационный номер пациента / Identifikačné číslo pacienta / ID broj pacijenta / Patientnummer / Hasta kimlik numarası / Идентификатор пациента / 患者标识号

 Fragile, Handle with Care / Чупливо, Роботете с необходимото внимание. / Křehké. Při manipulaci postupujte opatrně. / Forsiktig, kan gå i stykker. / Zerbrechlich, vorsichtig handhaben. / Εύθραστο. Χειριστείτε το με προσοχή. / Frágil. Manipular con cuidado. / Örn, käsitlege ettevaatlikult. / Fragile. Manipuler avec précaution. / Lomljivo, rukujte pažljivo. / Törékeny! Óvatosan kezelendő. / Fragile, maneggiare con cura. / Сынгыш, абайлап пайдаланыңыз. / 조심 깨지기 쉬운 처리 / Trapu, elkites atsargiai. / Trausls; rūkotiés uzmanigi / Breekbaar, voorzichtig behandelen. / Ømtålig, håndter forsigtig. / Krucha zawartość, przenosić ostrożnie. / Frágil, Manuseie com Cuidado. / Fragil, manipulați cu atenție. / Хрупкое! Обращаться с осторожностью. / Křehké, vyžaduje sa opatrná manipulácia. / Lomljivo - rukujte pažljivo. / Bräckligt. Hantera försiktigt. / Kolay Kırılır, Dikkatli Taşıyın. / Тендітна, звертатися з обережністю / 易碎，小心轻放

$R_x$  Only This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

 Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA

 EC REP

Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

FluMist is a registered trademark of MedImmune, LLC.

Viratrans is a trademark of Trinity Biotech, PLC.

© 2018 BD. BD, the BD Logo and BD Veritor are trademarks of Becton, Dickinson and Company.

# **Insert ends here**

DO NOT PRINT THIS PAGE

## **CLSI document should not be included in print file to vendor.**

**NOTES:**

The CLSI document is to be updated at the same time as the package insert and should be released as a separate PDF document to the web.

The package insert should also be released as a separate PDF document to the web.

The package insert and CLSI document should be released as one document to SAP.

Final art provided to Buyer/print vendor should not include this page or the CLSI portion of this document.

Proofreader: please proof the CLSI document along with the package insert.

**Procedure\* BD    Veritor™ System For Rapid Detection of Flu A+B**

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| Prepared by | Date Adopted | Supersedes Procedure # |
|-------------|--------------|------------------------|
|             |              |                        |

| Review Date | Revision Date | Signature |
|-------------|---------------|-----------|
|             |               |           |
|             |               |           |
|             |               |           |
|             |               |           |

| Distributed to | # of Copies | Distributed to | # of Copies |
|----------------|-------------|----------------|-------------|
|                |             |                |             |
|                |             |                |             |
|                |             |                |             |
|                |             |                |             |

---

\*Any modifications to this document are the sole responsibility of the facility. This “Sample Procedure” is not intended as a substitute for your facility procedure manual, instrument manual, or reagent labeling/package insert. This “Sample Procedure” is intended as a model for use by your facility to meet the needs of your laboratory.

## CLSI For Rapid Detection of Flu A+B

### INTENDED USE

The **BD Veritor** System for Rapid Detection of Flu A+B is a rapid chromatographic immunoassay for the direct and qualitative detection of influenza A and B viral nucleoprotein antigens from nasal and nasopharyngeal swabs of symptomatic patients.

The **BD Veritor** System for Rapid Detection of Flu A+B (also referred to as the **BD Veritor** System and **BD Veritor** System Flu A+B) is a differentiated test, such that influenza A viral antigens can be distinguished from influenza B viral antigens from a single processed sample using a single device. The test is to be used as an aid in the diagnosis of influenza A and B viral infections. A negative test is presumptive and it is recommended that these results be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay. Outside the U.S., a negative test is presumptive and it is recommended that these results be confirmed by viral culture or a molecular assay cleared for diagnostic use in the country of use. FDA has not cleared this device for use outside of the U.S. Negative test results do not preclude influenza viral infection and should not be used as the sole basis for treatment or other patient management decisions. The test is not intended to detect influenza C antigens.

Performance characteristics for influenza A and B were established during January through March of 2011 when influenza viruses A/2009 H1N1, A/H3N2, B/Victoria lineage, and B/Yamagata lineage were the predominant influenza viruses in circulation according to the *Morbidity and Mortality Weekly Report* from the CDC entitled "Update: Influenza Activity—United States, 2010-2011 Season, and Composition of the 2011-2012 Influenza Vaccine." Performance characteristics may vary against other emerging influenza viruses.

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Diagnostic tests available for influenza include rapid immunoassay, immunofluorescence assay, polymerase chain reaction (PCR), serology, and viral culture.<sup>4-11</sup> Immunofluorescence assays entail staining of specimens immobilized on microscope slides using fluorescent-labeled antibodies for observation by fluorescence microscopy.<sup>6,12,13</sup> Culture methods employ initial viral isolation in cell culture, followed by hemadsorption inhibition, immunofluorescence, or neutralization assays to confirm the presence of the influenza virus.<sup>13-15</sup>

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All **BD Veritor** System Flu A+B test devices are interpreted by a **BD Veritor** System Instrument, either a **BD Veritor** Reader or **BD Veritor** Plus Analyzer (the "Analyzer"). When using the **BD Veritor** Plus Analyzer, workflow steps depend on the selected operational mode and the Analyzer configuration settings. In **Analyze Now mode**, the instrument evaluates assay devices after manual timing of their development. In **Walk Away mode**, devices are inserted immediately after application of the specimen, and timing of assay development and analysis is automated. Additionally, connection of the Analyzer to a printer is possible if desired. Additional result documentation capabilities are possible with the integration of the **BD Veritor** InfoScan ("InfoScan") or **BD Veritor** InfoSync ("InfoSync") module. Please refer to the Analyzer *Instructions for Use* for details on how to implement these features. The **BD Veritor** InfoSync is not available in all regions.

### PRINCIPLES OF THE PROCEDURE

The **BD Veritor** System for Rapid Detection of Flu A+B is a qualitative, digital immunoassay for the detection of influenza A and B viral antigens in samples processed from respiratory specimens. When specimens are processed and added to the test device, influenza A or B viral antigens bind to anti-influenza antibodies conjugated to detector particles in the A+B test strip. The antigen-conjugate complex migrates across the test strip to the reaction area and is captured by the line of antibody on the membrane. A positive result for influenza A is determined by the **BD Veritor** System Instrument when antigen-conjugate is deposited at the Test "A" position and the Control "C" position on the **BD Veritor** System Flu A+B assay device. A positive result

for influenza B is determined by the **BD Veritor** System Instrument when antigen-conjugate is deposited at the Test "B" position and the Control "C" position in the **BD Veritor** System Flu A+B assay device. The instrument analyzes and corrects for non-specific binding and detects positives not recognized by the unaided eye to provide an objective digital result.

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## Warnings and Precautions:

### Warning

**H302** Harmful if swallowed. **H402** Harmful to aquatic life. **H412** Harmful to aquatic life with long lasting effects.

**P273** Avoid release to the environment. **P264** Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

1. For *in vitro* Diagnostic Use.
2. Test results are not meant to be visually determined. **All test results must be determined using the BD Veritor System Instrument.**
3. If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to the state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.
4. Pathogenic microorganisms, including hepatitis viruses, Human Immunodeficiency Virus and novel influenza viruses, may be present in clinical specimens. "Standard Precautions"<sup>16-19</sup> and institutional guidelines should be followed in handling, storing and disposing of all specimens and all items contaminated with blood and other body fluids.
5. Dispose of used **BD Veritor** System test devices as biohazardous waste in accordance with federal, state and local requirements.
6. Reagents contain sodium azide, which is harmful if inhaled, swallowed or exposed to skin. Contact with acids produces very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
7. For optimal results, use the flocked swabs provided with the kit for specimen collection.
8. Other than the flocked swabs that are used for specimen collection, kit components should not make contact with the patient.
9. Do not use kit components beyond the expiration date.
10. Do not reuse the device.
11. Do not use the kit if the Control A+/B- swab and Control B+/A- swab do not yield appropriate results.
12. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
13. To avoid erroneous results, swab specimens must be processed as indicated in the assay procedure section.
14. Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.
15. FluMist® is made from attenuated live flu virus and although the concentration tested (1%) was non-interfering, it is possible when tested with higher concentrations that an influenza A and/or influenza B false positive may occur.

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**SPECIMEN COLLECTION**

Acceptable specimens for testing with the **BD Veritor** System Flu A+B test include nasal swabs and nasopharyngeal (NP) swabs. Freshly collected specimens should be processed within 1 hour. It is essential that correct specimen collection and preparation methods be followed. Specimens obtained early in the course of the illness will contain the highest viral titers. Inadequate specimen collection, improper specimen handling and/or transport may yield a false negative result; therefore, specimen collection requires specific training and guidance due to the importance of specimen quality to accurate test results.

**Proper Nasal Swab Sample Collection**


- 1. The **BD Veritor** System Kit includes swabs with a flocked tip for nasal specimen collection.

- 2. Insert the swab into one nostril of the patient.

- 3. Rotate the swab two complete 360-degree turns; pressing firmly against the nasal mucosa to help ensure the swab obtains both cells and mucus.

- 4. Withdraw the swab from the nasal cavity. The sample is now ready for processing using the **BD Veritor** System Kit.

**Proper Nasopharyngeal Swab Sample Collection**


- 1. The **BD Veritor** System Kit includes swabs with a flocked tip for nasopharyngeal specimen collection.

2. Insert the swab into one nostril of the patient, reaching the surface of the posterior nasopharynx.

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##### Nasal and Nasopharyngeal Swab Test Procedure

###### NOTES:


- **Reagents, specimens and devices must be at room temperature (15–30 °C) prior to testing.**
- The CLIA-waived **BD Veritor** System for Rapid Detection of Flu A+B kit is only intended for nasal and nasopharyngeal swab specimens that are collected and tested directly (i.e., dry swabs that have **NOT** been placed in transport media). The kit includes a pre-diluted process reagent in a ready to use "unitized" tube. This CLIA-waived kit IS NOT INTENDED for testing liquid samples such as wash or aspirate samples or swabs in transport media as results can be compromised by over dilution.

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

Step 5:

- Press the attached tip firmly onto the **RV Reagent D** tube containing the processed sample (threading/twisting not required).
- Vortex or mix thoroughly by swirling or flicking the bottom of the tube before adding to assay device.
- Note: Do not use tips from any other product, including other products from BD or other manufacturers.


| After step 5, choose from the model and workflow option below before continuing to step 6: |                                                                                                                                                                                  |                                                                                              |                                                                                                                                                 |                                                                                              |
|--------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|
|                                                                                            | BD Veritor Reader or Analyzer in <b>Analyze Now</b> mode                                                                                                                         | BD Veritor Plus Analyzer in <b>Walk Away</b> mode                                            | BD Veritor Plus Analyzer with either The <b>BD Veritor</b> InfoScan or InfoSync module In <b>Analyze now</b> mode---or--- <b>Walk Away</b> mode |                                                                                              |
| Instructions in section:                                                                   |  <b>A</b>  | <b>B</b>  | <b>C</b>                                                     | <b>D</b>  |

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
**Step 9A: Record the Result**

- When analysis is complete, the test result appears in the display window. Record the result and discard the test device appropriately..

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
**Step 7B: Adding the specimen**

- When the display window reads "ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY"
  - Invert the tube, holding it vertically (approximately one inch above the **BD Veritor** System Flu A+B device sample well).
  - Gently squeeze the ridged portion of the tube, allowing three (3) **drops** of the processed specimen to dispense into the sample well of a labeled **BD Veritor** System Flu A+B device.

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
- **CAUTION:** A countdown timer displays the time remaining for test insertion. **Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.**

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**Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.****Step 9B: Record the Result**

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



The display window briefly shows “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the *Analyzer Instructions for Use* for configuration steps. Ignore this message and postpone this process when an assay is awaiting analysis.

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- Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.
- When analysis is complete a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.

If a printer is not connected, record the result before removing the assay device.

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
- When the display window reads: ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY:
  - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
  - **Gently** squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device. **NOTE: Squeezing the tube close to the tip may cause leakage.**
  - **CAUTION: A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.**

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1. Collect NP swab from the patient and follow Steps 1-5 of the test procedure above as instructed for the **BD Veritor** System for Rapid Detection of Flu A+B.
2. Using the sample from Step 5, continue the test procedure using the test device for RSV.
3. Refer to the product insert for **BD Veritor™** System for Rapid Detection of RSV, (Cat. # 256038) for the test procedure and full description of the **BD Veritor** RSV test. Follow the Instructions in the insert and the Instrument on-screen prompts to complete the test procedure and obtain results. Refer to the product insert for the **BD Veritor™** System RSV Kit (Cat. # 256038) for result interpretation.

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1. Insert the swab all the way into the appropriately labeled **RV Reagent D** tube and vigorously plunge the swab up and down in the fluid for a minimum of 15 seconds.
2. Continue processing the swab according to the Test Procedure for Nasal and Nasopharyngeal swabs above beginning at Step 4 "**Remove the swab.**"

**If the kit controls do not perform as expected, do not report patient results. Contact BD Technical Services at 1-800-638-8663.**

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## Analytical Studies

### Analytical Sensitivity (Limit of Detection)

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EID<sub>50</sub>/mL = 50% Egg Infectious Dose

### Strain Reactivity with Influenza A and B Viruses

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## Technical Support

For questions, or to report a problem, please call Technical Support at 1.800.638.8663. Test system problems may also be reported to the FDA using the MedWatch reporting system (phone: 1-800-FDA-1088; fax: 1-800-FDA-1078; or <http://www.fda.gov/medwatch>).

## AVAILABILITY

### Cat. No. Description

- 256045 **BD Veritor™** System for Rapid Detection of Flu A+B, 30 tests
- 256051 **BD Veritor™** System Flu A+B Control Swab Set, 10 pairs of swabs
- 220252 COPAN Flexible Minutip Flocked Swab, 100 swabs
- 256055 **BD Veritor™** System Reader
- 256066 **BD Veritor™** Plus Analyzer
- 256067 **BD Veritor™** InfoSync Module
- 256068 **BD Veritor™** InfoScan Module
- 443907 USB Printer Cable for **BD Veritor™** Analyzer
- 256038 **BD Veritor™** System for Rapid Detection of Respiratory Syncytial Virus (RSV)

## REFERENCES

1. Simonsen L., Fukuda K, Schonberger LB, Cox NJ. Impact of influenza epidemics on hospitalizations. *J. Infect. Dis.* 2000; 181:831-7
2. Thompson WW, Shay DK, Weintraub E, et al. Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA* 2003; 289:179-86
3. Treanor, J.J., Hayden, F.G., Vrooman, P.S., et al. 2000. Efficacy and safety of the oral neuraminidase inhibitor oseltamivir in treating acute influenza: a randomized controlled trial. *JAMA.* 283:1016-1024.
4. Kaiser, L., Couch, R.B., Galasso, G.J., Glezen, W.P., Webster, R.G., Wright, P.F., and Hayden, F.G. 1999. First international symposium on influenza and other respiratory viruses: summary and overview Kapalua, Maui, Hawaii, December 4-6, 1998. *Antiviral Res.*, 42:149-176
5. Cox, N.J., and Bender, C.A. 1995. The molecular epidemiology of influenza viruses. *Virology*, 6:359-370.
6. Todd, S.J., Minnich, L., and Waner, J.L. 1995. Comparison of rapid immunofluorescence procedure with TestPack RSV and Directigen Flu A for diagnosis of respiratory syncytial virus and influenza A virus. *J. Clin. Microbiol.*33:1650-1651.
7. Harris, P.O. 1989. Clinical relevance and efficient detection of seven major respiratory viruses. *ACL.* p. 15-19.
8. McElhaney, J.E., Gravenstein, S., Krause, P., Hooton, J.W., Upshaw, C.M., and Drinka, P. 1998. Assessment of markers of the cell-mediated immune response after influenza virus infection in frail older adults. *Clin. Diag. Lab. Immunol.* 5:840-844.
9. Fan, J., Henrickson, K.J., and Savatski, L.L. 1998. Rapid simultaneous diagnosis of infections with respiratory syncytial viruses A and B, influenza viruses A and B, and human parainfluenza virus types 1, 2, and 3 by multiplex quantitative reverse transcription-polymerase chain reaction-hybridization assay (hexaplex). *Clin. Infect. Disease* 26:1397-1402.
10. Wright, K.E., Wilson, G.A.R., Novosad, D., Dimock, C., Tan, D., and Weber, J.M. 1995. Typing and subtyping of influenza viruses in clinical samples by PCR. *J. Clin. Microbiol.*33:1180-1184.
11. Kendal, A.P. 1985. Influenza Viruses. p. 341-357. *Laboratory Diagnosis of Viral Infections*, In H. Lennette, (ed.) Marcel Dekker, Inc., New York.
12. McQuillen, J., Madeley, C.R., and Kendal, A.P. 1985. Monoclonal antibodies for the rapid diagnosis of influenza A and B virus infections by immunofluorescence. *Lancet.ii.* 911-914.
13. Guenther, S.H., and Linnemann, C.C., Jr. 1988. Indirect immunofluorescence assay for rapid diagnosis of influenza virus. *Laboratory Medicine.*19:581-583
14. Minnick, L.L., and Ray, C.G. 1986. Early testing of cell cultures for detection of hemadsorbing viruses. *J. Clin. Microbiol.* 25:421-422.
15. Schmidt, N.J., Ota, M., Gallo, D., and Fox, V.L. 1982. Monoclonal antibodies for rapid, strain specific identification of influenza virus isolates. *J. Clin. Microbiol.*16:763-765.
16. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
17. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.*17:53-80.
18. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S Government Printing Office, Washington, D.C.
19. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EECP). *Office Journal L262*, 17/10/2000, p.021-0045.

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