Veritor™ System

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CLIA-waived kit configured for testing nasopharyngeal swab samples freshly collected and dispensed directly onto assay test device.

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30

Determinations
Déterminations
Bestimmungen
Determinazioni
Determinaciones

For most current information regarding this Veritor System product
refer to www.bd.com/ds/veritorsystem

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1. The **BD Veritor** System RSV Kit includes swabs with a flocked tip for nasopharyngeal specimen collection.

2. Insert the swab into one nostril of the patient, reaching the surface of the posterior nasopharynx.

3. Rotate the swab over the surface of the posterior nasopharynx.

4. Withdraw the swab from the nasal cavity. The sample is now ready for processing using the **BD Veritor** System Kit.

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



After step 5, choose from the model and workflow option below before continuing to step 6:

	BD Veritor Reader or Analyzer in Analyze Now mode	BD Veritor Plus Analyzer in Walk Away mode	BD Veritor Plus Analyzer with either InfoScan or InfoSync modules In Analyze now mode---or--- Walk Away mode	
Instructions in section:	A	B	C	D

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
**Step 8C: Using the Analyzer**

During the incubation time, turn the BD Veritor Plus Analyzer on by pressing the blue button once.

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- When assay development time is complete and the Analyzer display window reads: “INSERT TEST DEVICE OR DOUBLE CLICK FOR WALK AWAY MODE”:
 - Insert the **BD Veritor** System RSV device into the **BD Veritor** Plus Analyzer.

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The results of the **BD Veritor** RSV test obtained on the same specimens by the intended users were also compared with the results obtained by viral cell culture. Among the 523 specimens tested, there were 134 specimens that tested positive and 389 that tested negative by the viral culture method. The sensitivity of the **BD Veritor** RSV assay for specimens that tested positive by the comparator viral cell culture was 91.8% (123/134) with the 95% CI of 85.9%– 95.4%, and the specificity was 93.3% (363/389) with the 95% CI of 90.4%– 95.4%). Of the 26 **BD Veritor** RSV positive, viral cell culture negative specimens there were 23 specimens that tested positive by the PCR comparator method.

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For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolni hranici teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Катэгорыя ніжняй тэмпературы / Limite inférieure de température / Alumine temperatur / Limite inferiore di temperatura / Temperaturgränzen / Limite inferior de temperatură / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температурныя межы рукост шэры / 제한 온도 / Žemiausia laikymo temperatūra / Temperatūros žemākā robeža / Laagste temperatuurilimiet / Nedre temperaturgrænse / Dolna granica temperature / Limite minimo de temperatura / Limittä minima di temperatura / Нижний предел температуры / Spodná hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sıcaklık alt sınırı / Миніміальна температура / 温度下限

Control / Контролно / Kontrola / Kontrol / Kontrolle / Μάρτυρας / Kontroll / Contrôle / Controllio / Бақылау / 컨트롤 / Kontrolè / Kontrolle / Controle /
 Control / Контроль / kontroll / Контроль / 对照

Positive control / Положительен контрол / Pozitivní kontrola / Positiv kontrol / Positive Kontrolle / Θετικός μάρτυρας / Control positivo / Positiivne kontroll / Contrôle positif / Pozitívna kontrola / Pozitiv kontroll / Controllo positivo / Он бакылай / 양성 컨트롤 / Teigiama kontrolė / Pozityva kontrolė / Posietive controle / Kontrola dodatnia / Controllo positivo / Control pozitiv / Положительный контроль / Pozitif kontrol / Позитивний контроль / 阳性对照试剂

Negative control / Отрицателен контрол / Negativní kontrola / Negativ kontrol / Negative Kontrolle / Αρνητικός μάρτυρας / Control negativo / Negativne kontroll / Contrôle négatif / Negativna kontrola / Negativ kontroll / Controllo negativo / Негативтік бақылау / 음성 컨트롤 / Neigiama kontrolė / Negativā kontrole / Negatieve controel / Kontrola ujemna / Controllo negativo / Control negativ / Отрицательный контроль / Negatif kontrol / Негативний контроль / 阴性对照试剂

Method of sterilization: ethylene oxide / Метод на стерилизација: етиленов оксид / Zrúboš sterilizácie: etylenoxid / Steriliseringensmetode: etylenoxid / Sterilisationsmethode: Etylenoxid / Μέθοδος αποστείρωσης: αιθυλένοξείδιο / Método de esterilización: óxido de etileno / Steriliseermiddeel: etyleenoksied / Méthode de stérilisation : oxyde d'éthylène / Metoda sterilizacji: etylenoksyd / Sterilizációs módszer: etilén-oxid / Metoda sterilizacije: oksid etilena / Sterilizacija ađici – tiljenje tototy / 소독 방법: 에틸렌옥사이드 / Sterilizavimo būdas: etileno oksidas / Sterilizācijas metode: etilēnoksīds / Gesteriliseerd met behulp van ethyleenoxide / Steriliseringensmetode: etylenoksid / Metoda sterilizacji: tlenek etylu / Método de esterilización: óxido de etileno / Metóda de sterilizace: óxido de etilén / Metóda sterilizácii: etilénoksid / Metóda sterilizácie: etylen oxid / Steriliseringensmetode: etenoxid / Sterilizasyon yöntemi: etilen oksit / Метод стерилизації: етиленоксид / 灭菌方法: 环氧乙烷

Method of sterilization: irradiation / Метод стерилизации: иррадиация / Způsob sterilizace: záření / Steriliseringsmetode: bestråling / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Sterileerimismeetod: kiirgus / Méthode de stérilisation : irradiation / Metoda sterilizacji: zraczenie / Sztereázás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Στερίλιση με ακτινοβολία – κύριου τύπου / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesterileerd met behulp van bestraling / Steriliseringsmetode: bestråling / Metoda sterilizacji: napromienianie / Metodo de esterilizație: iradiație / Metoda de sterilizare: iradiere / Метод стерилизации: облучение / Metoda sterilizácie: ožarovanie / Metoda sterilizacije: ozračevanje / Steriliseringsmetod: strålning / Sterilizasyon yöntemi: ışırayan / Metoda sterilizacji: opromienienie / 灭菌方法：照射

Biological Risks / Биологични рискове / Biologická rizika / Biologisk fare / Biogefährdung / Βιολογικοί κίνδυνοι / Riesgos biológicos / Bioloogilised riskid / Risques biologiques / Biološki rizik / Biológiallag veszélyes / Rischio biologico / Биологических тәуекелдер / 생물학적 위험 / Biologinis pavojus / Biologiskie riski / Biologisch risico / Biologisk risiko / Zagrożenia biologiczne / Perigo biológico / Riscuri biologice / Биологическая опасность / Biológické riziko / Biološki rizici / Biološki risk / Biološki Riskler / Биологічна небезпека / 生物学风险

Cautiun, consult accompanying documents / Вниманіе, направте справка в придружаващите документи / Pozor! Pristupajte si priloženou dokumentaci / Forsigtig, se ledsagende dokumenter / Achtung, Begleiddokumente beachten / Προσοχή, συμβουλευθείτε τα συνοδευτικά έγγραφα / Precuación, consulta la documentación adjunta / Ettevaattus! Luegda kaasvatud dokumentatsiooni / Attention, consultez les documents joints / Uprezorenje, konsultir pratecu dokumentaciju / Figeleml! Olvassa el a mellékelt tájékoztató / Attenzione: consultare la documentazione allegata / Абайланың, тиісті құжаттармен танысыңыз / 주의, 동봉된 설명서 참조 / Dmesio, zkontrolujte pridedavanu dokumentus / Piesardziba, skatit pavaddokumentus / Voorzichtig, raadpleeg begeleide documenten / Forsiktig, se vedlagte dokumentasjon / Naleyzy zapoznać się z załączonymi dokumentami / Vymahati, posmatrajte priložene dokumente / Vážení, prosíme, konzultujte přílohu / Vytřašha, posmatrajte priloženo dokumento / Vytřašha, posmatrajte priloženo dokumento / Vytřašha, posmatrajte priloženo dokumento / Obs! Se meðfarjandue dokumentu / hildite verilein þeleplere þessari / Увага, уважте придружаващите документи / 注意 添付書類を必ず確認してください

Upper limit of temperature / Горен лџит на температурата / Horní hranice teploty / Øvre temperaturgrænse / Temperaturöbergrenzen / Авърсто
 юр Върхуогрџцџ / Limite superio de temperatura / Ulemine temperaturiurip / Limite supérieure de température / Gornja dozvoljena temperatura
 / Felső hőmérsékleti határ / Limite superiore di temperatura / Температурни рџкст етџен џорџри џџ / 상한 온도 / Aukščiausia laikymo
 temperatūra / Aukščiausia temperatūros riba / Hoogste temperatuurlimiet / Øvre temperaturgrænse / Gornja granica temperature / Limite máximo
 de temperatura / Limită maximă de temperatură / Верхний предел температуры / Horní hranice teploty / Gornja granica temperature / Øvre
 temperaturgräns / Sicaklık üst sınırı / Максимальная температура / 温度上限

Keep dry / Pазетe cyxo / Skladujte v suchém prostředí / Opbevares tørt / Trocklagern / Φυλάξτε το στεγνό / Mantener seco / Hoida kuivas /
 Conserver au sec / Držati na suhom / Száraz helyen tartandó / Tenere all'asciutto / Κυβρακ κυλινδρε ψετα / 건조 상태 유지 / Laiykite sausiai
 Uzglabāt sauc / Droog houden / Holdes tørt / Przechowywać w stanie suchym / Manter seco / Å se feri de umezzala / Не допускать пона,
 enari / Uchovávať v suchu / Držiti na suvom mestu / Förråras tørt / Kuru bir sekilde muhafaza edin / Бeрeгитe вaшe вoздух / 請保持干燥

Collection time / Време на събиране / Čas odběru / Opsamlingstidspunkt / Entnahmeuhrzeit / Ωρα συλλογής / Hora de recogida / Kogumisaie /
Heure de prélèvement / Sati prikupljanja / Mintavétel időpontja / Ora di raccolta / Жинау уакыты / 수집 시간 / Paemimo laikas / Savākšanas
laiks / Verzameltijd / Tid prøvetaking / Godzina pobrania / Hora de coleta / Ora colectării / Время сбора / Doba odberu / Vreme prikupljanja /
Uopsamlingstid / Toolama zaman / Час забору / 采集时间

Peel / Обелете / Otevfete zde / Άβη / Abziehen / Αποκολλήστε / Desprendre / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Υτρίνη / қабатын алып таста / 벗기기 / Plešti cia / Atīmēt / Schillen / Trekk av / Oderwać / Destacar / Se dezilpește / Отрпнеуть / Odrhñite / Oljušiti / Dra isår / Avırma / Baktneit / 撕下

Perforation / Перфорация / Perforace / Perforering / Διτρήση / Perforación / Perforatsioon / Perforacija / Perforálás / Perforazione / Tesik tescy /
 절취선 / Perforacija / Perforăția / Perforatie / Perforacja / Perfuração / Perforare / Перфорация / Perforácia / Perforasyon / Перфорация / 穿孔

Do not use if package damaged / Не используйте, ако опаковката е повредена / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballagen er beskadiget / Inihal beschädigter Packungnicht verwenden / Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά. / No usar si el paquete está dañado / Mitte kasutada, kui pakend on kahjustatud / Ne pas l'utiliser si l'emballage est endommagé / Ne koristiti ako je oštećeno pakiranje / Ne használja, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Erep naket bűzyltán bõnca, palidalanba / पैकिज का 손상जन क्वाउ सारु क्मज / Jei pakuoatė pažeista, nenaudoti / Nelietoti, ja iepakojums bojāts / Niet gebruiken indien de verpakking beschadigd is / Må ikke brukes hvis pakke er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / A nu se folosi dacă pachetul este deteriorat / Не использовать при повреждении упаковки / Nepoužívajte, ak je obal poškodený / Ne koristiti ako je pakovanje oštećeno / Använd ej om förpackningen är skadad / Ambalaj hasar görmüşse kullanmayın / Не використовувати за пошкодженої упаковки / 如果包裝破損, 請勿使用

Keep away from heat / Пазете от топлина / Nevystavujte přilísnému teplu / Må ikke udsættes for varme / Vor Wärme schützen / Κρατήστε το μακριά από τη θερμότητα / Mantener alejado de fuentes de calor / Hoida eemal valgusest / Protéger de la chaleur / Držati dalje od izvora topline / Övja a melegtől / Tenere lontano dal calore / Салқын жерде сақта / Ялған көкейге қам / Laikyti atokiau nuo šilumos šaltinių / Sargāt no karstuma / Beschermen tegen warmte / Må ikke utsettes for varme / Przechowywać z dala od źródeł ciepła / Manter ao abrigo do calor / A se feri de căldură / Не нагревать / Uchovávať mimo zdroja tepla / Držite dalje od toplote / Får ej utsättas för värme / Isidan uzak tutun / Беретти від дії тепла / 请远离热源 / 高热源

Cut / Срежете / Odstřihněte / Klip / Schneiden / Κόψτε / Cortar / Lõigata / Découper / Reži / Vágja ki / Tagliare / Keciňja / 잘라내기 / Kirpti / Nogriez / Knippen / Kutt / Odciać / Cortar / Decupați / Отрезать / Odstřihněte / Iseći / Klipp / Kesme / Розрізати / 剪下

Collection date / Дата на събиране / Datum odběru / Opsamlingsdato / Entnahmedatum / Ημερομηνία συλλογής / Fecha de recogida / Kogumiskuupäev / Date de prélèvement / Dani prikupljanja / Mintavétel dátuma / Data di raccolta / Жинаған тізбекүні / 수집 날짜 / Paëmio data / Savākšanas datums / Verzameldatum / Dato pobrania / Data de colecta / Data colectării / Дата сбора / Datum odboru / Datum prikupljanja / Uppsamlingsdatum / Toplama tarihi / Дата забору / 采集日期 / μL/test / μL/тест / μL/Test / μL/тэст / μL/prueba / μL/тест / μL/테스트 / мкн/тест / μL/tyrimas / μL/pārbaude / μL/teste / мкн/анализ / μL/檢測

Keep away from light / Пазете от светлина / Nevystavujte světlu / Må ikke udsættes for lys / Vor Licht schützen / Κρατήστε το μακριά από το φως / Mantener alejado de la luz / Hoida eemal valgusest / Conserver à l'abri de la lumière / Držati dalje od svetla / Fény nem érheti / Tenere al riparo dalla luce / Қараңғыланған жерде ұста / 빛을 피해야 함 / Laikyti atokiau nuo šilumos šaltinių / Sargāt no gaismas / Niet blootstellen aan zonlicht / Må ikke utsettes for lys / Przechowywać z dala od źródeł światła / Manter ao abrigo da luz / Ferir de lumină / Хранить в темноте / Uchovávať mimo dosahu svetla / Držite dalje od svetlosti / Får ej utsättas för ljus / Işıktan uzak tutun / Беретти від дії світла / 请远离光线

Hydrogen gas generated / Образован е водород газ / Možnost úniku plyného vodíku / Frembringer hydrogengas / Wasserstoffgas erzeugt / Δημιουργία αερίου υδρογόνου / Producción de gas de hidrógeno / Vesinikgaasi tekitatud / Produit de l'hydrogène gazeux / Sadržaji hydrogen vodik / Hidrogén gázt fejleszt / Produzione di gas idrogeno / Газтөрескө сытери пайда болды / 수소 가스 생성됨 / Išiskiria vandenilio dujas / Rodas ūdegradis / Waterstofgas gegenereerd / Hydrogengas generert / Powoduje powstawanie wodoru / Produção de gás de hidrogénio / Generare gaz de hidrogen / Выделение водорода / Vyrobené použitím vodíka / Oslobađa se vodonik / Genererad vätegas / Açığa çıkan hidrojen gazı / Реакция з виділенням водню / 会产生氢气

Patient ID number / ИД номер на пациента / ID pacienta / Patientens ID-nummer / Patienten-ID / Αριθμός αναγνώρισης ασθενούς / Número de ID del paciente / Patsienti ID / No d'identification du patient / Identifikacijski broj pacijenta / Beteg azonosító száma / Numero ID paciente / Пациенттин идентификациялык нөмүрү / 환자 ID 번호 / Paciento identifikavimo numeris / Pacienta ID numurs / Identificationnummer van de patiënt / Patientens ID-nummer / Numer ID pacienta / Número da ID do doente / Număr ID pacient / Идентификационный номер пациента / Identifikačné číslo pacienta / ID broj pacijenta / Patientnummer / Hasta kimlik numarası / Идентификатор пациента / 患者标识号

Fragile, Handle with Care / Чупливо, Работете с необходимо внимание. / Kéhekké. Pfi manipuláci posturupje opatrné. / Forsigtig, kan gå i stykker. / Zerbrechlich, vorsichtig handhaben. / Ευφροιστο. Χειριστείτε το με προσοχή. / Frágil. Manipular con cuidado. / Őrn, kásitsege etevaatikult. / Fragile. Manipuler avec précaution. / Lomljivo, rukujte pažljivo. / Törékeny! Óvatosan kezelendó. / Fragile, maneggiare con cura. / Сыныш, абайлап пайдаланыңыз. / 조심 깨지기 쉬운 처리. / Trapu, elkites atsargiai. / Trausis; rikotes uzmanigi / Breekbaar, voorzichtig behandelen. / Ömtålig, håndter försiktigt. / Krucha zawartość, przenosić ostrożnie. / Frágil, Manuseie con Cuidado. / Frágil, manipulați cu atenție. / Хрупкое! Обращаться с осторожностью. / Krehké, využaduje sa opatrná manipulácia. / Lomljivo - rukujte pažljivo. / Bräckligt. Hantera försiktigt. / Kolay Kırılır, Dikkatli Taşın. / Тендітна, звертатися з обережністю / 易碎, 小心轻放

R_x Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Solo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

FluMist is a registered trademark of MedImmune, LLC.

© 2017 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.

Insert ends here

DO NOT PRINT THIS PAGE

CLSI document should not be included in print file to vendor.

NOTES:

The CLSI document is to be updated at the same time as the package insert and should be released as a separate PDF document to the web.

The package insert should also be released as a separate PDF document to the web.

The package insert and CLSI document should be released as one document to SAP.

Final art provided to Buyer/print vendor should not include this page or the CLSI portion of this document.

Proofreader: please proof the CLSI document along with the package insert.

Procedure* BD Veritor™ System For Rapid Detection of Group A Strep

For use with throat swab specimens.

Prepared by	Date Adopted	Supersedes Procedure #

Review Date	Revision Date	Signature

Distributed to	# of Copies	Distributed to	# of Copies

***Any modifications to this document are the sole responsibility of the facility. This “Sample Procedure” is not intended as a substitute for your facility procedure manual, instrument manual, or reagent labeling/package insert. This “Sample Procedure” is intended as a model for use by your facility to meet the needs of your laboratory.**

For Rapid Detection of Respiratory Syncytial Virus (RSV)

For use with nasopharyngeal swab specimens.

CLIA Complexity-WAIVED

For *in vitro* diagnostic use only.

Rx Only

A Certificate of Waiver is required to perform this test in a CLIA waived setting. To obtain a Certificate of Waiver, please contact your state health department.

Additional CLIA waiver information is available at the Centers for Medicare and Medicaid website at www.cms.hhs.gov/CLIA or from your state health department.

Failure to follow the instructions or modification to the test system instructions will result in the test no longer meeting the requirements for waived category.

INTENDED USE

The **BD Veritor™** System for Rapid Detection of Respiratory Syncytial Virus (RSV) is a chromatographic immunoassay with an instrumented read for the direct and qualitative detection of RSV fusion protein from a direct nasopharyngeal swab from patients suspected of having a viral respiratory infection. This test is intended for *in vitro* diagnostic use to aid in the diagnosis of RSV infections in infants and pediatric patients under the age of 6 years. Negative results do not preclude RSV infection and should not be used as the sole basis for treatment or for other management decisions. A negative test is presumptive. It is recommended that negative test results be confirmed by viral cell culture or an alternative method, such as a FDA-cleared molecular assay. Outside the U.S., a negative test is presumptive and it is recommended that these results be confirmed by viral culture or a molecular assay cleared for diagnostic use in the country of use. FDA has not cleared this device for use outside of the U.S. Negative test results do not preclude viral infection and should not be used as the sole basis for treatment or other patient management decisions. The test is intended for professional and laboratory use. It is to be used in conjunction with the **BD Veritor** System Instrument.

SUMMARY AND EXPLANATION

Viral respiratory tract infections are responsible for widespread disease. Respiratory syncytial virus is a leading cause of lower respiratory tract infections (LRI) in young children in both the developed and developing worlds. Worldwide, it is estimated that RSV is responsible for greater than 30 million cases of LRI in children under 5 years of age each year.^{1,2}

Diagnostic methods for detection of respiratory viruses include viral cell culture, direct fluorescent antibody (DFA), rapid immunoassays, and nucleic acid amplification assays such as the polymerase chain reaction (PCR).^{3,4} Each has been demonstrated to have clinical utility for the detection of respiratory viruses including RSV. Rapid immunoassays available for specific viruses such as influenza A/B and RSV allow a quick diagnosis so that patients may be appropriately isolated and treated to prevent the nosocomial spread of infections to fellow patients with compromised cardiac, respiratory or immune functions.⁵ In addition, rapid tests assist with the selection of appropriate antiviral therapy. The most common specimen types collected for RSV testing include nasopharyngeal washes, nasopharyngeal aspirates, nasal swabs and nasopharyngeal swabs.

The **BD Veritor** System for Rapid Detection of RSV (also referred to as the **BD Veritor** System and **BD Veritor** System RSV) is a chromatographic immunoassay to detect RSV fusion protein extracted from nasopharyngeal swab specimens from symptomatic patients. All **BD Veritor** System RSV test devices are interpreted by a **BD Veritor** System Instrument, either a **BD Veritor** Reader or **BD Veritor** Plus Analyzer (the "Analyzer"). When using an Analyzer, procedures to evaluate test devices depend on the workflow configuration chosen. In Analyze Now mode, the instrument evaluates assay devices after manual timing of their development. In Walk Away mode, devices are inserted immediately after application of the specimen, and timing of assay development and analysis is automated. Additionally, connection of an Analyzer to a printer is possible if desired. Additional result documentation capabilities are possible with the integration of a **BD Veritor** InfoScan ("InfoScan") or **BD Veritor** InfoSync ("InfoSync") module. Please refer to the Analyzer *Instructions for Use* for details on how to implement these features. InfoSync is not available in all regions.

PRINCIPLES OF THE PROCEDURE

The **BD Veritor** System for Rapid Detection of RSV is a qualitative, digital immunoassay for the detection of RSV fusion protein in samples processed from nasopharyngeal specimens. When specimens are processed and added to the test device, RSV antigen binds to anti-RSV antibodies conjugated to detector particles in the RSV test strip. The antigen-conjugate complex migrates across the test strip to the reaction area and is captured by the line of RSV antibody on the membrane. A positive result for RSV is determined by the **BD Veritor** System Instrument (purchased separately) when antigen-conjugate is deposited at the Test "T" position and the Control "C" position on the **BD Veritor** System RSV assay device. The Instrument analyzes and corrects for non-specific binding and detects positives not recognized by the unaided eye to provide an objective digital result.

REAGENTS

The following components are included in the **BD Veritor** System for Rapid Detection of RSV kit:

BD Veritor System RSV Devices	30 devices	Foil pouched device containing one reactive strip. Each strip has one test line of monoclonal antibody specific to RSV viral antigen and murine monoclonal control line antibodies.
RV Reagent D	30 tubes with 400 µL reagent	Detergent with < 0.1% sodium azide (preservative).
Flexible minitip flocked swab	30 each	Swab for nasopharyngeal collection
RSV Positive Control Swab	1 each	RSV Positive Control Swab, RSV antigen (noninfectious cell lysate) with < 0.1% sodium azide (preservative)
RSV Negative Control Swab	1 each	RSV Negative Control Swab, (detergent-treated non-infected cells) with < 0.1% sodium azide (preservative)

Materials Required But Not Provided: **BD Veritor™** System Reader (Cat. No. 256055) or **BD Veritor™** Plus Analyzer (Cat. No. 256066), Timer, Tube Rack for specimen testing.

Optional Equipment: **BD Veritor™** InfoScan Module (Cat. No. 256068), **BD Veritor™** InfoSync Module (Cat. No. 256067), USB Printer cable for **BD Veritor™** Analyzer (Cat. No. 443907), Epson Printer model TM-T20 II.

Warnings and Precautions:

Warning

H302 Harmful if swallowed. **H402** Harmful to aquatic life. **H412** Harmful to aquatic life with long lasting effects.

P273 Avoid release to the environment. **P264** Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

1. For *in vitro* Diagnostic Use.
2. Test results are not meant to be visually determined. **All test results must be determined using the BD Veritor System Instrument.**
3. The RSV Positive Control Swab and the positive control line on the **BD Veritor** System for Rapid Detection of RSV device have been prepared from RSV-infected tissue culture cells which have been inactivated by detergent treatment and sonication then subsequently tested by bioassay procedures.
4. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁶⁻⁹ and institutional guidelines should be followed in handling, storing and disposing of all specimens and all items contaminated with blood and other body fluids.
5. Dispose of used **BD Veritor** System test devices as biohazardous waste in accordance with federal, state and local requirements.
6. Reagents contain sodium azide, which is harmful if inhaled, swallowed or exposed to skin. Contact with acids produces very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
7. Other than the flocked swabs that are used for specimen collection, kit components should not come in contact with the patient.
8. Do not use kit components beyond the expiration date.
9. Do not reuse the **BD Veritor** System test device.
10. Do not use the kit if the Control RSV Positive Swab and Control RSV Negative Swab do not yield appropriate results.
11. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
12. To avoid erroneous results, specimens must be processed as indicated in the assay procedure section. The addition of excess sample may give invalid test results.
13. Proper specimen collection, storage and transport are critical to the performance of this test.
14. Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.

Storage and Handling: Kits may be stored at 2–30 °C. **DO NOT FREEZE.** Reagents and devices must be at room temperature (15–30 °C) when used for testing.

SPECIMEN COLLECTION AND HANDLING

Specimen Collection and Preparation: The acceptable specimen for testing with the **BD Veritor** System for Rapid Detection of RSV is a nasopharyngeal (NP) swab. It is essential that correct specimen collection and preparation methods be followed. Specimens obtained early in the course of the illness will contain the highest viral titers.

Inadequate specimen collection or improper specimen handling and/or transport may yield a false negative result; therefore, training in specimen collection is highly recommended due to the importance of specimen quality for generating accurate test results.

Specimen Transport and Storage:

Freshly collected specimens should be processed and tested within one hour. It is essential that correct specimen collection and preparation methods be followed.

1. The **BD Veritor** System RSV Kit includes swabs with a flocked tip for nasopharyngeal specimen collection.
2. Insert the swab into one nostril of the patient, reaching the surface of the posterior nasopharynx.
3. Rotate the swab over the surface of the posterior nasopharynx.
4. Withdraw the swab from the nasal cavity. The sample is now ready for processing using the **BD Veritor** System Kit.

Dos and Don'ts of Sample Collection

- Do collect sample as soon as possible after onset of symptoms
- Do test sample immediately
- BD recommends flocked swabs which are provided in the **BD Veritor** System RSV Kit
- Do not use cotton tips and wood shafts
- Do not use calcium alginate swabs

TEST PROCEDURE FOR NASOPHARYNGEAL SWABS

NOTES:


- Reagents, specimens and devices must be at room temperature (15–30 °C) for testing.
- The CLIA-waived **BD Veritor** System for Rapid detection of RSV kit is only intended for nasopharyngeal swab specimens that are collected and tested directly (i.e dry swabs that have not been placed in transport media). The kit includes a pre-diluted processing reagent in a ready to use "unitized" tube. This CLIA-waived kit is NOT INTENDED for testing liquid samples such as wash or aspirate samples or swabs in transport media as results can be compromised by over dilution.

Prepare for testing

The following steps assume that users of a **BD Veritor Plus Analyzer** have chosen and set all configuration options, and that the Analyzer is ready to use. To choose or change these settings, see the **BD Veritor Plus Analyzer Instructions for Use**, section 4.7. A printer is not necessary to display results. However, if your facility has chosen to connect the **BD Veritor Plus Analyzer** to a printer, check that the Analyzer is plugged into a power source, paper supply is adequate and any necessary network connections are enabled before testing.

Step 1: For each patient specimen:


- Remove one **RV Reagent D** tube/tip and one **BD Veritor** System RSV device from its foil pouch immediately before testing.
- Label with patient's name or ID number
- Place the labeled **RV Reagent D** tube(s) in the designated area of the tube rack.

Prepare the Sample


Step 2:

- Remove and discard the cap from the **RV Reagent D** tube corresponding to the sample to be tested.


Step 3:

- Insert the patient sample swab all the way into the **RV Reagent D** tube and swirl it against the inside wall three (3) times.


Step 4:

- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Properly dispose of the swab.

Step 5:

- Press the attached tip firmly onto the **RV Reagent D** tube containing the processed sample (threading/twisting not required).
- Vortex or mix thoroughly by swirling or flicking the bottom of the tube.
- Do not use tips from any other product, including other products from BD or other manufacturers.


After step 5, choose from the model and workflow option below before continuing to step 6:				
	BD Veritor Reader or Analyzer in Analyze Now mode	BD Veritor Plus Analyzer in Walk Away mode	BD Veritor Plus Analyzer with either InfoScan or InfoSync modules In Analyze now mode---or--- Walk Away mode	
Instructions in section:	A	B	C	D

Step 6A: Adding the specimen

- Invert the **RV Reagent D** tube and hold the tube vertically (approximately one inch above the labeled **BD Veritor** System RSV device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled **BD Veritor** System RSV device.


NOTE: Squeezing the tube too close to the tip may cause leakage.

**Step 7A: Timing development**

- After adding the sample, allow the test to run for 10 minutes before inserting into the **BD Veritor** Instrument.
- **NOTE:** If running test under laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.

**Step 8A: Using the BD Veritor Instrument:**

- During incubation time, turn the **BD Veritor** Instrument on by pressing the power button once.
- Insert assay device when the 10 minute assay development time is complete. Follow the on-screen prompts to complete the procedure.
- The status of the assay analysis process appears in the display window.

**Step 9A: Record the Result**

- When analysis is complete, the test result appears in the display window.

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).

To use Walk Away mode - Connect the AC power adapter to the Analyzer and a power source


Step 6B: Starting Walk Away mode

- Turn the Analyzer on by pressing the blue power button once.
- When the display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE,”
 - **Double-click** the blue power button.

Step 7B: Adding the specimen:


- When the display window reads “ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY”:
 - Invert the tube, holding it vertically (approximately one inch above the **BD Veritor** System RSV device sample well).
 - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled **BD Veritor** System RSV device.

NOTE: Squeezing the tube close to the tip may cause leakage.

Step 8B: Start the development and reading sequence

- Immediately insert the test device into the slot on the right side of the Analyzer.
The test device must remain horizontal to prevent spilling the specimen out of the sample well.
- “DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- A countdown timer in the display window shows the remaining analysis time.

Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.


Step 9B: Record the Result

- When analysis is complete, the test result appears in the display window.

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (if the AC power adapter is connected).

Step 6C: Adding the specimen

- Invert the tube, holding it vertically (approximately one inch above the **BD Veritor** System RSV device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) **drops** of the processed specimen into the sample well of a labeled **BD Veritor** System RSV device. **NOTE: Squeezing the tube close to the tip may cause leakage.**

**Step 7C: Timing development**


- Allow the test to develop for **10 minutes**.
- If running the test in a laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.

**Step 8C: Using the Analyzer**

During the incubation time, turn the **BD Veritor Plus Analyzer** on by pressing the blue button once.


The display window briefly shows “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the *Analyzer Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.

- When assay development time is complete and the Analyzer display window reads: “INSERT TEST DEVICE OR DOUBLE CLICK FOR WALK AWAY MODE”:
 - Insert the **BD Veritor** System RSV device into the **BD Veritor Plus Analyzer**.

**Step 9C: Using the Bar Code scanner**

- Follow the prompts on the display window to complete any required barcode scans of:
 - OPERATOR ID
 - SPECIMEN ID and/or
 - KIT LOT NUMBER.

According to site requirements and Analyzer settings.

- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 8C. To restart this step, remove and reinsert the test device to initiate a new test sequence.
- Move the barcode slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the kit Lot Number in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

After required scans are completed, the Analyzer displays a countdown timer and test analysis begins.

- **Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**
- When analysis is complete, a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.

If a printer is not connected, record the result before removing the assay device.

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).

Step 10C: Remove the test device

- Pull the device out. The display will show INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE to indicate the Analyzer is ready to perform another test.


If an InfoSync module is installed the ENVELOPE symbol will appear to indicate that results are transmitting.

- In the event that the **BD Veritor** Plus Analyzer does not detect adequate cellular network strength while the ENVELOPE symbol is still displayed, it will queue all results to be transmitted and continuously attempt to transmit them. If it is powered off during this time, it will attempt to transmit as soon as power is restored.

To use Walk Away mode - Connect the AC power adapter to the Analyzer and a power source

Step 6D: Starting Walk Away mode


- Turn the Analyzer on by pressing the blue power button once.
The display window will briefly show “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.
- When the display window reads:
“INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”
– **Double-click the blue power button.**

Step 7D: Using the Bar Code scanner

- Follow the prompts on the display window to complete any required barcode scans of:
 - OPERATOR ID
 - SPECIMEN ID and/or
 - KIT LOT NUMBER


According to site requirements and Analyzer settings.

- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 6D. To restart this step, double-click the power button.
- Move the barcode slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the kit Lot Number in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

Step 8D: Add the specimen to the test device:

- When the display window reads: “ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY”:
 - Invert the tube, holding it vertically (approximately one inch above the **BD Veritor** System RSV device sample well).
 - Gently squeeze the ridged portion of the tube away from the tip, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled **BD Veritor** System RSV device. **NOTE: Squeezing the tube close to the tip may cause leakage.**

Step 9D: Start the development and reading sequence


- Immediately insert the test device into the slot on the right side of the Analyzer.

The test device must remain horizontal to prevent spilling the specimen out of the sample well.

- “DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- A countdown timer in the display window shows the remaining analysis time.

Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.

- When analysis is complete, a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed. **If a printer is not connected, record the result before removing the assay device.**

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).

Step 10D: Remove the test device

- Pull the device out. The display will show “INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE” to indicate the Analyzer is ready to perform another test. Note that the Analyzer returns to Analyze Now mode at the conclusion of each read sequence.

If an InfoSync module is installed the ENVELOPE symbol will appear to indicate that results are transmitting.

- In the event that the **BD Veritor** Plus Analyzer does not detect adequate cellular network strength while the ENVELOPE symbol is still displayed, it will queue all results to be transmitted and continuously attempt to transmit them. If it is powered off during this time, it will attempt to transmit as soon as power is restored.

OPTIONAL TEST PROCEDURE: Testing for RSV and INFLUENZA A+B using a single NP swab. For patients under the age of 6 years.

Note: The BD Veritor™ System for Rapid Detection of Flu A+B (Catalog # 256045) is required for this procedure in addition to the BD Veritor™ System for Rapid Detection of RSV (Catalog # 256038).

This procedure allows for use of the remaining processed sample from Step 5 above to test additionally for Influenza A+B. When using this optional test procedure, the sample may be used up to 15 minutes after initial processing.

1. Collect NP swab from the patient and follow Steps 1–5 of the test procedure above as instructed for RSV.
2. Using the sample from Step 5, continue the test procedure using the test device for Flu A+B.
3. Refer to the product insert for **BD Veritor™ System for Rapid Detection of Flu A+B**, (Cat. # 256045) for the test procedure and full description of the **BD Veritor** Flu A+B test. Follow the Instrument on-screen prompts to complete the procedure and obtain test results. Refer to the product insert for the **BD Veritor™ System Flu A+B CLIA waived Kit** (Cat. # 256045) for result interpretation.

INTERPRETATION OF RESULTS:

The **BD Veritor** System Instrument (purchased separately) must be used for interpretation of all test results. Operators should not attempt to interpret assay results directly from the test strip contained within the **BD Veritor** System RSV assay device.

Display	Interpretation
RSV: +	Positive Test for RSV (RSV antigen present)
RSV: -	Negative Test for RSV (no RSV antigen detected)
CONTROL INVALID	Control line error. Repeat the test.

Invalid Test – If the test is invalid, the **BD Veritor** System Instrument will display “CONTROL INVALID” and the test or control must then be repeated. If the “CONTROL INVALID” reading recurs, contact BD Technical Support.

REPORTING OF RESULTS

- Positive Test** Positive for the presence of RSV antigen. A positive result may occur in the absence of viable virus.
- Negative Test** Negative for the presence of RSV antigen. Infection due to RSV cannot be ruled out because the antigen present in the sample may be below the detection limit of the test. In the U.S., a negative test is presumptive and it is recommended that these results be confirmed by viral cell culture or an FDA-cleared RSV molecular assay.
- Control Invalid** Do not report results. Repeat the test.

QUALITY CONTROL:

To utilize an Analyzer’s QC documentation capability, the Analyzer must be equipped with either an InfoScan or InfoSync module and specimen barcode scanning must be enabled. Please refer to the Analyzer *Instructions for Use*, section 4, to choose or change this configuration, and for specific QC procedure steps.

Each **BD Veritor** RSV device contains both positive and negative internal/procedural controls:

1. The internal positive control validates the immunological integrity of the device, proper reagent function, and assures correct test procedure.
2. The membrane area surrounding test lines functions as a background check on the assay device.

The **BD Veritor** System Instrument evaluates the positive and negative internal/procedural controls after insertion of each **BD Veritor** System test device. The **BD Veritor** System Instrument prompts the operator if a quality issue occurs during assay analysis. Failure of the internal/procedural controls will generate an invalid test result. **NOTE: The internal controls do not assess proper sample collection technique**

External Positive and Negative Controls:

RSV + and RSV - control swabs are supplied with each kit. These controls provide additional quality control material to assess that the test reagents and the **BD Veritor** System Instrument perform as expected. Prepare kit control swabs and test using the same procedure (either Analyze Now or Walk Away mode) as used for patient specimen swabs. When using the barcode scanning feature to document QC procedures, scan the barcode on the control swab packaging when prompted for a Specimen ID.

Your laboratory’s standard Quality Control procedures and applicable local, state and/or federal regulations or accreditation requirements dictate the performance of external quality control procedures.

ATTENTION: To utilize the Analyzer’s QC documentation capability, specimen barcode scanning must be enabled on an Analyzer equipped with either an InfoScan or InfoSync module. Please refer to the Analyzer *Instructions for Use*, section 4, to choose or change this configuration.

BD recommends controls be run once for:

- each new kit lot,
- each new operator,
- each new shipment of test kits,
- as required by internal quality control procedures and in accordance with local, state and federal regulations or accreditation requirements.

Test Procedure For Kit Swab Controls:

1. Remove and discard the cap from the **RV Reagent D** tube corresponding to the sample to be tested.
 2. Insert the control swab into the tube and vigorously plunge the swab up and down in the fluid for a minimum of 15 seconds.
 3. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
 4. Continue processing the swab according to the Test Procedure for Nasopharyngeal Swabs above, beginning at Step 5.
- If the kit controls do not perform as expected, do not test patient specimens. Contact BD Technical Services at 1.800.638.8663.

LIMITATIONS OF THE PROCEDURE

- Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- The contents of this kit are to be used for the qualitative detection of RSV antigens from nasopharyngeal swabs.
- The **BD Veritor** System for Rapid Detection of RSV is capable of detecting both viable and non-viable RSV particles. The **BD Veritor** System for Rapid Detection of RSV performance depends on antigen load and may not correlate with other diagnostic methods performed on the same specimen.
- Results from the **BD Veritor** System for Rapid Detection of RSV test should be correlated with the clinical history, epidemiological data and other data available to the clinician evaluating the patient.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly; therefore, a negative test result does not eliminate the possibility of RSV infection, and should be confirmed by viral cell culture or in the U.S., an FDA-cleared RSV molecular assay.
- Positive test results do not rule out co-infections with other pathogens.
- Negative test results are not intended to rule in other non-RSV viral or bacterial infections.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false positive results during periods of little/no RSV activity when disease prevalence is low. False negative test results are more likely during peak RSV activity when prevalence of disease is high.
- This device has been evaluated for use with human specimen material only.
- Monoclonal antibodies may fail to detect, or detect with less sensitivity, RSV viruses that have undergone minor amino acid changes in the target epitope region.
- The performance of this test has not been evaluated for use in patients without signs and symptoms of respiratory infection.
- The validity of the **BD Veritor** System for Rapid Detection of RSV test has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.
- Therapeutic anti-RSV monoclonal antibodies may interfere with the **BD Veritor** System for Rapid Detection of RSV.
- Performance characteristics have not been established for use with patients older than 5 years of age or for immunocompromised patients.

EXPECTED VALUES

The rate of positivity observed in RSV testing will vary depending on the method of specimen collection, handling/transport system employed, detection method utilized, time of year, age of the patient, geographic location and most importantly, local disease prevalence. In the 2012/2013 clinical trial the overall prevalence of RSV as determined by viral cell culture positive specimens was 25.6% (range of 7.7 to 65.2%). Overall prevalence of RSV as determined by PCR positive specimens tested was 34.2% (range of 15.4% to 69.6%).

PERFORMANCE CHARACTERISTICS

Clinical Performance

Performance characteristics for the **BD Veritor** System for Rapid Detection of RSV test were established in a prospective multi-center clinical study conducted at eight Point of Care (POC) U.S. testing sites during the 2012–2013 respiratory season. The performance of the **BD Veritor** RSV test was compared to a commercially available PCR method as well as viral culture. A total of 540 specimens were enrolled in the clinical trial. A total of 523 specimens had evaluable results for all three testing methods: PCR, viral cell culture and **BD Veritor** RSV. The study population was 42.6 % females and 57.4% males. The table below shows age distribution of the study population.

Demographics Summary – Age Group		
Age Group	Number	Percentage
<2	305	58.3
2 – 5	218	41.7
Total	523	100

The table below summarizes the performance obtained with clinical specimens using the **BD Veritor** System RSV test in comparison to a commercially available PCR. The overall positive percent agreement (PPA) and negative percent agreement (NPA) of the **BD Veritor** System RSV with a PCR comparator, based on these 523 specimens, are 81.6% (146/179) and 99.1% (341/344), respectively.

BD Veritor RSV Compared to PCR			
BD Veritor RSV	PCR		
	P	N	Total
P	146	3	149
N	33	341	374
Total	179	344	523
Reference Method: PCR			
PPA: 81.6% (95% CI: 75.2%, 86.6%)			
NPA: 99.1% (95% CI: 97.5%, 99.7%)			

BD Veritor System RSV test performance compared to viral cell culture was also evaluated in this study. For the same 523 specimens 91.8% (123/134) were positive by both **BD Veritor** RSV and cultures, 93.3% (363/389) were negative by both **BD Veritor** RSV and culture. There were 26 **BD Veritor** RSV positive, viral cell culture negative specimens of which 23 were demonstrated to be RSV positive by an FDA cleared molecular assay.

Invalid rates for the **BD Veritor** System for Rapid Detection of RSV were calculated by dividing the number of invalids by the total number of evaluable specimens tested by the **BD Veritor** System. The overall invalid rate for the **BD Veritor** System for RSV based on the 523 specimens was determined to be 0.2% (1/523, 95% CI: 0.0%, 1.1%).

Reproducibility

The reproducibility of the **BD Veritor** System for Rapid Detection of RSV test was evaluated at two POC sites (P-1, P-2) and one clinical laboratory site (S-1). The reproducibility panel was composed of 12 simulated RSV samples. These included moderate positive samples, low positive samples (near the assay limit of detection), high negative samples (i.e., containing very low concentrations of virus) and negative samples. The panel was tested by two operators at each site for five consecutive days. The results are summarized below.

BD Veritor RSV Reproducibility (% RSV Positive Results)								
Sample	P-1		P-2		S-1		Total	
	% Pos	95% CI	% Pos	95% CI	% Pos	95% CI	% Pos	95% CI
High negative RSV	6.7% (2/30)	1.8%, 21.3%	6.7% (2/30)	1.8%, 21.3%	13.3% (4/30)	5.3%, 29.7%	8.9% (8/90)	4.6%, 16.6%
Low positive RSV	90.0% (27/30)	74.4%, 96.5%	76.7% (23/30)	59.1%, 88.2%	80.0% (24/30)	62.7%, 90.5%	82.2% (74/90)	73.1%, 88.8%
Moderate positive RSV	100% (30/30)	88.6%, 100%	100% (30/30)	88.6%, 100%	100% (30/30)	88.6%, 100%	100% (90/90)	95.9%, 100%
Negative	0% (0/30)	0%, 11.3%	0% (0/30)	0%, 11.3%	0% (0/30)	0%, 11.3%	0% (0/90)	0%, 4.1%

Analytical Studies

Analytical Sensitivity (Limit of Detection)

The limit of detection (LOD) for the **BD Veritor** System for Rapid Detection of RSV test was established for the following RSV strains. The LOD for each strain represents the lowest concentration producing a positivity rate of $\geq 95\%$ based on testing 60 to 80 replicates.

RSV Viral Strain	Calculated LOD (TCID ₅₀ /mL)	No. Positive / Total	% Positive
VR-26 (Long subgroup A)	1.43×10^5	57/60	95.0
VR-955 (9320 subgroup B)	3.98×10^4	57/60	95.0
VR-1540 (A-2)	1.94×10^3	59/60	98.3
VR-1580 (Washington subgroup B)	1.08×10^4	58/60	96.7
VR-1400 (Wild Type subgroup B)	2.96×10^3	76/80	95.0

TCID₅₀/mL = 50% Tissue Culture Infectious Dose

Analytical Specificity (Cross Reactivity)

The **BD Veritor** System for Rapid Detection of RSV test was evaluated with bacteria and yeast at a target concentration of approximately 5×10^6 CFU/mL (CFU – Colony Forming Units). The viruses were evaluated at concentrations of 10^4 TCID₅₀/mL or greater. Of the microorganisms tested, none showed cross-reactivity in the RSV test.

<i>Bacteroides fragilis</i>	<i>Neisseria</i> sp. (<i>Neisseria perflava</i>)	Adenovirus, Type 1
<i>Bordetella pertussis</i>	<i>Neisseria subflava</i>	Adenovirus, Type 7
<i>Candida albicans</i>	<i>Peptostreptococcus anaerobius</i>	Cytomegalovirus
<i>Chlamydia pneumoniae</i>	<i>Porphyromonas asaccharolyticus</i>	Enterovirus
<i>Corynebacterium diphtherium</i>	<i>Prevotella oralis</i>	HSV Type 1
<i>Escherichia coli</i>	<i>Propionibacterium acnes</i>	Human Coronavirus OC43
<i>Fusobacterium nucleatum</i>	<i>Proteus mirabilis</i>	Human metapneumovirus (hMPV-27 A2)
<i>Haemophilus influenzae</i>	<i>Pseudomonas aeruginosa</i>	Human Parainfluenza
<i>Haemophilus parainfluenzae</i>	<i>Serratia marcescens</i>	Influenza A/California/7/2009 H1N1
<i>Kingella kingae</i>	<i>Staphylococcus aureus</i>	Influenza A/Brisbane/10/2007 H3N2
<i>Klebsiella pneumoniae</i>	<i>Staphylococcus epidermidis</i>	Influenza A/Victoria/3/75 H3N2
<i>Lactobacillus</i> sp.	<i>Streptococcus mutans</i>	Influenza B/Brisbane/60/2008
<i>Legionella</i> sp.	<i>Streptococcus pneumoniae</i>	Influenza B/Florida/4/2006
<i>Moraxella catarrhalis</i>	<i>Streptococcus pyogenes</i>	Influenza B/Lee/40
<i>Mycobacterium tuberculosis</i>	<i>Streptococcus</i> sp. Group C	Measles virus
<i>Mycoplasma pneumoniae</i>	<i>Streptococcus</i> sp. Group G	Mumps virus
<i>Neisseria gonorrhoeae</i>	<i>Streptococcus salivarius</i>	Rhinovirus
<i>Neisseria meningitidis</i>	<i>Veillonella parvula</i>	
<i>Neisseria mucosa</i>		

Interfering Substances

Various substances were evaluated with the **BD Veritor™** System for Rapid Detection of RSV test. These substances included whole blood (2%) and various medications. No interference was noted with this assay for any of the substances at the concentrations tested.

Substance	Concentration
Ayr Saline Nasal Gel	10 mg/mL
4-Acetamidophenol	10 mg/mL
Acetylsalicylic acid	20 mg/mL
Albuterol	0.083 mg/mL
Amantadine Hydrochloride	500 ng/mL
Beclomethasone	500 ng/mL
Budesonide	500 ng/mL
Chlorpheniramine maleate	5 mg/mL
Dexamethasone	10 mg/mL
Dextromethorphan	10 mg/mL
Diphenhydramine HCl	5 mg/mL
Fexofenadine	500 ng/mL
FluMist™	1%
Flunisolide	500 ng/mL
Fluticasone	500 ng/mL
Guaiacol Glyceryl Ether	20 mg/mL
Homeopathic Allergy Medicine	10 mg/mL
Ibuprofen	10 mg/mL
Loratidine	100 ng/mL

Substance	Concentration
Menthol Throat Lozenges	10 mg/mL
Mometasone	500 ng/mL
Mupirocin	500 ng/mL
Oseteltamivir	500 ng/mL
Oxymetazoline	0.05 mg/mL
Phenylephrine	1 mg/mL
Pseudoephedrine HCl	20 mg/mL
Purified Mucin Protein	1 mg/mL
Ribavirin	500 ng/mL
Rimantadine	500 ng/mL
Synagis	4 µg/mL
Tobramycin	500 ng/mL
Triamcinolone	500 ng/mL
Zanamivir	1 mg/mL
Four OTC nasal sprays	10%
Four OTC throat drops	12.5%
Two OTC mouthwashes	5%
Whole Blood	2%

CLIA WAIVER STUDY

The accuracy of the **BD Veritor** System for Rapid Detection of RSV test was evaluated at eight non-laboratory intended use Point of Care (POC) sites. A total of 22 operators representative of CLIA waived site personnel (intended users) participated in the study. No training on the use of the test was provided. The **BD Veritor** RSV test results obtained by the intended users were compared with results obtained by a commercially available PCR method. There were 523 prospectively collected specimens evaluated in this study. Among those there were 179 specimens that tested positive and 344 that tested negative by the PCR method. The positive percent agreement (PPA) of the **BD Veritor** RSV assay for specimens that tested positive by the comparator PCR method was 81.6% (146/179) with the 95% confidence interval (CI) of 75.2%–86.6%, and the negative percent agreement (NPA) was 99.1% (341/344) with the 95% CI of 97.5%– 99.7%.

The results of the **BD Veritor** RSV test obtained on the same specimens by the intended users were also compared with the results obtained by viral cell culture. Among the 523 specimens tested, there were 134 specimens that tested positive and 389 that tested negative by the viral culture method. The sensitivity of the **BD Veritor** RSV assay for specimens that tested positive by the comparator viral cell culture was 91.8% (123/134) with the 95% CI of 85.9%– 95.4%, and the specificity was 93.3% (363/389) with the 95% CI of 90.4%– 94.4%. Of the 26 **BD Veritor** RSV positive, viral cell culture negative specimens there were 23 specimens that tested positive by the PCR comparator method.

The summary of the results is presented below.

BD Veritor RSV Compared to PCR				BD Veritor RSV Compared to Viral Culture		
BD Veritor RSV	PCR			Viral Culture		
	P	N	Total	P	N	Total
P	146	3	149	123	26*	149
N	33	341	374	11	363	374
Total	179	344	523	134	389	523
PPA: 81.6% (95% CI: 75.2%, 86.6%) NPA: 99.1% (95% CI: 97.5%, 99.7%)				Sensitivity: 91.8% (95% CI: 85.9%, 95.4%) Specificity: 93.3% (95% CI: 90.4%, 94.4%)		
				*Of the 26 BD Veritor RSV positive, viral cell culture negative specimens, there were 23 specimens that tested positive by the PCR comparator method.		

The overall invalid rate for the **BD Veritor** System for RSV based on the 523 specimens tested in the study was 0.2% (1/523) with the 95% CI of 0.0%–1.1%.

Another study was designed to assess the capability of untrained users to test weakly reactive samples and deliver results with accuracy. This study was conducted at three CLIA waived intended use sites using simulated swab samples. The samples were spiked with RSV virus targeting three concentrations (high negative ~5% positivity, low positive ~95% positivity and moderate positive ~100% positivity). The swabs were provided to the operators in panels that were masked and randomized prior to shipment to the sites. The samples were also tested by trained laboratorians at one clinical site. Each site had two operators who each tested the panel on each of ten days. Each swab was processed and tested in a single device according to the test procedure.

The table below shows the rate of RSV detection for weakly reactive samples when the test was used by untrained intended users. The rate of RSV detection for weakly reactive samples when the test was used by trained laboratorians at one clinical laboratory site is also provided.

Sample	Untrained Users						Trained Users	
	POC 1		POC 2		POC 3		Clinical Site	
	% Pos	95% CI	% Pos	95% CI	% Pos	95% CI	% Pos	95% CI
High negative RSV	5.0% (1/20)	0.9%, 23.6%	0% (0/20)	0%, 16.1%	5.0% (1/20)	0.9%, 23.6%	15.0% (3/20)	5.2%, 36.0%
Low positive RSV	75.0% (15/20)	53.1%, 88.8%	75.0% (15/20)	53.1%, 88.8%	70.0% (14/20)	48.1%, 85.5%	75.0% (15/20)	53.1%, 88.8%
Moderate positive RSV	100% (20/20)	83.9%, 100%	100% (20/20)	83.9%, 100%	100% (20/20)	83.9%, 100%	100% (20/20)	83.9%, 100%

Using risk analysis as a guide, analytical flex studies were conducted. The studies demonstrated that the test is insensitive to stresses due to environmental conditions or from potential user errors.

Technical Support

For questions, or to report a problem please call **Technical Support at 1.800.638.8663**. Test system problems may also be reported to the FDA using the MedWatch reporting system (phone: 1.800.FDA.1088; fax: 1.800.FDA.1078; or <http://www.fda.gov/medwatch>).

AVAILABILITY

Cat. No. Description

256038	BD Veritor™ System for the Rapid Detection of RSV, 30 test kit
256045	BD Veritor™ System for Rapid Detection of Flu A+B, 30 test kit
256061	BD Veritor™ System RSV Control Swab Set, 10 pairs of swabs
220252	COPAN Flexible Minutip Flocked Swab, 100 swabs
256055	BD Veritor™ System Reader
256066	BD Veritor™ Plus Analyzer
256068	BD Veritor™ InfoScan Module
256067	BD Veritor™ InfoSync Module
443907	USB Printer Cable for BD Veritor™ Analyzer

REFERENCES

- Hall CB, Weinberg GA, Iwane MK, et al., The Burden of Respiratory Syncytial Virus Infection in Young Children. *N Engl J Med* 2009;360:588–98.
- Nair H, Nokes DJ, Gessner BD, et al., Global burden of acute lower respiratory infections due to respiratory syncytial virus in young children: a systematic review and meta-analysis. *Lancet* 2010;375:1545–55.
- Falsey AR and EE Walsh, Respiratory syncytial virus infection in adults. *Clin Microbiol Rev* 2000;13:371–84.
- Murata Y and AR Falsey, RSV Infection in Elderly Adults, In: Patricia Cane, Editor(s), *Perspectives in Medical Virology*, 2006, Elsevier, Volume 14, Pages 163–82.
- Crowcroft NS, Cutts F and Zambon MC, Respiratory syncytial virus: an underestimated cause of respiratory infection, with prospects for a vaccine. *Commun Dis Public Health*. 1999;2:234–41.
- Thompson WW, Shay DK, Weintraub E, et al., Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA*. 2003;289:179–86.
- Henrickson KJ and CB Hall, Diagnostic assays for respiratory syncytial virus disease. *Pediatr Infect Dis J* 2007;26 (Suppl):S36–40.
- Popow-Kraupp, T and JH Aberle, Diagnosis of Respiratory Syncytial Virus Infection. *Open Microbiol J* 2011;5:128–34. Published online 2011 December 30.
- Barenfanger J, Drake C, Leon N, et al., Clinical and financial benefits of rapid detection of respiratory viruses: an outcomes study. *J Clin Microbiol* 2000;38:2824–8.
- Clinical and Laboratory Standards Institute. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI 2005. Wayne, PA.

11. Garner JS, Hospital Infection Control Practices Advisory Committee, U. S. Department of Health and Human Services, Center for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect Control Hospit Epidemiol* 1996;17:53–80.
12. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
13. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021–0045.

Technical Information: In the United States, contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.