



URINE DRUG SCREENING DEVICE

# iCUP<sup>®</sup> AND iCUP<sup>®</sup> A.D.

Our most widely used  
rapid drug test cup

- Integrated test cup screens for up to 14 prescription and illicit drugs (3-13 panel options available)
- Self-contained collection cup and test help minimize exposure to urine
- Fast results in 5 minutes
- Try the iCup<sup>®</sup> A.D. Drug Screen for expanded specimen validity testing: CRE, GLUT, NIT, OX, pH, SG<sup>1</sup>



*For professional in vitro diagnostic use.*

# iCup® and iCup® A.D. Urine Drug Screen


## Drug screen options and cutoff levels

iCUP® AND iCUP® A.D.	CALIBRATOR	CUTOFF LEVEL (ng/mL)
Amphetamine (AMP 1000)	d-Amphetamine	1,000
Amphetamine (300)	d-Amphetamine	300
Barbiturates (BAR)	Secobarbital	300
Benzodiazepines (BZO)	Oxazepam	300
Buprenorphine (BUP)	Buprenorphine	10
Cocaine (COC 300)	Benzoyllecgonine	300
Cocaine (COC 150)	Benzoyllecgonine	150
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9-COOH	50
Methadone (MTD)	Methadone	300
Methamphetamine (mAMP 1000)	d-Methamphetamine	1,000
Methamphetamine (mAMP 500)	d-Methamphetamine	500
Methylenedioxymethamphetamine (MDMA - Ecstasy)	d,l-Methylenedioxymethamphetamine	500
Opiate (MOP 300)	Morphine	300
Opiate (OPI 2000)	Morphine	2,000
Oxycodone (OXY)	Oxycodone	100
Phencyclidine (PCP)	Phencyclidine	25
Propoxyphene (PPX)	Propoxyphene	300
Tricyclic Antidepressants (TCA)	Nortriptyline	1,000

1. **SPECIMEN VALIDITY TESTS:** Creatinine (CRE), Glutaraldehyde (GLUT), Nitrite (NIT), Oxidants/PCC (OX), Potential Hydrogren (pH), Specific Gravity (SG)

## SIMPLE COLLECTION AND TESTING PROCEDURE

**iCup® Device Components**



**1** Collect specimen and secure the lid.



**2** Peel label to reveal drug test and specimen validity (if applicable).



**3** Interpret drug test results at 5 minutes.



PRODUCT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings. For healthcare professionals including professionals at point of care sites. Immunoassay for in vitro diagnostic use. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography-mass spectrometry (GC-MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

**TO LEARN MORE, CALL 800-340-4029  
OR VISIT ALERETOXICOLOGY.COM**

© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. MKT50266 REV2 2/19

