Fentanyl Urine Control Set

Instructions for Use for 3009-5 Fentanyl Urine Control Set

Contents:

English	2
Français	4
Italiano	6
Deutsch	8
Español	10

English

FENTANYL URINE CONTROL SET

IVD

REF

3009-5

INTENDED USE

The Immunalysis Fentanyl Urine Control Set is used as control material in the Immunalysis SEFRIA™ Fentanyl Urine Enzyme Immunoassay.

Refer to the SEFRIA™ Fentanyl Urine Immunoassay Instructions for Use for additional information.

MATERIALS PROVIDED

Contents

Description: SEFRIA™ Fentanyl Urine Control Set

- Low Control: One 5 mL bottle of synthetic urine containing 0.5 ng/mL of fentanyl with stabilizers and sodium azide as a preservative.
- High Control: One 5 mL bottle of synthetic urine containing 1.5 ng/mL of fentanyl with stabilizers and sodium azide as a preservative.
- One Instructions for Use.

The insert range for the control is not lot specific and represents the total range of values which may be generated throughout the life of the product. It is recommended that each laboratory establish its own acceptable range which should fall within the instructions for use range.

ADDITIONAL MATERIALS REQUIRED (SOLD SEPARATELY)

REF	Description
NOTE: Some kit sizes are not available in all countries. Please contact your local distributor.	
618UR-0025 618UR-0100 618UR-0500 618UR-0060W	Immunalysis SEFRIA™ Fentanyl Urine Reagent Kit
NEG-600-5	Immunalysis Negative Urine Calibrator, 1 x 5 mL
10019-5	Immunalysis Fentanyl Urine Calibrator 1 containing 1 ng/mL of fentanyl, 1 x 5 mL
10020-5	Immunalysis Fentanyl Urine Calibrator 2 containing 2 ng/mL of fentanyl, 1 x 5 mL
10021-5	Immunalysis Fentanyl Urine Calibrator 3 containing 4 ng/mL of fentanyl, 1 x 5 mL
Clinical Chemistry Analyzer system suitable for running the SEFRIA™ Fentanyl Urine reagents. Refer to SEFRIA™ Fentanyl Urine reagent kit Instructions for Use for further details.	

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Exercise the normal precautions required for handling all laboratory reagents. Please refer to the Instructions for Use for the associated assay.
- It is the responsibility of the operator to ensure correct result interpretations for controls, calibrators and specimen.
- Dispose of contents and container, in accordance with local, regional, national and/or international regulation.

- Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.
- Safety Data Sheets are available at www.globalpointofcare.eifu.abbott or contact your local representative.
- No substances contained within the Fentanyl Urine Control Set are classified as hazardous in accordance with Regulation (EC) No. 1272/2008 [CLP] and US Hazcom 2012. No hazard pictograms, hazard statements or precautionary statements required. Refer to Safety Data Sheet for further details.

STORAGE, STABILITY AND HANDLING

The controls are provided ready to use and may be used directly from the refrigerator.

The control bottles should remain closed when not in use with screw caps tightly closed.

Store refrigerated at 2-8°C. When stored and handled as directed, the unopened controls are stable until the expiration date listed on the label. The opened bottles are stable for 60 days or the printed expiration date, whichever comes first. Do not use the controls beyond their expiration dates.

DO NOT FREEZE

Ensure all controls are stored in accordance with the specified conditions. Storage conditions are critical to the performance of the assay.

Do not interchange or mix controls from different lots.

Do not interchange control caps.

PREPARATIONS FOR USE

The controls are ready to use and do not require reconstitution. Mix the contents of the bottle before each use by gently inverting the bottle 3-4 times. Record the date on which the control was opened.

Use the controls as specified in the SEFRIA™ Fentanyl Urine Instructions for Use or Application Sheet.

TRACEABILITY

The SEFRIA™ Fentanyl Urine Control Set is standardized from a commercially available fentanyl standard.

INSTRUMENT PROCEDURE

Before performing the assay, refer to the assay-specific reagent Instructions for Use for additional instructions on performing quality control for qualitative analysis.

Verify that the correct control values have been entered into the control file on the analyzer.

INDICATIONS OF INSTABILITY OR DETERIORATION


Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, or if control is not within the established ranges.

English

TECHNICAL ASSISTANCE

If technical assistance is required when using this product, please contact Immunalysis Corporation at the contact details below.

Region	Group Emails
Europe & Middle East	Phone: +44 161 483 9032 Option 3 Email: EME.TechSupport@abbott.com
Asia Pacific	Phone: +61 7 3363 7100 Email: AP.TechSupport@abbott.com
Africa	Phone: +27 10 500 9700 Option 3 Email: arcis.techsupport@abbott.com
Russia, and the Commonwealth of Independent States (RCIS)	Phone: +44 161 483 9032 Option 3 Email: arcis.techsupport@abbott.com
Latin America	Phone: +57 1 482 4033 Email: LA.TechSupport@Abbott.com
Canada	Phone: 1-800-818-8335 Email: Canproductsupport@abbott.com
US	Phone: +1 877 441 7440, Option 2 Email: USproductsupport@abbott.com

CEpartner4U BV
Esdoornlaan 13, 3951DB Maarn,
The Netherlands
P: +31.343.444.875
CEPartner4U.eu

IMMUNALYSIS Corporation
829 Towne Center Drive, Pomona, CA 91767
+1 (909) 482-0840
www.immunalysis.com

Additional copies of the Instructions for Use are available from the Abbott Rapid Diagnostics website: www.globalpointofcare.eifu.abbott

© 2022 Immunalysis Corporation. All rights reserved. All trademarks referenced are trademarks of Immunalysis Corporation or their respective owners. SEFRIA™ is a trademark of Advanced Immunoassay, LLC.
KI-3009-IMM vAG 09/22

GLOSSARY OF SYMBOLS

Symbol	Description	ISO 15223-1 ⁽¹⁾ Reference #
	Authorized representative in the European Community	5.1.2
	Batch code	5.1.5
	Catalogue number	5.1.6
	CE Marking of Conformity	Not Applicable
	Consult Instructions for Use	5.4.3
	In vitro diagnostic medical device	5.5.1
	Manufacturer	5.1.1
	QR Code for Unique Device Identifier (UDI)	Not Applicable
	Temperature limit	5.3.7
	Use by date	5.1.4

⁽¹⁾ ISO 15223-1:2016 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements

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