Liquid QC® Immunoassay Control

REF

94102

LOT

1706091B

2021-02

INTENDED USE

FOR IN VITRO DIAGNOSTIC USE

CLINIQA Liquid QC Immunoassay Control is intended for use as an assayed quality control material for the constituents listed in this package insert. CLINIQA Liquid QC Immunoassay Control is not intended for use as a standard.

SUMMARY

The use of quality control materials to objectively monitor the accuracy and precision of procedures is well established in clinical laboratories. CLINIQA Liquid QC Immunoassay Control is provided at three levels to assist in the monitoring of analytical systems within the clinical range.

PRODUCT DESCRIPTION

CLINIQA Liquid QC Immunoassay Control is human serum-based containing constituents of human and animal origin as well as purified chemicals and therapeutic drugs. Preservatives and stabilizers have been added to maintain product integrity. CLINIQA Liquid QC Immunoassay Control is a ready-to-use liquid control requiring no reconstitution.

STORAGE AND STABILITY

CLINIQA Liquid QC Immunoassay Control is stable until the expiration date on the vial label when stored unopened at $\leq -20^{\circ}\text{C}$. Once opened, vials of control are stable for 30 days when stored tightly capped at $2-8^{\circ}\text{C}$.

PROCEDURE

CLINIQA Liquid QC Immunoassay Control should be treated the same as patient specimens and run in accordance with the instructions accompanying the test system being used. Gently mix the contents of each vial before sampling to ensure homogeneity. Replace cap immediately and store at $2-8^{\circ}\text{C}$.

QC materials should be used in accordance with local, state, and/or federal regulations or accreditation requirements.

LIMITATIONS OF PROCEDURE

CLINIQA Liquid QC Immunoassay Control should not be used past the expiration date on the vial label.

CLINIQA Liquid QC Immunoassay Control is a stabilized liquid product.

To obtain consistent assay values, CLINIQA Liquid QC Immunoassay Control requires storage and handling as detailed in STORAGE and STABILITY.

Accurate and reproducible results are dependent upon properly functioning instruments and reagents.

ASSIGNMENT OF VALUES

The assigned mean values were derived from analyses of vials representative of the entire lot. Analyte values were obtained from laboratories of the instrument manufacturer, manufacturers of instrument specific reagents, or from reference laboratories.

The Expected Range of the Mean is provided to assist the laboratory until it has established its own mean and standard deviation. It is considered good laboratory practice for each laboratory to establish its own mean and standard deviation for its test methods. The indicated Mean and Expected Range of the Mean should serve as a guide in assessing the performance of each test method.

The values are usually method dependent. The variations which can occur over time and between laboratories may be attributed to differences in laboratory technique, instrumentation, reagent lot, method modifications, and other systematic errors including random errors. Analyte values for some methods may have been obtained using dilutions according to the manufacturer's instructions. The Expected Range of the Mean for each analyte value is calculated using the appropriate factors and may exceed the linearity for that method.

Refer to www.cliniqa.com/support/documents.aspx for lot specific value assignment. To request a faxed or emailed copy of the value assignment, contact your local distributor or csgroup@cliniqua.com.

CAUTION

Human source material. Treat as potentially infectious.

Each serum/plasma donor unit used in the manufacture of this product has been tested by FDA accepted methods and found non-reactive for the presence of HBsAg and antibody to HIV-1/2, HCV and HIV-1 Ag. While these methods are highly accurate, they do not guarantee that all infected units will be detected. Because no known test method can offer complete assurance the hepatitis B virus, hepatitis C virus, human immunodeficiency virus (HIV) or other infectious agents are absent, all products containing human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

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© 2008 Cliniqa Corporation. All rights reserved.