

---

# TRIPLE SUGAR IRON (TSI) AGAR

---

## INTENDED USE

Remel Triple Sugar Iron (TSI) Agar is a solid medium recommended for use in qualitative procedures for differentiation of microorganisms on the basis of dextrose, lactose, and sucrose fermentation and hydrogen sulfide production.

## SUMMARY AND EXPLANATION

In 1911, Russell described a medium containing dextrose and lactose for differentiation of enteric gram-negative bacilli.<sup>1</sup> Krumwiede and Kohn modified Russell's double sugar medium by adding sucrose which allowed for detection of sucrose-fermenting gram-negative bacilli.<sup>2</sup> In 1940, Sulkin and Willet further modified the medium by adding ferrous sulfate, an indicator of H<sub>2</sub>S production.<sup>3</sup> Hajna developed the recent formulation for Triple Sugar Iron Agar by adding sucrose to Kligler Iron Agar.<sup>4</sup>

## PRINCIPLE

Casein and meat peptones provide nitrogenous compounds, amino acids, and peptides necessary for bacterial growth. Dextrose, lactose, and sucrose are fermentable carbohydrates. Phenol red is an indicator of carbohydrate fermentation. Fermentation reactions are read on the slant and in the butt, with a color change from red (alkaline) to yellow (acid). The dextrose concentration in TSI Agar is one-tenth the concentration of lactose and sucrose which serves to distinguish dextrose-only fermenting organisms from those which ferment lactose and/or sucrose. The small amount of acid produced in the slant during dextrose fermentation oxidizes rapidly, causing the slant to revert to alkaline (red) while the yellow acid reaction is maintained in the butt due to the absence of oxygen. If the lactose or sucrose is also fermented, sufficient acid is produced to retain the yellow color on the slant and in the butt. Ferric ammonium citrate and sodium thiosulfate are indicators of H<sub>2</sub>S production which form a black precipitate (ferrous sulfate) in the butt of the tube. Gas production is indicated by bubbles, splitting of the agar, or displacement of the agar in the tube.

## REAGENTS (CLASSICAL FORMULA)\*

Casein Peptone.....	10.0 g	Dextrose .....	1.0 g
Lactose.....	10.0 g	Sodium Thiosulfate.....	0.2 g
Meat Peptone.....	10.0 g	Ferric Ammonium Citrate.....	0.2 g
Sucrose.....	10.0 g	Phenol Red.....	25.0 mg
Sodium Chloride.....	5.0 g	Agar.....	13.0 g
		Demineralized Water.....	1000.0 ml

pH 7.3 ± 0.2 @ 25°C

\*Adjusted as required to meet performance standards.

## PRECAUTIONS

This product is For Laboratory Use only. It is not intended for use in the diagnosis of disease or other conditions.

## PREPARATION OF DEHYDRATED CULTURE MEDIUM

1. Suspend 59.5 grams of medium in 1000 ml of demineralized water.
2. Heat to boiling with agitation to completely dissolve.
3. Dispense into tubes and sterilize by autoclaving at 121C for 15 minutes.
4. Cool in a slanted position so that deep butts are formed.

## PROCEDURE

1. Consult current editions of appropriate references for the recommended procedure for sample preparation, inoculation, testing, and interpretation.

## INTERPRETATION OF THE TEST

### Carbohydrate Fermentation:

Positive Test for Slant Reaction - Yellow (acid = A)  
Negative Test for Slant Reaction - Red (alkaline = ALK)

Positive Test for Butt Reaction - Yellow (acid = A)  
Negative Test for Butt Reaction - Red (alkaline = ALK)

ALK/A - Dextrose only fermented  
A/A - Dextrose and/or sucrose and/or lactose fermented  
ALK/ALK - No sugar fermented

### Hydrogen Sulfide Production:

Positive Test - Black color throughout the entire butt (may mask acidity), a black ring at the juncture of the slant and butt, or a black precipitate  
Negative Test - No blackening of medium

### Gas Production (CO<sub>2</sub> and H<sub>2</sub>):

Positive Test - Bubbles in the medium, cracking and displacement of the medium, or separation of the medium from the side and bottom of the tube  
Negative Test - No bubbles and no separation or displacement of the medium

## QUALITY CONTROL

Each lot number of Triple Sugar Iron Agar has been manufactured, packaged, and processed in accordance with current Good Manufacturing Practice regulations. All lot numbers have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, sample results should not be reported.

CONTROL	INCUBATION	RESULTS			
		GAS	SLANT	BUTT	H <sub>2</sub> S
<i>Escherichia coli</i> ATCC® 25922	Ambient, 18-24 h @ 33-37°C	+	A	A	-
<i>Salmonella enterica</i> serovar Typhimurium ATCC® 14028	Ambient, 18-24 h @ 33-37°C	+	ALK	A	+
<i>Pseudomonas aeruginosa</i> ATCC® 27853	Ambient, 18-24 h @ 33-37°C	-	ALK	ALK	-

## LIMITATIONS

1. To enhance the alkaline condition in the slant a free exchange of air must be permitted. If TSI tubes are tightly capped, an acid reaction caused solely by dextrose fermentation will also involve the slant. Therefore, tubes must have loosened caps during incubation.<sup>5</sup>
2. Before inoculation, a slight precipitate may be present on the slant. This will not effect the performance of the medium.<sup>6</sup>
3. Studies have demonstrated that sucrose utilization may suppress the enzyme mechanisms responsible for H<sub>2</sub>S production.<sup>7</sup>
4. SIM (sulfide-indole-motility) agar has been reported to be more sensitive for the detection of H<sub>2</sub>S than TSI.<sup>5</sup>
5. Some organisms, such as *Proteus* spp., may produce reactions similar to *Salmonella* and *Shigella*. Additional biochemical testing is required for definitive identification. Consult appropriate references for further instructions.<sup>8,9</sup>

## BIBLIOGRAPHY

1. Russell, F.F. 1911. J. Med. Res. 25:217-229.
2. Krumwiede, C. and L. Kohn. 1917-1918. J. Med. Res. 37:225-227.
3. Sulkin, S.E. and J.C. Willett. 1940. J. Lab. Clin. Med. 25:649-653.
4. Hajna, A.A. 1945. J. Bacteriol. 49:516-517.
5. MacFaddin, J.F. 1985. Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria. Vol. 1. Williams & Wilkins, Baltimore, MD.
6. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1995. Manual of Clinical Microbiology. 6<sup>th</sup> ed. ASM, Washington, D.C.
7. Bulmash, M.J. and M. Fulton. 1964. J. Bacteriol. 88:1813.
8. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10<sup>th</sup> ed. ASM Press, Washington, D.C.
9. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12<sup>th</sup> ed. Mosby Elsevier, St. Louis, MO.

Refer to the front of Remel *Technical Manual of Microbiological Media* for **General Information** regarding precautions, product storage and deterioration, sample collection, storage and transportation, materials required, quality control, and limitations.

ATCC® is a registered trademark of American Type Culture Collection.  
IFU 454982, Revised December 1, 2011

Printed in U.S.A.

**remel**

12076 Santa Fe Drive, Lenexa, KS 66215, USA  
General Information: (800) 255-6730 Website: [www.remel.com](http://www.remel.com) Email: [remel@remel.com](mailto:remel@remel.com)  
Local/International Phone: (913) 888-0939 International Fax: (913) 895-4128