



MedSchenker® STM®: Smart Transport Medium®

For collection, transport, maintenance and long term freeze storage of clinical specimens containing viruses

Product Information and How-to-Use

INTENDED USE

The MedSchenker STM is intended for the collection and transport of clinical specimens containing viruses from the collection site to the testing laboratory. The STM can be processed using standard clinical laboratory operating procedures for viral cultures.

SUMMARY AND EXPLANATION

One of the routine procedures in the diagnosis of infections caused by viruses, involves the collection and safe transportation of biological samples. This can be accomplished using the MedSchenker STM. The STM includes a transport medium that is room temperature stable which can sustain viability (and infectivity) of a plurality of organisms that include clinically important viruses during transit to the testing laboratory. The formulation of the STM includes proteins for stabilization, antibiotics to minimize bacterial and fungal contamination and a buffer to maintain a neutral pH.

The STM is provided in a labeled, screw-cap tube designed for transport of the clinical sample. The STM is also supplied as a sample collection kit that comprises a package which contains one screw-cap tube of STM and one or two sterile specimen collection CavSwab™ flocked swabs.

Once a CavSwab flocked swab sample is collected it should be placed immediately into the transport tube where it comes into contact with the transport medium. Swab specimens for virus isolation should be submitted to the laboratory as quickly as possible after collection. Although the STM can maintain fragile organisms for long periods of time at room temperature, it is recommended that specimens be refrigerated at 2 – 8°C or kept on wet ice following collection and while in transit. If there will be a long delay before processing, specimens should be frozen at – 70°C or colder and transported on dry ice. Storage at –20°C is less satisfactory than storage at 4°C or –70°C and can result in the loss of infectivity.

PRINCIPLE

The STM consists of modified Hank's balanced salt solution supplemented with fetal bovine serum, cysteine, gelatin, sucrose and glutamic acid. The pH is buffered with HEPES buffer. Phenol red is used to indicate pH. Vancomycin, amphotericin B, and colistin are incorporated in the medium to inhibit growth of competing bacteria and yeast. The medium is isotonic and non-toxic to mammalian host cells. The presence of sucrose acts as a cryoprotectant which aids in the preservation of viruses if specimens are frozen (–70°C) for prolonged storage.

PRECAUTIONS

- This product is For In Vitro Diagnostic Use.
- Observe approved bio-hazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
- All specimens and materials used to process them should be considered potentially infectious and handled in a manner which prevents infection of laboratory personnel. Sterilize all bio-hazard waste including specimens, containers and media after their use.
- Directions should be read and followed carefully.

STORAGE

This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 2 – 25°C until used. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual sterile pouch and the specimen transport tube label.

SHELF LIFE

The MedSchenker STM has a shelf life of 18 months after the date of manufacture. Do not use after the listed expiration date.

472.441

PRODUCT DETERIORATION

Unless expressly approved by the manufacturer, the STM should not be used if (1) there is evidence of damage or contamination to the product (2) there is evidence of leakage (3) there is evidence of breakage, tear, or damage to the packaging container (4) the color of the medium has changed from reddish pink - magenta (5) the expiration date has passed (6) the swab pouch is open (7) there are other signs of deterioration.

SPECIMEN COLLECTION, STORAGE AND TRANSPORTATION

Specimens for virus investigation should be collected and handled following published manuals and guidelines 2,3,4,7,9,10,11. To maintain optimum viability, transport the specimen to the laboratory as soon as possible. Best recovery is obtained when specimens are refrigerated at 2 – 8°C or kept on wet ice following collection and while in transit. If there will be a long delay before processing, specimens should be frozen at –70°C or colder and transported on dry ice. Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations 8,11,12 Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All Specimens should be processed as soon as they are received in the laboratory.

MATERIALS SUPPLIED

MedSchenker STM comes in a screw-cap tube containing 1.5 / 2 / 3 ml of transport medium.

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for isolating, differentiating and culturing viruses. These materials include tissue culture cell lines, tissue culture medium, incubation systems and reading equipment. Refer to appropriate references for recommended protocols for isolation and identification of viruses 2,3,4,7,10.

DIRECTIONS FOR USE

MedSchenker STM is available in the product configurations indicated in the table below

SKU	STM Tube Description	Pack Size	Sampling Sites*
STM15-A1	1.5 mL Screw Cap with Tube	50 qty	Culture: viruses
STM20-A	2.0 mL Screw Cap with Tube	50 qty	Culture: viruses
STM30-A1	3.0 mL Screw Cap with Tube	50 qty	Culture: viruses

Table 1. Comparison Results of MedSchenker STM and COPAN UTM for detection of Influenza A H1N1

	Influenza A H1N1		Influenza A		Influenza A H1N1		Influenza A		Influenza A H1N1		Influenza A		Influenza A H1N1		Influenza A	
	STM	UTM	STM	UTM	STM	UTM	STM	UTM	STM	UTM	STM	UTM	STM	UTM	STM	UTM
01	P	P	P	P	N	N	N	N	N	N	P	P	N	N	P	P
02	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N
03	N	N	N	N	N	N	N	N	N	N	P	P	N	N	P	P
04	P	P	P	P	N	N	N	N	N	N	N	N	N	N	N	N
05	N	N	N	N	P	P	P	P	N	N	N	N	N	N	N	N
06	N	N	N	N	N	N	N	N	N	N	P	P	N	N	P	P
07	P (32)	N	P	N	N	N	N	N	N	N	N	N	N	N	N	N
08	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N
09	N	N	N	N	P	P	P	P	P	P	N	N	P	P	N	N
10	N	N	N	N	N	N	N	N	N	N	P	P	N	N	P	P
11	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N
12	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N
13	N	N	N	N	N	N	N	N	N	N	P	P	N	N	P	P
14	P	P	P	P												

7th specimen was positive after nucleic sequence diagnosis. Agreement (%) of MedSchenker STM and COPAN UTM for detection of Influenza A H1N1 was 98.11% (52/23).

* Performance testing with the MedSchenker STM was conducted using laboratory strains spiked onto a swab and not using human specimens.

Specimens should be collected as soon as possible after the clinical onset of disease. Highest viral titers are present during the acute illness For STM15-A / STM20-A / STM30-A Tubes

- Aseptically remove cap from tube
- Aseptically place vesicle aspirates, corneal or conjunctival scrapings, small pieces of tissue or stool samples into the tube with STM
- Replace cap to tube and close tightly
- Label with appropriate patient information
- Send to the laboratory for immediate analysis

QUALITY CONTROL

All lot numbers of the STM are tested for microbial contamination, toxicity to host cells and the ability to maintain viability of desired agents. Procedures for quality control of STM and virus culture media are described in several publications by the American Society for Microbiology 3,7,10 and by NCCLS 5,6. If aberrant quality control results are noted, patient results should not be reported.

Table 2. Effect of storage time for the cycle of threshold (Ct) of Influenza A H1N1 Positive specimens

Ct of Influenza A H1N1			Ct of Influenza A		
24 hrs	3 Days	7 Days	24 hrs	3 Days	7 Days
21.67	21.99	22.74	22.20	21.76	22.38
23.10	23.55	24.89	21.67	22.20	23.60
27.60	24.66	24.68	27.63	25.45	Within
19.39	19.85	19.96	19.44	20.20	20.22
32.36	32.38	32.74	32.22	31.86	32.89

Table 3. Effect of storage time for delta RN of Influenza A H1N1 positive specimens

Ct of Influenza A H1N1			Ct of Influenza A		
24 hrs	3 Days	7 Days	24 hrs	3 Days	7 Days
4336	14,204	18,788	13,620	13,410	13,998
9837	10,899	13,851	12,693	12,030	13,269
8972	9,626	15,809	9,733	9,731	13,703
6587	15,614	18,731	14,365	13,621	14,940
6705	4,491	6,750	4,515	6,023	5,589

LIMITATIONS

- Specimens should be handled aseptically.
- Condition, timing and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.1,2,3,4,7,10.
- Repeated freezing and thawing of specimens may reduce the recovery of viable organisms.
- The STM is intended for use as a collection and transport medium for viral agents only. This medium can serve as a cryoprotectant for clinical viruses.
- Because calcium alginate swabs are toxic for many enveloped viruses and may interfere with fluorescent antibody tests, they should not be used for specimen collection. Wooden shaft swabs may contain toxins and formaldehydes and should not be used. Polyester (Dacron) tipped swabs and flocked swabs are suitable when specimen collection by a swab is appropriate.
- The STM is intended to be used with MedSchenker CavSwab flocked swabs. The use of swabs from any other source could affect the performance of the product.

WARNINGS

- Do not re-pack
- Not suitable to collect and transport microorganisms other than viruses
- Not suitable for any other application than intended use
- The use of this product in association with a rapid diagnostic kit or with diagnostic instrumentation should be previously validated by the user
- Do not ingest the STM
- To be handled by trained personnel only
- Do not use the STM for pre-moistening or pre-wetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites

DISCLAIMER

The manufacturer of MedSchenker STM shall not be liable for any losses and/or damages caused by the user's sole negligence or intentional tort on using STM or non-compliance with all manuals and guidelines of STM, including but not limited to, laboratory errors, test errors, damage to samples or specimens, laboratory contamination, infection, or other similar losses or damages.

RESULTS

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

- Currently undergoing further testing and analysis to provide a full spectrum of virus transportation data.

The MedSchenker STM was able to maintain the viability of the following organisms for at least 48 hours at both room temperature (20 - 25°C) and in the refrigerator (2 - 8°C) under the test conditions described above.

- STM is a transport medium that serves as a culture media non-propagating transport. It has not been reviewed by the FDA
- Specimen stability for this media was not validated for recovery of viral infectious particles using a culture-based assay

Manufacturer
Korea Standard
Jeollanamdo, South Korea

MedSchenker
493 Old Trenton Road
East Windsor NJ 08512 USA
Exclusive Distributor



STM-UNV-0002-Q122