

HemoTrol

Superior QC for HemoCue B-Hemoglobin & Hb 201 Systems



No other quality control offers so many benefits:

- **18 month refrigerated shelf-life** – ensures long term QC continuity and virtually eliminates product waste due to outdating
- **30-day open vial stability at room temperature** – makes QC more convenient and accessible, eliminates waste due to accidental room temperature exposure
- **Hemoglobin matrix** – closely mimics whole blood, outperforms ‘pickled’ or stabilized RBC formulations
- **Optimized fill volume** – 20 tests per vial reduces risk of using expired product
- **Recommended and distributed by HemoCue**

Best of all, HemoTrol is competitively priced and available from stock for immediate delivery!

<u>Level</u>	<u>Target</u>	<u>Packaging</u>
Low	8.0 g/dL	2 vials/box
Normal	12.0 g/dL	2 vials/box
High	16.0 g/dL	2 vials/box