

Hyperbaric QC

For Validation of elevated pO_2 Blood Gas Measurements



Ideal for Analytical Measurement Range (AMR) Validation:

- Access accuracy and precision of elevated pO_2 measurements
- Extend your lab's reportable range
- Allows reporting of actual results--simplify recordkeeping
- Satisfy federal and accrediting agency AMR requirements
- Ideal for cardiac catheterization labs & OR suites
- Three year shelf-life at room temperature (15 - 30°C)

<u>Level</u>	<u>Target</u>	<u>Packaging</u>
High	460mmHg	4 amps/box
X-High	710mmHg	4 amps/box