

# BinaxNOW® RSV Test Kit

Distributed by

CLIAwaived™ Inc.



The BinaxNOW® RSV Test Kit is a rapid immunochromatographic assay for the qualitative detection of respiratory syncytial virus (RSV) fusion protein antigen in nasal wash and nasopharyngeal (NP) swab specimens from symptomatic patients. The test is intended for in vitro diagnostic use to aid in the diagnosis of RSV infections in neonatal and pediatric patients under the age of five. It is recommended that negative test results be confirmed by cell culture.

#### BENEFITS:


- Onboard extraction: no extra steps
- Sample types: nasopharyngeal swab and nasal wash
- Shelf life up to 24 months
- One step test—no additional reagents needed
- FDA-Cleared / CLIA Waived

Ordering Information:	
Product Description	Catalog #
BinaxNOW® RSV Test Kit (22 tests)	INV-430-122
BinaxNOW® RSV Test Kit (10 tests)	INV-430-100
Reimbursement Information:	
CPT Code #	National Limit Amount
87807QW - Rapid immunoassay for the qualitative detection for respiratory syncytial virus antigen.	\$17.08

DISCLAIMER – The CPT information contained in this document is provided as representative examples of reimbursement in this category. It is intended to assist providers in accurately obtaining reimbursement for health care services. It is not intended to increase or maximize reimbursement by any payor. Providers should consult their payor organizations with regard to local reimbursement policies. The information provided in this document is for information purposes only and represents no statement, promise or guarantee by CLIAwaived™ Inc. All CPT codes are supplied for information purposes only and represent no statement; promise or guarantee by CLIAwaived™ Inc. that these codes will be appropriate or that reimbursement will be made. CPT codes and descriptions are copyrights of the American Medical Associations. CPT does not include fee schedules, relative values or related ratings. The source for this information is the Center for Medicare and Medicaid Services. The content provided by the Center for Medicare and Medicaid Services is updated frequently. It is the responsibility of the service provider to confirm the appropriate coding required by their local Medicare carriers, fiscal intermediaries and commercial payors. All product information and prices are subject to change without notice at any time. Please call our office before ordering to verify current prices.

# BinaxNOW® RSV Test Kit Instructions

## Test Procedure:



Rev. 1.24.11