

UrineCheck 7 Drug Adulteration Test Strips

Distributed by

CLIAwaived™ Inc.



Note: The UrineCheck 7 Drug Adulteration Test Strips are not for use in diagnosing disease or illness to effect treatment, and /or cure, it is for forensic/ toxicology use only.

The UrineCheck 7 Drug Adulteration Test Strips are a fast dip-and-read test for the determination of diluted or adulterated urine specimens. It is an important pre-screening test for any drug-testing program. The UrineCheck 7 Drug Adulteration Test Strips include assays for Creatinine, pH, Specific Gravity, Nitrite, Bleach, Pyridinium Chlorochromate and Glutaraldehyde. *(The test strip is used to determine specimen urine integrity and is for forensic use only.)*

BENEFITS:

- Accurate results in 30 seconds to 2 minutes
- Room temperature storage
- Ergonomic easy to read color chart
- Available in 50 or 100 strips/bottle sizes
- Ready-to-use and disposable
- Fast dip-and-read test

| Ordering Information: | |
|--|-----------------------|
| Product Description | Catalog # |
| UrineCHECK 7 Adulterant Tests (25 test strips) | TC-D700-25 |
| Amtox UrineCheck 7 Adulterant Controls (2 bottles - 25ml each) | AMTX-URINECHECK7C |
| Reimbursement Information: | |
| CPT Code # | National Limit Amount |
| N/A | N/A |

All product information and prices are subject to change without notice at any time. Please call our office before ordering to verify current prices.

Rev. 03.16.10