

Product Information

Product Name:	COVID 19 Negative Control Set (2x3mL)
Product Number:	COV-19-NEG-CTRL-SET
Date of Manufacture:	07/20
Lot Number:	C5652
Expiration Date:	2022-01-31
Intended Use:	In vitro diagnostic
Homogeneity:	Sufficiently homogeneous for analytical testing

Precautions

Matrix:	SMx Serum
Form:	Frozen
Anticoagulant:	N/A
Preservative:	N/A
pH:	6.81 – 6.96
Virology:	Non-reactive for HBsAg, HCV, and HIV
Microbial Growth:	<100 colony-forming units per milliliter

Storage and Stability

Storage Condition:	-10°C (14°F)
Open Vial Stability:	30 days*
Vial Volume:	2 mL

*The matrix may be thawed, aliquoted, and refrozen; although, this product should not be exposed to multiple freeze/thaw cycles.

Origin of Raw Material

UTAK Laboratories

Analyte	Supplier	Lot#	Purity
L-Alanine	Sigma-Aldrich	BCBG7593V	>98%
Recombinant Albumin-Cellastim (rHSA)	InVitria	P0689	98%
L-Arginine	Sigma-Aldrich	MKCG4020	100%
Calcium Chloride	Sigma-Aldrich	110M0005V	100%
Magnesium Chloride	Sigma-Aldrich	021M01071V	98%
Sodium Chloride	Sigma-Aldrich	SLBZ8529	99%
Potassium Chloride	Sigma-Aldrich	SLBR9127V	99%
Citric Acid	Spectrum Chemical	UD0899	98%
Creatinine	Sigma-Aldrich	SLBW4791	99%
D-(+)-Glucose	Sigma-Aldrich	SLBG4983V	99%
L-Glycine	Sigma-Aldrich	SLBF4675V	100%
L-Lysine	Sigma-Aldrich	BCBH1450V	100%
Potassium Phosphate Monobasic	Sigma-Aldrich	SLBD0340V	99%
L-Proline	Sigma-Aldrich	SLBL4766V	100%
Pyruvic Acid	Sigma-Aldrich	MKBF6056V	98%
Sodium Bicarbonate	Sigma-Aldrich	SLBL0544V	99%
Sodium Phosphate Dibasic	Sigma-Aldrich	SLBT4651	99.8%
Sodium Sulfate	Sigma-Aldrich	MKBG8705V	97%
L-Threonine	Sigma-Aldrich	SLBK2255V	100%
Urea	Sigma-Aldrich	SLBT1918	100%
Uric Acid	Sigma-Aldrich	STBB7849V	98%
L-Valine	Sigma-Aldrich	SLBZ5498	100%

Approved by: Kevin Kopp
Quality Assurance

