

**NOW  
DETECTS  
BUP!**

## CLIAwaived,™ Inc. Instant Drug Test Card "IDTC"

The CLIAwaived, Inc IDTC is a one-step dip card for detection for drugs of abuse urine sample. When quality and cost is your goal, this test card is the answer. The IDTC test card is designed for simplicity of use and economical cost.



- One step qualitative assay for over 13 unique targets
- Read negative results as soon as 1 minute
- Stable results up-to 1 hour
- Room temperature storage

### SINGLE DIP DRUG TEST

1-PANEL SDDT		
Tests for:	Sensitivity (ng/mL)	Catalog #
<b>Cocaine (COC)</b>	150	CLIA-IDTC-II-11
<b>Opiates (OPI)</b>	300	CLIA-IDTC-II-13
<b>Methamphetamine (MET)</b>	500	CLIA-IDTC-II-16
<b>Benzodiazepines (BZO)</b>	300	CLIA-IDTC-II-18
<b>Barbiturates (BAR)</b>	300	CLIA-IDTC-II-19
<b>Methadone (MTD)</b>	300	CLIA-IDTC-II-20
<b>Oxycodone (OXY)</b>	100	CLIA-IDTC-II-23
<b>Buprenorphine (BUP)</b>	10	CLIA-IDTC-II-25
<b>Methadone Metabolite (EDDP)</b>	300	CLIA-IDTC-II-30



**CLFS National Limit: \$14.86 - Code 80305 - Drug test(s) presumptive, any number of drug classes, any number of devices or procedures capable of being read by direct optical observation only (eg, utilizing immunoassay [eg, dipsticks, cups, cards, or cartridges]), includes sample validation when performed, per date of service. The code description includes "per date of service," therefore it should only be billed with 1 unit per date of service.**

DISCLAIMER: The CPT information contained in this document is provided as representative examples of reimbursement in this category. It is intended to assist providers in accurately obtaining reimbursement for healthcare services. It is not intended to increase or maximize reimbursement by any payor. Providers should consult their payor organizations with regard to local reimbursement policies. The information provided in this document is for information purposes only and represents no statement, promise or guarantee by CLIAwaived,™ Inc. All CPT codes are supplied for information purposes only and represent no statement, promise or guarantee by CLIAwaived,™ Inc. that these codes will be appropriate or that reimbursement will be made. CPT codes and descriptions are copyrights of the American Medical Association. CPT does not include fee schedules, relative values or related ratings. The source for this information is the Center for Medicare and Medicaid Services. The content provided by the Center for Medicare and Medicaid Services is updated frequently. It is the responsibility of the service provider to confirm the appropriate coding required by their local Medicare carriers, fiscal intermediaries and commercial payors. All product information and prices are subject to change without notice at any time. Please call our office before ordering to verify current prices.

\*NOTE: There are no uniformly recognized cutoff levels for drugs of abuse testing for Benzodiazepine, Oxycodone, Tricyclic Antidepressant and Barbiturates. The test is not intended to screen or monitor individuals who are prescribed these drugs by a physician; the test may yield positive results for individuals taking such drugs, as prescribed.

To order, call 888.882.7739 or visit [CLIAwaived.com](http://CLIAwaived.com)