

# On·Call® Pro

## Blood Glucose Monitoring System

Cleared for multi-patient use in LTC Facilities

### Fast and Reliable

- Clinically proven accuracy
- Auto calibration circuitry in every test strip
- Audible fill detection and redosing to avoid wasted strips
- 24/7 Support: 1-800-838-9502

### Less Pain


- Tiny 0.4 µL blood sample
- Alternate site testing (Palm, Forearm)

### Easy

- 4 second testing
- No coding required
- User friendly design
- Large display with bold numbers
- Convenient test strip ejector reduces risk of contamination
- Strip port faces the patient

### Advanced Features

- Daily Control Test Reminder option
- Hypo and hyper alarms
- Ketone test alert



| Feature                     | Specification  |
|-----------------------------|--|
| Measurement Range           | 20 - 600 mg/dl   |
| Result Calibration          | Plasma-equivalent  |
| Minimum Sample Size         | 0.4 µL   |
| Test Time                   | 4 seconds  |
| Battery Life                | 1,000 tests  |
| Glucose Units of Measure    | The meter is pre-set to milligrams per deciliter (mg/dl) |
| Memory                      | Up to 300 records with time and date                     |
| Automatic shutoff           | 2 minutes after last action or 5 sec after strip removal |
| Meter Size                  | 85.0 mm x 49.0 mm x 16.4 mm                              |
| Display Size                | 35 mm x 32 mm  |
| Operating Temperature       | 41-113°F (5-45°C)  |
| Operating Relative Humidity | 10 - 90% (non-condensing)                                |
| Hematocrit Range            | 30-55%   |
| Data Port                   | USB  |

CLIAwaived™ Inc.