



HANGZHOU GENESIS BIODETECTION & BIOCONTROL CO.,LTD.

No.139, St.10th (East), Hangzhou Economic & Technology
Development Zone Zhejiang Province, China

Tel : +86-571-86736936 FAX : +86-571-87824695

PRODUCT SHELF LIFE EXTENSION

June 30, 2021

Dear Valued Customer :

Hangzhou GENESIS Biodetection and Biocontrol CO.,LTD. determines the shelf life of our products after a carefully assessment. The KaiBiLi™ Extended ViralTrans 3mL Viral Transport Medium (Cat. No. M221006) was tested to be stabled for 14 months with appropriate transportation and storage condition. For the best product performance, we set 12 months as the labeled shelf life. For product beyond 12 months to 14 months, we would suggest users to perform QC test before applying the goods to general practice to validate the product quality.

For more information, please contact Hangzhou GENESIS Biodetection and Biocontrol CO.,LTD. at customer@hgb.com.cn. Thank you very much.

Most sincerely,

A handwritten signature in black ink, reading "Reuben Weng", written over a horizontal line.

Reuben Weng

Vice President, International Business