



Biochemical Diagnostics, Inc.

180 Heartland Blvd, Edgewood, NY 11717 • Phone (800) 223-4835

Fax (631) 595-9204 • www.biochemicaldiagnostics.com

SALIVABUSE® LIQUID CONTROL, ORAL-SKREEN



Constituents (ng/mL)	Negative	3X Positive
Cannabinoids (parent)	0	120
Cocaine (parent)	0	90
Benzoylcegonine	0	90
Phencyclidine (PCP)	0	30
Opiates/morphine	0	120
Oxycodone	0	90
Amphetamines	0	150
MDMA	0	150
Methamphetamines	0	150
Barbiturates (Secobarbital)	0	150
Benzodiazepines (Oxazepam)	0	150
Methadone	0	90
Cotinine	0	90
Ethanol (mg/dL)	0	75
Buprenorphine	0	30

Kitbox part numbers for distribution

702030	Oral Fluid Control Kitbox, Neg/3X POSITIVE, 2x20mL
20200101-6pk	Negative, 6x20mL
20200109-6pk	Oral Fluid, 3X Cutoff, 6x20mL

The 3X Positive and Negative controls are designed to work with most DOA devices on the market as an external control for customer convenience. Oral Fluid controls can be ordered in +/-50% or other custom levels as needed.

Please contact Colleen Gang to discuss our new Salivabuse Controls, they are available for distribution to work with your oral fluid devices in stock or custom formulations
Email: colleen@biochemicaldiagnostics.com phone: 631-637-2024.