



Nitrile Chemo Exam Gloves



Features and Benefits

- Medical Grade
- Powder Free
- Chemotherapy rated
- Ambidextrous
- Eliminates exposure to latex
- Comfortable fit and feel
- Textured surface for improved grip
- Fatigue free movement and unparalleled dexterity

Specifications



Ordering Information

ITEM #	SIZE	INNER BOX	OUTER CASE	QUANTITY
A18A11	SM	8-10052-97055-4	008-10052-97056-1	100bx/10cs
A18A12	MD	8-10052-97058-5	008-10052-97059-2	100bx/10cs
A18A13	LG	8-10052-97061-5	008-10052-97062-2	100bx/10cs
A18A14	XL	8-10052-97064-6	008-10052-97065-3	100bx/10cs

INNER BOX DIMENSIONS (mm)	L 240 x W 120 x H 63
OUTER CASE DIMENSIONS (mm)	L 345 x W 258 x H 235

Finding the right hand protection doesn't have to be complicated. Consider VGuard® gloves as an extension of our promise to you: High quality products, great value and excellence in service. If ever you have a question regarding which glove is best for your application, we're here to help. VGuard® is protection you can count on!

Vanguard International Solutions, Inc.

Cumming, GA 30041 • 1-866-VGUARD-6 (1-866-848-2736)
info@vanguardsourcing.com • www.vanguardsourcing.com

Usage



Dental



Healthcare



Laboratory



Home



Foodservice



Manufacturing



Vet



Janitorial



Chemotherapy



Nitrile Chemo Exam Gloves

Composition

Material:	Nitrile synthetic rubber
Allergen content:	Contains no proteins found in latex
Powdered:	No, as per ASTM D6124

Quality

Quality Control:	ASTM D6319, ASTM F1671 ASTM D6978-05 (Chemo)
Quality Assurance:	ISO:9001, GMP, HACCP
Regulatory Compliance:	FDA approved

Technical Data

THICKNESS	MIL
Palm:	3.0 ± 0.5
Finger:	4.0 ± 0.5

Weight*:	4.0 gm ± 0.2
-----------------	--------------

*Based on size M

PHYSICAL PROPERTIES	BEFORE AGING	AFTER AGING
Tensile Strength: (MPa)	Min. 14	Min. 14
Ultimate Elongation (%):	Min. 500	Min. 400
Freedom from holes:	AQL 1.5	N/A

DIMENSIONS			
SIZE	ITEM#	Palm Width (mm)	Length (mm)
Small	A18A11	85 ± 5	Min. 240
Medium	A18A12	95 ± 5	Min. 240
Large	A18A13	105 ± 5	Min. 240
Extra Large	A18A14	115 ± 5	Min. 240

Vanguard International Solutions, Inc.

Cumming, GA 30041 • 1-866-VGUARD-6 (1-866-848-2736)
info@vanguardsourcing.com • www.vanguardsourcing.com