

# N300 Series Nitrile Exam Gloves

## Features and Benefits

- Nitrile
- Thickness (Finger): 5 mil
- Powder Free
- ASTM D6319, D6978
- FDA 21 CFR 177.2600
- Exam Grade
- Color: Cobalt
- Texture: Finger and Palm
- Length: Standard (240mm)
- Cuff: Beaded
- Non-Sterile
- Chemotherapy Tested<sup>1</sup>

<sup>1</sup>Chemotherapy drug testing information on back of box.



## Usage

Dependable for house-wide applications within Healthcare, Labs, EMS and other settings where high-use and disposable exam glove protection is important. Full palm texture provides excellent grip.

### Specifications



### Usage



## Ordering Information

ITEM #	SIZE	GLOVES PER BOX	BOXES PER CASE	GLOVES PER CASE
N300S	Small	100	10	1000
N300M	Medium	100	10	1000
N300L	Large	100	10	1000
N300XL	X-Large	100	10	1000

# N300 Series Nitrile Exam Gloves

100 CT



## Spec Information

### OVERVIEW

Material	Nitrile
Powder Content	Powder-Free
Freedom From Holes	1.5 AQL
Surface	Full hand textured (Finger and Palm)
Color	Cobalt Blue
Cuff Length	240mm, Standard
Not Made With Natural Rubber Latex	Yes
Chemo Rated (Yes/No)	Yes, in accordance with FDA ASTM D6978
Food Rated (Yes/No)	Yes, in accordance with FDA 21 CFR 177.2600
Weight	4.2g +/- 0.2 (size S) 4.7g +/- 0.2 (size M) 5.2g +/- 0.2 (size L) 5.7g +/- 0.2 (size XL)

### THICKNESS MEASUREMENTS

Finger	0.10mm - 0.14mm
Palm	0.08mm - 0.12mm
Cuff	0.07mm - 0.10mm

### PHYSICAL PROPERTIES

Criteria	Before Aging	After Aging
Tensile Strength	>= 14	>= 14
Elongation	500 (min.)	400 (min.)

### DIMENSIONS

Item#	Size	Palm Width (mm)	Length
N300S	Small	85 ± 5	240mm / 9.5 inches
N300M	Medium	95 ± 5	240mm / 9.5 inches
N300L	Large	105 ± 5	240mm / 9.5 inches
N300XL	Extra-Large	115 ± 5	240mm / 9.5 inches