



## Vinyl Exam Gloves



### Features and Benefits

- Clear Color
- Powder Free
- Medical Grade
- Single Use Only
- Smooth Surface
- Approved for Food Handling
- Prop 65 Compliant

### Specifications



### Ordering Information

ITEM #	SIZE	INNER BOX	OUTER CASE	QUANTITY
A21A11	SM	8-50001-98101-9	008-50001-98150-7	100bx/10cs
A21A12	MD	8-50001-98103-3	008-50001-98104-0	100bx/10cs
A21A13	LG	8-50001-98106-4	008-50001-98107-1	100bx/10cs
A21A14	XL	8-50001-98109-5	008-50001-98110-1	100bx/10cs

INNER BOX DIMENSIONS (mm)	L 230 x W 120 x H 65
OUTER CASE DIMENSIONS (mm)	L 340 x W 248 x H 244

Finding the right hand protection doesn't have to be complicated. Consider VGuard® gloves as an extension of our promise to you: High quality products, great value and excellence in service. If ever you have a question regarding which glove is best for your application, we're here to help. VGuard® is protection you can count on!

## Vanguard International Solutions, Inc.

Cumming, GA 30041 • 1-866-VGUARD-6 (1-866-848-2736)  
info@vanguardsourcing.com • www.vanguardsourcing.com

### Usage



Dental



Healthcare



Laboratory



Home



Foodservice



Manufacturing



# Vinyl Clear Exam Gloves

## Composition

<b>Material:</b>	Poly Vinyl chloride
<b>Allergen content:</b>	Contains no proteins found in latex
<b>Powdered:</b>	No, as per ASTM D6124

## Quality

<b>Quality Control:</b>	ASTM D5250, ASTM F1671
<b>Quality Assurance:</b>	ISO:9001, GMP, HACCP
<b>Regulatory Compliance:</b>	FDA approved

## Technical Data

THICKNESS	MIL
<b>Palm:</b>	3.0 ± 0.5
<b>Finger:</b>	4.0 ± 0.5

<b>Weight*:</b>	5.0 gm ± 0.2
-----------------	--------------

*\*Based on size M*

PHYSICAL PROPERTIES	BEFORE AGING	AFTER AGING
<b>Tensile Strength: (MPa)</b>	Min. 11	Min. 11
<b>Ultimate Elongation (%):</b>	Min. 300	Min. 300
<b>Freedom from holes:</b>	AQL 2.5	N/A

DIMENSIONS			
SIZE	ITEM#	Palm Width (mm)	Length (mm)
Small	A21A11	80 ± 10	Min. 220
Medium	A21A12	95 ± 10	Min. 230
Large	A21A13	110 ± 10	Min. 230
Extra Large	A21A14	120 ± 10	Min. 230

## Vanguard International Solutions, Inc.

Cumming, GA 30041 • 1-866-VGUARD-6 (1-866-848-2736)  
info@vanguardsourcing.com • www.vanguardsourcing.com