


LUCIRA™ CHECK IT COVID-19 TEST KIT



- ✓ Easy to use
 - 3 simple steps: Swab, Stir, Detect
 - 100% of users were able to run the test¹
- ✓ Results in 30 minutes or less
 - Detects positives in as fast as 11 minutes and confirms negatives in 30 minutes.
- ✓ Clinically Proven
 - 98% accurate compared to one of the best lab PCR tests²

Lucira Community Testing Study Positive Percent Agreement (PPA) Summary


PPA = 92% across all data (n=132) PPA = 97% on all CT Values ≤ 37.5 (n=122)



Samples shown by decreasing amounts of virus

The # of cycles (CTs) required to detect virus increases when the amount of virus in the sample decreases.

- Lucira’s EUA-authorized LOD is comparable to high sensitivity PCR assays where LOD is 1000 cps/ml or less



✓ Able to detect the Alpha, Beta, Delta, Gamma, Lambda, Mu and Omicron variants³

✓ Current Authorizations:

- USA
- Canada
- Israel
- Singapore
- Taiwan

✓ Some organizations we have partnered with:



¹ Lucira Usability Clinical Study 07A-CLI-004/005 conducted Summer 2020 (n=352)

² Excluding samples with very low levels of virus (>37.5 Ct), Lucira Community Testing Study 07A-CLI 006 Fall 2020 and Lucira Community Testing study 07A-CLI 007 Winter 2021 (n=404)

³ <https://www.lucirahealth.com/technical-brief/>