



SAFE life™

T-Cup® T-Cup® Compact Multi-Drug Urine Test Cup

Wondfo has tested our T-Cup and C-cup and has determined they can withstand a 25.1 inHg (inches mercury) vacuum, the maximum obtainable value in our vacuum chamber, for over 10 minutes without any sign of specimen leakage.

According to FedEx, packages can experience as much as 12inHg vacuum in the cargo compartment of their planes. Generally, this can result in a partial loss or total loss of specimen during overnight air transit to the laboratory. With our SAFElife urine cups you can rest assured there will be no leakage on the way to the lab.

