

PCR 8-Strip Tubes Caps

Features

1. Made from USP grade VI polypropylene.
2. Compatible with both NEST PCR 8-strip tube and PCR plate.
3. Domed and flat caps are available
4. Tight sealing, easy to open
5. Flat cap is ideal for qPCR
6. The design of domed cap can better prevent liquid evaporation


Product Range

Name	Cat.No.	Cap Type	Color	/Pack	/Case
PCR 8-Strip	406012	Flat	Transparent	125	1250
Tubes Caps	406112	Domed	Transparent	125	1250

Description

Temperature Range : Approx. -80°C to +121°C
 Shelf Life : 5 Years after date of production (ensure package is in good)
 Sterile : No

Technical Drawing



**NEST
China**

Tel : 0510-6878 8698 / 6878 8718
 Email : info@nest-wuxi.com
 Website : www.cell-nest.com
 Add : No.530, xida road, new district, Wuxi, Jiangsu, China

**NEST
USA**

Tel : +1-732-381-0268
 Email : info@nestscientificusa.com
 Website : www.nestscientificusa.com
 Add : 1592 Hart St., Rahway, NJ07065, USA