

# hemochroma<sup>PLUS</sup>

The hemochroma PLUS<sup>®</sup> Analyzer provides an accurate hemoglobin and hematocrit measurement.



- Reagent-free microcuvettes
- Microcuvette shelf life of up to 2 years, even after the canister has been opened
- Cleared for use on patients 6 months or older
- Measurement range 5.0 to 25.6 g/dL
- Small 15  $\mu$ l sample volume
- 3.25" Color LCD display
- CLIA Waived



#### Reimbursement

Test	Hemoglobin
CPT Code	85018-QW
Blood count; hemoglobin (Hgb)	

