



EMPLOYMENT & INSURANCE

Oral Fluid Drug Screen (EUO)

MRK220009

Oral fluid on-site drug testing allows for the delivery of fast and accurate results. A streamlined specimen collection and test procedure facilitates drug screening for employment and insurance purposes*.

Testing oral fluid for drugs is a non-invasive method to allow for direct observation of specimen collection which reduces risk of specimen tampering. In addition, the collection swab features a saturation indicator that ensures ample specimen volume is collected prior to testing.

FEATURES

- Highly accurate, sensitive and specific
- Simultaneous detection of multiple drugs
- Qualitative results within 10 minutes
- Lightweight and portable
- Non-invasive specimen collection
- Swab saturation indicator ensures ample specimen volume
- Lid closure ensures sample integrity
- Sample port for confirmation testing
- 16 strip capacity
- Customizable panels



BENEFITS

- On-site drug testing
- Direct observation of specimen collection
- Rapid turnaround of results
- Elimination of shy bladder issues
- Self-contained sample allows for shipment to laboratory for confirmation
- Easily administered test procedure

TEST MENU

Test	Cut-off (ng/mL)	Test	Cut-off (ng/mL)
Amphetamine (AMP)	50	Marijuana (THC)	40
Barbiturates (BAR)	50	Methadone (MTD)	30
Benzodiazepines (BZO)	10	Methamphetamine (MET)	50
Buprenorphine (BUP)	5	Methaqualone (MQL)	100
Cocaine (COC)	20	Opiates (OPI)	40
Ecstasy (MDMA)	50	Oxycodone (OXY)	20
LSD	10	Propoxyphene (PPX)	50
Heroin (6-MAM)	10	Tricyclic Antidepressants (TCA)	100

*Exemption is limited to test systems intended for employment and insurance testing that include a statement in their labeling that the device is intended solely for use in employment and insurance testing, and does not include test systems intended for Federal drug testing programs (e.g., programs run by the Substance Abuse and Mental Health Services Administration (SAMHSA), the Department of Transportation (DOT), and the U.S. military.)

AVAILABLE PANELS

Oral Fluid Drug Screen (EUO)

MRK220009

Initial offerings include panel sizes ranging from 4 to 11 tests.

PANEL SIZE	SKU GBDSA-XXXX	TEST											
4	3245 ESI-EUO	AMP50	COC20	MET50	OPI40								
	9344 ESI-EUO	AMP50	COC20	OPI40	THC40								
	9444 ESI-EUO	COC20	MET50	OPI40	THC40								
5	31155 ESI-EUO	AMP50	COC20	MET50	OPI40	OXY20							
	9154 ESI-EUO	AMP50	COC20	MET50	OPI40	THC40							
6	31965 ESI-EUO	AMP50	BZO10	COC20	MET50	OPI40	OXY20						
	91864 ESI-1-EUO	AMP50	COC20	MET50	OPI40	OXY20	THC40						
	9264 ESI-2-EUO	AMP50	BZO10	COC20	MET50	OPI40	THC40						
7	91174 ESI-EUO	AMP50	COC20	MET50	MTD30	OPI40	OXY20	THC40					
8	3985 ESI-EUO	AMP50	BAR50	COC20	MDMA50	MET50	MTD30	OPI40	OXY20				
9	91094 ESI-EUO	AMP50	BUP5	BZO10	COC20	MET50	MTD30	OPI40	OXY20	THC40			
	9894 ESI-EUO	AMP50	BAR50	BZO10	COC20	MET50	MTD30	OPI40	OXY20	THC40			
	9994 ESI-EUO	AMP50	BAR50	COC20	MDMA50	MET50	MTD30	OPI40	OXY20	THC40			
10	318105 ESI-EUO	AMP50	BAR50	BUP5	BZO10	COC20	MDMA50	MET50	MTD30	OPI40	OXY20		
	91104 ESI-1-EUO	AMP50	BAR50	BUP5	BZO10	COC20	MET50	MTD30	OPI40	OXY20	THC40		
11	95114 ESI-EUO	AMP50	BAR50	BUP5	BZO10	COC20	MDMA50	MET50	MTD30	OPI40	OXY20	THC40	

TEST PROCEDURE

