

90X53MM, ID LABEL W/TEMP.STRIP

# Urine Container

The urine cups are made of 100% virgin polypropylene resin, attached screw cap and tamper tab patient ID label meets JACHO standards. Stackable design helps lab space-saving. Pre-assembled graduated sterile cups have a wide mouth 53mm diameters and With temp.strip ranges from 90F°-100F° / 32°C -38°C, ideal for leakproof specimen transport.



Available in different color caps for option




**Model#**  
JYB-SPUC-90mL

**Sterilization**  
EO

**Stackable**  
Space - Saving

**Package**  
Individual Pack  
Case of 400



Popular Volume



Wide-Mouth



Leak-Resistant

120ML URINE CUP W/NEEDLE & VACCUM TUBE

# Transfer Device

This urine collection cup with the integrated transfer device, a full set includes a urine container (needle built-in) and vacuum tubes. Due to the special design of transfer lids, it's convenient for laboratory staff and donors to collect urine, pour urine into a tube, and processing. Each cup is sterilization and can be traced by lot number, provide temp.strip attached for option.




Available in different color caps for option



● Plain tube, 5/10ml



● Collection Cup / 4oz



**Model#**  
JYB-VAC-CUP-T

**Sterilization**  
Yes

**Temp.Strip**  
Available

**Package**  
Individual Pack  
Case of 250