

# BACtrack

## S35 Breathalyzer



The BACtrack S35 is one of the smallest, sleekest breathalyzers on the market. Simply blow into the replaceable mouthpiece and a Blood Alcohol Content (BAC) estimate is displayed on the cool blue screen in just seconds.

### Key Features

- Powerful proprietary sensing technology delivers highest possible accuracy in a unit this size.
- Compact size and sleek design fits into a pocket or purse for on-the-go lifestyles.
- Crisp three digit multifunction display.

Replaceable Hygienic Mouthpieces

Three Digit LCD Display



Ultra-portable Design

## Accurate

The BACtrack S35 comes equipped with an Advanced Semiconductor sensor and an algorithm specially tuned for its intended use. This enables it to deliver accuracy never before seen in a breathalyzer this small and affordable.

## Insightful

- Discover how your body is affected by alcohol.
- Quantify the effect of each drink on your Blood Alcohol Content.
- Make smart, more informed decisions.



### Product Specifications

Dimensions	2.5 x 2.25 x 0.66 in.
Weight	1.7 oz (47g)
Sensor Technology	Advanced Semiconductor
Detection Range	0.00-0.400% BAC
Mouthpieces	3 included
Power Supply	Two AAA batteries
Warm up time	15 seconds
Blowing time	5 seconds
Operating Temperature	50-104 °F
Warranty	1 year

## Portable

- Designed in California with ease of use and style in mind.
- Fits easily into a pocket or purse for testing family or friends anywhere you go.
- Smooth palm-of-hand feel enables discreet self-testing.
- Caribbean-blue LCD screen is easy on the eyes—especially in low light conditions.