

microINR[®]



LIMITED TIME OFFER



microINR the most efficient
INR Point-of-Care device
in the market

BUNDLE OFFER
SKA0001AD
METER + 6 BOXES



**friendly
technology**

SMART MONITORING OF WARFARING THERAPY

- ✓ Smallest blood sample required (3µL vs. 8-12µL)
- ✓ Direct application to test chip (No pipettes needed)
- ✓ Automatic Strip Identification (No need to match lot numbers)
- ✓ Auto Calibration (Parameters encoded within each individual chip)
- ✓ On-board and Internal Controls (in both chip AND monitor)
- ✓ No need of external liquid controls

www.ilinemicrosystems.com