

microINR[®]



**LIMITED TIME
OFFER**



BUNDLE OFFER
FMA0001AD
Buy 12 boxes of chips
and get a FREE meter!

**microINR,
the new CLIA waived INR
point-of-care device!**



**friendly
technology**

SMART MONITORING OF WARFARING THERAPY

- ✓ Small blood sample required (3µL vs. 8-12µL)
- ✓ Direct application to test chip (No pipettes needed)
- ✓ Automatic Strip Identification (No need to match lot numbers)
- ✓ Auto Calibration (Parameters encoded within each individual chip)
- ✓ On-board and Internal Controls (in both chip AND monitor)
- ✓ No need of external liquid controls