

RAPID DIAGNOSTICS


osom[®]

BVBLUE[®] Test

FOR THE DETECTION OF SIALIDASE ACTIVITY FROM GARDNERELLA, BACTEROIDES, PREVOTELLA, & MOBILUNCUS



COMPARISON OF BV
DETECTION METHODS*



www.osomtests.com



FEATURES:

- ✓ **ACCURATE;** 92.8% sensitive, 98% specific vs. gram stain
- ✓ **FAST;** Less than 1 minute hands-on time; result in 10 minutes
- ✓ **SAMPLE;** Vaginal swab
- ✓ **OBJECTIVE;** Instant color change provides clear, easy-to-read results
- ✓ **TRAINING;** Online training and certification; www.osomtraining.com
- ✓ **CUSTOM CARE;** Technical support staffed with certified medical technologists



RAPID DIAGNOSTICS

www.osomtests.com

osom[®]

BVBLUE[®] Test

Sekisui Diagnostics, LLC (AMERICAS)
4 Hartwell Place Lexington, MA 02421
Phone: 800 332 1042
Fax: 800 762 6311
Email: osomrapidtests@sekisui-dx.com
ISO 13485 Certified

Sekisui Diagnostics (UK) Limited (INTERNATIONAL)
50 Gibson Drive, Kings Hill, West Malling, Kent, ME19 4AF, UK
Phone: +44 (0) 1732 220022 Fax: +44 (0) 1732 220024
Email: info@sekisui-dx.com
ISO 9001 Certified, ISO 13485 Certified

TEST PROCEDURE / RESULTS



DESCRIPTION	CONFIGURATION	CATALOG #	STORAGE	CPT CODE
OSOM [®] BVBLUE [®] TEST	25 TESTS	183	REFRIGERATED	87905QW
OSOM [®] BVBLUE [®] CONTROL KIT		184	REFRIGERATED	
OSOM [®] VAGINITIS TRANSPORT SYSTEM (TRICH & BV)	50 SWAB SETS	185	ROOM TEMP	
OSOM [®] BVBLUE [®] TRANSPORT SWAB	30 SWABS	186	ROOM TEMP	

OSOM[®] is a registered U.S. trademark of Sekisui Diagnostics, LLC.
BVBLUE[®] is a registered trademark of Gryphus Diagnostics, LLC.
©2013 Sekisui Diagnostics, LLC – All rights reserved.

* See package insert.

