

QUICKVUE[®]+

One-Step
hCG Combo **TEST**

A Simple Test You Can Depend On For The Early Detection Of Pregnancy

CLIAwaived[™].com

For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121
Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739
E-mail: info@cliawaived.com
Visit us on the web at: www.cliawaived.com
Distributed by: CLIAwaived, Inc.



Rapid Diagnosis

- Results in 5 minutes (serum), 3 minutes (urine).
- Strong positive may become visible sooner.

Excellent Performance

- >99% sensitivity; specificity; accuracy.
- 10 mIU/mL sensitivity in serum or 20 mIU/mL sensitivity in urine.

Convenient

- Room temperature storage.
- 21-month shelf life.

Easy to Use

- Simple, one-step procedure uses only 4 drops of sample.
- Two-color plus / minus (+ / -) endpoint for easy-to-read and interpret results.

Reliable

- Highly sensitive assay with built-in controls.

Lot-to-Lot Consistency


- Highly controlled and monitored manufacturing environment.

QUIDEL[®]
CORPORATION

For more information on the QuickVue+ One-Step hCG Combo test kit, visit our website at www.quidel.com, or call your Quidel Account Manager at 1.800.874.1517 (in the U.S.) or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • www.quidel.com

Procedure:



CLIAwaived.comTM



To order: Call your Quidel Account Manager at 1.800.874.1517

CLIAwaived.comTM



For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121
Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739
E-mail: info@cliawaived.com
Visit us on the web at: www.cliawaived.com
Distributed by: CLIAwaived, Inc.

QUIDEL[®]
CORPORATION

For more information on the QuickVue+ One-Step hCG Combo test kit, visit our website at www.quidel.com, or call your Quidel Account Manager at 1.800.874.1517 (in the U.S.) or 858.552.1100 (outside U.S.).

Quidel Corporation ▶ Worldwide Headquarters ▶ 10165 McKellar Court ▶ San Diego, CA 92121 USA ▶ www.quidel.com

Distributor Information: