

CardioChek® Plus

Faster. Easier. Smarter...It All Adds Up.

eGLU™

[electrochemical glucose test strip]

- + Small blood sample; apply to eGLU glucose test strip directly from a fingerstick blood drop.
- + Glucose result in approximately 10 seconds as a stand-alone test.
- + Can be run as a stand-alone test or in conjunction with existing PTS Panels® Lipid Panel test strips.

3 QUICK STEPS

Lipid+eGLU test

Faster.

1


Insert MEMO Chip®. Insert test strips.



Easier.

2



Apply blood to eGLU & Lipid test strips.



Smarter.

3

Read results.



CardioChek® Plus eGLU™ Test Strip Ordering Information

Product Description*	Total Units	REF #
CardioChek® Plus Analyzer	1 analyzer	2700
eGLU™ Test Strips** (electrochemical glucose)	50 tests per vial	2713

*All existing PTS Panels® test strips and test strip panels are compatible with the CardioChek® Plus analyzer. Ordering information for the other available test strips not listed can be found by contacting your local sales representative, or by visiting our website at: www.cardiochek.com.

**The stand-alone eGLU has a single yellow stripe on the test strip box. The stand-alone test is designated with a Lot prefix of "X."

CardioChek® Plus
Faster. Easier. Smarter...It All Adds Up.

Lipid+eGLU™ Smart Bundle™ pack

+ Lipid + eGLU bundles contain identical Lipid Panel and eGLU test strips found in the stand-alone package configurations.

+ Unique bundle design and Lot # designation allow Lipid Panel and eGLU tests to be performed simultaneously.

CardioChek® Plus Analyzer Ordering Information

Product Description*	Total Units	REF #
CardioChek® Plus Analyzer	1 analyzer	2700

*All existing PTS Panels® test strips are compatible with the CardioChek® Plus analyzer. Ordering information for the other available test strips not listed can be found by contacting your local sales representative, or by visiting our website at: www.cardiochek.com.

CardioChek® Plus Lipid+eGLU™ Smart Bundle™ Pack

Product Description	Total Units	REF #
CardioChek® Plus Lipid+eGLU™ Smart Bundle™ pack	1 Smart Bundle™*	2727

*Each Smart Bundle includes (3) boxes of 15 ct. ea. Lipid Panel test strips (bundle version) and (1) box of 50 ct. ea. eGLU test strips (bundle version). There are (4) MEMo Chips® total that are included with each bundle set. These are calibrated together and have the same Lot #, which is required to run the Lipid Panel+eGLU tests simultaneously.

STAND-ALONE



(1)

SMART BUNDLE



+

(3)

STAND-ALONE



SMART BUNDLE



=

Lipid+eGLU™ test strips
Smart Bundle™ pack

(4)



Distributed by:

CLIAwaived™ Inc.

858.481.5031 | cliawaived.com