

CDS NextGeneration™

CDS Hematology Calibrator

KIT LOT: 1008-731 EXP: 19-Oct-10

Open Vial Stability: 5 days

Analyzer	WBC x10 ⁹ /L	RBC x10 ¹² /L	Hgb g/dL	Hct %	MCV fL	PLT x10 ⁹ /L	MPV fL
MINDRAY 3000 Series	7.4	4.52	12.7	38.5	85.1	242	10.5
Range Limits:	± 0.3	± 0.10	± 0.2	± 1.0	± 2.0	± 15	1.5

Range Limits are provided for evaluating analyzer data, and for determination of analyzer reproducibility and operator skill.
All trademarked analyzer names belong to the respective trademark holders.



Distributed by: **CLIAwaived.com™**
San Diego, CA 92121
tel 858-481-5031
toll free 888-882-7739
www.cliawaived.com

PN 201058A R01.15.08

