

Medonic *M-SERIES* Hematology Analyzers

Open Vial Sampling System

Cap Piercing Sampling System

Choose The Analyzer
That Best Fits
Your Needs

Simple

Accurate

Reliable

Economical



What Differentiates The *M-SERIES*

Accurate and Precise — Delivers excellent patient results.

Bar Coding — Standard with all systems giving you positive patient identification and reduces transcription errors.

Zero Daily Maintenance — No daily shutdown, engineered to eliminate lengthy cleaning procedures.

Service — Highly reliable, excellent uptime record.



Patent Finger Stick Technology

20µ **draw is all it takes**

Eliminates the need for costly
drawing devices

Automate Your Records & Data

Integrate a Data Management Package
with your system.

- Easily interfaced to Data Management Systems
- Presently compatible with several systems on the market





Clinical Diagnostic Solutions, Inc.

1800 NW 65th Avenue
Plantation FL, 33313

Local: 954.791.1773
Toll-Free: 800.453.3328
Fax: 954.791.7118
www.cdsolinc.com

Specifications

16-Parameter, 3 Part Diff
Up to 60 Samples Per Hour
Results in 45 Seconds
Micro Sampling
LCD Color Touch Screen Display
Histogram on Screen Display
On Board Step By Step Operation Instructions
Bar Code Reader
Sample Memory
Control Memory
Program Files
Clot Removal
On Board QC Including LJ, XB
Blood Detectors
Auto Hgb Blank
Hgb Correction

Reagent Level Sensors
Warning Flags
Parameter Flags
Interface Ready
Zero Daily Maintenance

Options:

Key Board

Small Footprint

16.15 x 11.43 x 18.12

Weight

39.6 Pounds - Open / Closed Vial
Power Consumption 100VA
Power Consumption (stand-by) 20VA
Frequency 50/60HZ



Boule Medical AB
Box 42056
SE-126 13 Stockholm,
Sweden

Copyright © 2013 Clinical Diagnostic Solutions (CDS). All rights reserved.
Clinical Diagnostic Solutions and the CDS logo are properties of Clinical Diagnostic Solutions, Inc.
Medonic and the Medonic logo are properties of Boule Medical AB, Stockholm, Sweden.

P/N 202325A R02.22.13