

BACtrack

Mobile Breathalyzer

Introducing the world's most advanced breathalyzer.



Get more out of every test.

BACtrack Mobile makes it easy to be smart about drinking. Test, track, and share BAC results, to gain a better understanding of how alcohol affects you.

Technical Specs:

Weight	1.75 oz (47g)
Mouthpieces	3 included
Battery	Rechargeable with included micro USB cable
Sensor Technology	Xtend™ Electro-chemical fuel cell
Warm-up Time	10-20 Seconds
Blowing Time	5 Seconds
BAC Range	0.000 - 0.400% BAC (Blood Alcohol Content)
Warranty	1 Year

*Actual Size

Compact, lightweight design

Sleek shape with stunning LED indicator



2.75"


Quick pairing with Bluetooth 4.0

Long-lasting rechargeable lithium ion battery

Test

Fuel cell sensor provides professional-grade accuracy.

Internal air pump ensures precise, consistent results every time.




Track

See your BAC rise and fall over time.

- ▶ View results from the last 6 hours to the last 14 days and beyond.
- ▶ Track what you drank and where.

Learn

- ▶ Gain insight into your drinking habits and learn how your body processes alcohol.
- ▶ With ZeroLine™, estimate when your BAC will return to 0.00%.
- ▶ Your BAC results remain private, or chose to share them via text message, Facebook, or Twitter.



DISTRIBUTED BY:

CLIAwaived™ Inc.

11578 Sorrento Valley Road #26 | San Diego, CA 92121 | 858.431.5031 | cliawaived.com