

CDS M-SERIES Hematology Analyzers

Open Vial Sampling System

Cap Piercing Sampling System

Choose The Analyzer
That Best Fits
Your Needs

Simple

Accurate

Reliable

Economical



What Differentiates The M-SERIES

Accurate and Precise — Delivers excellent patient results.

Bar Coding — Standard with all systems giving you positive patient identification and reduces transcription errors.

Zero Daily Maintenance — No daily shutdown, engineered to eliminate lengthy cleaning procedures.

Service — Highly reliable, excellent uptime record.



Patent Finger Stick Technology

20µ draw is all it takes

Eliminates the need for costly
drawing devices

Automate Your Records & Data

Integrate a Data Management Package
with your system.

- Ensures compliance with all patient regulations regarding patient data, patient records, quality control and general record keeping.
- Easy to use.
- Affordable solution for the physician's office.
- Easy interface to EMR.





Clinical Diagnostic Solutions, Inc.

CLIAwaivedTM.com

For more information please contact CLIAwaivedTM Inc. at:

TEL: (858) 481-5031 TOLL FREE: (888) 882-7739

E-MAIL: info@cliawaived.com WEB: www.cliawaived.com

Distributed by: CLIAwaivedTM Inc., San Diego, CA

Specifications

16-Parameter, 3 Part Diff
Up to 60 Samples Per Hour
Results in 45 Seconds
Micro Sampling
LCD Color Touch Screen Display
Histogram on Screen Display
On Board Step By Step Operation Instructions
Bar Code Reader
Sample Memory
Control Memory
Program Files
Clot Removal
On Board QC Including LJ, XB
Blood Detectors
Auto Hgb Blank
Hgb Correction
Reagent Level Sensors
Warning Flags
Parameter Flags
Interface Ready
Zero Daily Maintenance

Options:

Key Board
Data Management Package

Small Footprint

16.15 x 11.43 x 18.12

Weight

39.6 Pounds - Open / Closed Vial
Power Consumption 100VA
Power Consumption (stand-by) 20VA
Frequency 50/60HZ