

# CDS M<sup>c</sup>SERIES Hematology Analyzers

Open Vial Sampling System

Cap Piercing Sampling System

Choose The Analyzer  
That Best Fits  
Your Needs

Simple

Accurate

Reliable

Economical



## What Differentiates The M<sup>c</sup>SERIES

**Accurate and Precise** — Delivers excellent patient results.

**Bar Coding** — Standard with all systems giving you positive patient identification and reduces transcription errors.

**Zero Daily Maintenance** — No daily shutdown, engineered to eliminate lengthy cleaning procedures.

**Service** — Highly reliable, excellent uptime record.



### Patent Finger Stick Technology

20µ draw is all it takes

Eliminates the need for costly  
drawing devices

### Automate Your Records & Data

Integrate a Data Management Package with your system.

- Ensures compliance with all patient regulations regarding patient data, patient records, quality control and general record keeping.
- Easy to use.
- Affordable solution for the physician's office.
- Easy interface to EMR.





Clinical Diagnostic Solutions, Inc.

Distributed by: **CLIAwaived.com™**  
San Diego, CA 92121  
tel 858-481-5031  
toll free 888-882-7739  
www.cliawaived.com

### Specifications

16-Parameter, 3 Part Diff  
Up to 60 Samples Per Hour  
Results in 45 Seconds  
Micro Sampling  
LCD Color Touch Screen Display  
Histogram on Screen Display  
On Board Step By Step Operation Instructions  
Bar Code Reader  
Sample Memory  
Control Memory  
Program Files  
Clot Removal  
On Board QC Including LJ, XB  
Blood Detectors  
Auto Hgb Blank  
Hgb Correction  
Reagent Level Sensors  
Warning Flags  
Parameter Flags  
Interface Ready

Zero Daily Maintenance

### Options:

Key Board  
Data Management Package

### Small Footprint

16.15 x 11.43 x 18.12

### Weight

39.6 Pounds - Open / Closed Vial  
Power Consumption 100VA  
Power Consumption (stand-by) 20VA  
Frequency 50/60HZ