



the GE100 advantage

The Top 5 Reasons

Why The GE 100 Is Well Suited For Your Diabetes Management




Reason 1: Large Screen Size

- One of the largest screen sizes on the market measuring at **1674mm²** aiding those with visual impairment.

Reason 2: Innovative Strip Design

- Larger, thicker test strips help patients with dexterity issues and are less likely to bend.
- Unique handle allows handling the strip without touching the sample area.



Reason 3: Accuracy

- Code embossed on every strip and utilizes gold electrodes for increased stability

Reason 4: Price

- The GE100 offers your patients both quality and accuracy at a price at or below most co-pay insurance programs.



Reason 5: A Brand You Can Trust

- GE consistently ranks as one of the world's leading corporations as well as one of the world's most recognizable brands.



Why Co-Pay?.....Try Low Pay!

Are your patients using Co-Pay for their Diabetes Management?
With the GE100, they can now use "Low-Pay".

GE100 Blood Glucose Monitoring System



True Auto Coding

Testing made simple.

Accurate

Reliable results in about 5 seconds.

Large Screen

Know your results in a glance



Clean and No Mess

Unique strip handle reduces unnecessary sample handling.

Easy to Handle

Easy to hold, use & test.

Tiny Blood Sample

Small sampling with quick & accurate results.



The GE100 offers your patients quality, accuracy and a brand you can trust at a price at or below most co-pay insurance programs.

Description	NDC #	UPC #
GE100 Meter Kit	08539-7100-02	883489000781
GE100 Blood Glucose Test Strip(25CT)	08539-8100-02	883489000767
GE100 Blood Glucose Test Strip(50CT)	08539-8100-01	883489000750
GE100 Blood Glucose Test Strip(100CT)	08539-8100-03	883489000774
GE100 CONTROL SOLUTION (NORMAL- 4ML)	08539-9001-00	883489000729
GE LANCING DEVICE	08539-9004-00	883489000743