

*INSTANT-VIEW®*

# Drugs of Abuse Tests

ONE-STEP DRUGS OF ABUSE IMMUNOASSAY DEVICE



**FDA 510(K)-CLEARED** FOR PROFESSIONAL USE  
**CE MARKED** FOR SALE IN EUROPEAN UNION  
**MADE IN USA**

**SIMPLE • FAST • SENSITIVE**  
**SPECIFIC • CONVENIENT**

**SIMPLE**  
**One-Step Test**

**FAST**  
**Results in 4-7 Minutes**

**CONVENIENT**

**Room Temperature Storage**  
**Built-in Procedural Control**  
**Test up to 12 Drugs of Abuse at One Time**  
**Pipette or Dip-and-Read Procedure**  
**Adulteration Controls Available**

**SENSITIVE**

**SAMHSA Defined Levels**

**SPECIFIC**

**92.2%-99.8%**



Product Configuration	Sensitivity	Kit Size	INSTANT-VIEW® Kit Catalog	Bulk Catalog
Amphetamine Urine Cassette Test (AMP)	1,000ng/ml	25	03-3012	08-3022
Amphetamine Urine Dip Strip Test (AMP)	1,000ng/ml	50	03-3017	08-3027
Barbiturate Urine Cassette Test (BAR)	200ng/ml	25	03-3152	08-3162
Barbiturate Urine Dip Strip Test (BAR)	200ng/ml	50	03-3157	08-3167
Benzodiazepine Urine Cassette Test (BZD)	300ng/ml	25	03-3132	08-3142
Benzodiazepine Urine Dip Strip Test (BZD)	300ng/ml	50	03-3137	08-3147
Cocaine Urine Cassette Test (COC)	300ng/ml	25	03-3052	08-3062
Cocaine Urine Dip Strip Test (COC)	300ng/ml	50	03-3057	08-3067
MDMA (Ecstasy) Urine Cassette Test (XTC)	500ng/ml	25	03-2952	08-2962
MDMA (Ecstasy) Urine Dip Strip Test (XTC)	500ng/ml	50	03-2957	08-2967
Methamphetamine Urine Cassette Test (MET)	1,000ng/ml	25	03-3032	08-3042
Methamphetamine Urine Dip Strip Test (MET)	1,000ng/ml	50	03-3037	08-3047
Methamphetamine Urine Cassette Test (MET)	500ng/ml	25	03-3033	08-3043
Methamphetamine Urine Dip Strip Test (MET)	500ng/ml	50	03-3038	08-3048
Morphine Urine Cassette Test (MOR)	2,000ng/ml	25	03-3092	08-3102
Morphine Urine Dip Strip Test (MOR)	2,000ng/ml	50	03-3097	08-3107
Morphine Urine Cassette Test (MOR)	300ng/ml	25	03-3093	08-3103
Morphine Urine Dip Strip Test (MOR)	300ng/ml	50	03-3098	08-3108
Methadone Urine Cassette Test (MTD)	300ng/ml	25	03-3172	08-3182
Methadone Urine Dip Strip Test (MTD)	300ng/ml	50	03-3177	08-3187
Phencyclidine Urine Cassette Test (PCP)	25ng/ml	25	03-3112	08-3122
Phencyclidine Urine Dip Strip Test (PCP)	25ng/ml	50	03-3117	08-3127
Propoxyphene Urine Cassette Test (PPX)	300ng/ml	25	03-2932	08-2942
Propoxyphene Urine Dip Strip Test (PPX)	300ng/ml	50	03-2937	08-2947
Tricyclids Urine Cassette Test (TCA)	1,000ng/ml	25	03-2972	08-2982
Tricyclids Urine Dip Strip Test (TCA)	1,000ng/ml	50	03-2977	08-2987
Marijuana Urine Cassette Test (THC)	50ng/ml	25	03-3072	08-3082
Marijuana Urine Dip Strip Test (THC)	50ng/ml	50	03-3077	08-3087

Product Configuration	Panel Size	Kit Size	INSTANT-VIEW® Kit Catalog
MET/THC	2-Panel	25	03-3226
COC/THC	2-Panel	25	03-3270
COC/MET/THC	3-Panel	25	03-3371
COC/MOR/THC	3-Panel	25	03-3373
COC/MET500/MOR300/THC	4-Panel	25	03-3445
BZD/COC/MOR/THC	4-Panel	25	03-3451
AMP/COC/MOR/THC	4-Panel	25	03-3458
COC/MET/MOR/THC	4-Panel	25	03-3474
COC/MET/MOR/PCP/THC w/adulteration	5-Panel	25	03-3502
COC/MET500/MOR300/PCP/THC	5-Panel	25	03-3548
COC/MET/MOR/THC w/UID	5-Panel	25	03-3570
BZD/COC/MET/MOR/THC	5-Panel	25	03-3571
AMP/COC/MOR/PCP/THC	5-Panel	25	03-3576
AMP/COC/MET/MOR/THC	5-Panel	25	03-3577
COC/MET/MOR/PCP/THC	5-Panel	25	03-3578
AMP/BZD/COC/MET/MOR/THC	6-Panel	25	03-3635
AMP/COC/MET/MOR/THC w/UID	6-Panel	25	03-3677
AMP/COC/MET/MOR/PCP/THC	6-Panel	25	03-3679
AMP/BAR/BZD/COC/MET/MOR/THC w/UID	8-Panel	20	03-3870
AMP/BAR/BZD/COC/MET/MOR/MTD/PCP/THC	9-Panel	20	03-3997
AMP/BAR/BZD/COC/MET/MOR/MTD/PCP/TCA/THC	10-Panel	20	03-3910
AMP/BAR/BZD/COC/MET/MOR/MTD/PCP/PPX/TCA/THC/XTC	12-Panel	20	03-3951