



Clearly different.
Clearly better.



Women's Health | HIV / STDs | Respiratory | Cardiovascular | Enterics & GI | General Health

Clearview® hCG

Clearview hCG Testing

One-step pregnancy tests for early diagnosis.

Clearly different

Clearview's patented innovative technology provides highly sensitive and specific tests for pregnancy.

- **Rapid Response:** Single-step; positive results in as little as 1 minute.
- **Simplicity:** CLIA waived (urine), easy-to-use test.
- **Reliability:** Highly sensitive assay with built-in controls ensures accuracy.
- **Flexibility:** Use with multiple sample types — urine or serum.

Clearly better

Clearview provides better patient care with early diagnosis for proper planning, reducing potential risks to mother and fetus.

- Earlier detection allows for improved care.
- Rapid results help prevent infection or contraindicated treatment.



Point of Care Diagnostics

CLIAwaived™.com

For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121
Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739

E-mail: info@cliawaived.com

Visit us on the web at: www.cliawaived.com

Distributed by: CLIAwaived, Inc.

**Clearview hCG Combo II**

40 Cassettes

Order No. 92215

Urine - CLIA Waived

Serum - Non-waived

Clearview 25 hCG Combo

40 Cassettes

Order No. 92225

Urine - CLIA Waived

Serum - Non-waived

Clearview hCG Urine

50 Dipsticks

Order No. 92211

CLIA Waived

Clearview hCG Urine

40 Cassettes

Order No. 92217

CLIA Waived

hCG Controls

1 Positive, 1 Negative


Order No. 92218

Clearview® hCG Combo Tests Comparison Grid


Test Name	Clearview hCG Combo II & 25 hCG Combo	Quidel QuickVue® One-Step hCG Combo	Beckman Coulter ICON® 25 hCG	Genzyme OSOM® hCG Combo Test
Format	40 Cassettes	50 Cassettes	25 Cassettes	25 Cassettes
CLIA Status: Sample Types	Waived: Urine Non-waived: Serum	Waived: Urine Non-waived: Serum	Waived: Urine Non-waived: Serum	Waived: Urine Non-waived: Serum
Sensitivity	hCG Combo II 20 mIU/mL Urine 10 mIU/mL Serum 25 hCG Combo 25 mIU/mL	25 mIU/mL	25 mIU/mL	20 mIU/mL Urine 10 mIU/mL Serum
Specificity	>99%	>99%	>99%	>99%
Test Steps	1	1	1	1
Minutes to Result	(+) = 1* (-) Urine = 3 (-) Serum = 5	Urine = 3 Serum = 5	Urine = 3 Serum = 5	Urine = 3 Serum = 5
Shelf Life	24 months	18 months	24 months	24 months
Storage	Room Temperature	Room Temperature	Room Temperature	Room Temperature

Clearview hCG Urine Dipstick ProcedureDip to line,
5 secondsRead in
3 minutes


Interpret results



Positive




Negative



Invalid

Clearview hCG Combo II / 25 Combo / Urine Cassette Procedure

Draw sample

Add 3 drops
to sample wellRead in
3 minutes - Urine
5 minutes - Serum

Positive



Negative



Invalid

Interpret results

visit

www.invernessmedicalpd.com

for more information

CLIAwaived.com

For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121
Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739
E-mail: info@cliawaived.com
Visit us on the web at: www.cliawaived.com
Distributed by: CLIAwaived, Inc.

CLIAwaived.com

For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121

Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739

E-mail: info@cliawaived.comVisit us on the web at: www.cliawaived.com

Distributed by: CLIAwaived, Inc.