



Respiratory | Women's Health & STD | Enterics & GI | General Health | Instrumentation

RUBELLA-plus[®]

Clearly Different

Wampole Rubella-plus provides simple, accurate detection of rubella antibodies.

- **Fast: results available in just 15 minutes**
- **Easy: no special training or rotator needed**
- **Reliable: detects antibodies at a level greater than or equal to 10 IU/mL**
- **Flexible: can be used with either serum or plasma**

Clearly better

The rapid, reliable results provided by Rubella-plus can lead to improved patient care.

- **Early detection in women of child-bearing age may prevent the transmission of the disease to unborn children**

Clearview

RUBELLA-plus

For the qualitative detection of Rubella antibody (IgG) in human serum or plasma.



Point of Care Diagnostics

CLIAwaived.com[™]

For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121
Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739

E-mail: info@cliawaived.com

Visit us on the web at: www.cliawaived.com

Distributed by: CLIAwaived, Inc.

Clearview
RUBELLA-plus
 30 Test Kit
 Order No. 84R2
 CLIA Non-waived

About Rubella

Rubella (also called German measles or 3-day measles) is a disease caused by the rubella virus. Rubella spreads from person to person through the air (often through close contact such as talking, coughing, or sneezing) and is moderately contagious.

Although rubella can strike people of all ages, **it poses the greatest danger to unborn babies.** Congenital rubella syndrome (CRS) occurs when the rubella virus attacks a developing fetus. Up to 85% of infants infected during the first trimester will be born with birth defects, including deafness, blindness, heart defects, and mental retardation.¹ Miscarriages are also common.²

Fortunately, due to widespread immunization with a rubella vaccine, rubella and CRS are now uncommon in the U.S.¹ Rubella vaccine is often given together with measles and mumps vaccines as part of an MMR shot.

RUBELLA-plus Test Procedure

Add 10 μ L of serum or plasma to the Sample well (S)
 Follow with 3 to 4 drops Developer Solution (D)
 Allow 15 minutes
 Interpret results

CLIAwaived.com
 For more information please contact CLIAwaived Inc. at
 11578 Sorrento Valley Road ■ San Diego, CA 92121
 Tel: (858) 481-5031 ■ Toll Free: (888) 882-7739
 E-mail: info@cliawaived.com
 Visit us on the web at: www.cliawaived.com
 Distributed by: CLIAwaived, Inc.

1. Immunization Action Coalition. Vaccine Information for the Public and Health Professionals: Rubella Questions and Answers. Available at: <http://www.vaccineinformation.org/rubella/qandadis.asp>. Accessed April 21, 2005.
2. Centers for Disease Control. Parents Guide to Childhood Immunization: Rubella (German Measles). Available at: <http://www.cdc.gov/nip/publications/Parents-Guide/pg-rubella.pdf>. Accessed April 21, 2005.

RUBELLA-plus[®] is a registered trademark of the Inverness Medical family of companies.

331 S. 104th Street 1.800.637.3717
 Louisville, CO 80027 USA P: 303.530.3888
www.invernessmedicalpd.com F: 800.637.3827

© 2006 Inverness Medical - BioStar Inc.
 02831.00

visit
www.invernessmedicalpd.com
 for more information