

MiniPet® Specifications

Volumes

5 L	10 L	11 L	14 L	20 L	25 L	30 L
35 L	40 L	50 L	75 L	100 L	200 L	

Performance (with generic Eppendorf® style 0-200 L Tip)

Volume	±% Accuracy	±% Precision (CV)
100 L and 200 L (full volume)	3	2
75 L	4	2
50 L (full volume)	3	2
35-45 L	4	2
20-30 L	5	3
10 L, 11 L and 14 L	10	4
5 L	20	5

Physical

Dimensions (instrument only)

3.75 in. or 9.53 cm length

.540 in. or 1.37 cm handle diameter

Weight

.20 oz or 6.25 gm (instrument only)

Appearance

Red	White	Pink
5 L	10 L	20 L

Orange	Plum	Light Blue	Light Green
25 L	30 L	35 L	40 L

Yellow	Green	Blue	Teal
50 L	75 L	100 L	200 L

Fluid Contact Material

Tip only (disposable)

Accessories

0-200 L Yellow Tips

DeTipper Bottle

CLIAwaived™

For more information please contact CLIAwaived™ Inc. at:

TEL: (858) 481-5031 TOLL FREE: (888) 882-7739

E-MAIL: info@cliawaived.com WEB: www.cliawaived.com

Distributed by: CLIAwaived™ Inc., San Diego, CA