

Hypoxic QC

For Validation of low pO_2 Blood Gas Measurements



Revolutionary New QC for ABG Labs:

- Assess accuracy and precision of low pO_2 measurements
- True hemoglobin buffering – 10 minute open ampule stability!
- Ideal for Method Comparisons and Analytical Measuring Range (AMR) validation
- Comparable to whole blood tonometry – w/o the tonometer!

Catalog #	Level	Target	Packaging
166.001.004	X-Low	15 mmHg	4 amps/box

Note: Due to its blood-like matrix, Hypoxic QC has a one month shelf life once shipped

Distributed by: **CLIAwaived.com™**
San Diego, CA 92121
tel 858-481-5031
toll free 888-882-7739
www.cliawaived.com