

The ART of Flu Testing

- Accurate, objective and automated results in **as few as 3 minutes**
- Secure, remote instrument management through the Virena ecosystem
- **Excellent performance** compared to molecular and viral culture methods
- Integrated data management
 automatically stores test and user history
- Flexible, dual mode testing and simple, CLIA-waived procedures for a variety of laboratory environments
- User and patient ID captured with **onboard barcode scanner**


Sofia 2 with **Advance Result Technology**

Delivering accurate, objective and automated results in as few as three minutes!






Sofia Influenza A+B FIA - Procedures



Sofia 2 Development Modes



READ NOW Mode – Batch multiple samples per hour



Sofia Influenza A+B FIA

- Meets reclassification requirements for rapid influenza antigen tests
- Includes Nasal swabs, NP swabs and NP aspirate/wash samples types
- Single specimen may be used for both Sofia RSV and Sofia Influenza
- Compatible with viral transport media
- Room temperature storage
- 21-month shelf life from date of manufacture

Sofia Influenza A+B FIA Results vs. Culture*

Samples		CLIA Waiver Studies	
		Sens	Spec
Influenza A	NS/NPS	93%	95%
	NPA/NPW	99%	96%
Influenza B	NS/NPW	90%	96%
	NPA/NPW	88%	96%

^{*}Refer to the Package Insert for additional performance claims.







