



---

## Clearview® Strep A Exact II Dipstick Quiz

Name: \_\_\_\_\_

Date: \_\_\_\_\_

**Circle T (True) or F (False) for each Question:**

- |  |   |   |
|--|---|---|
| 1. The Clearview® Strep A Exact II Dipstick must be refrigerated at 2-8°C.   | T | F |
| 2. After opening the foil pouch, you must use the test dipstick as soon as possible.   | T | F |
| 3. Throat swab is the only acceptable sample type for testing on the Clearview® Strep A Exact II Dipstick.                                   | T | F |
| 4. Throat swab specimens may be held up to 72 hours at refrigerated temperatures before performing the Clearview® Strep A Exact II Dipstick. | T | F |
| 5. The Clearview® Strep A Exact II Dipstick does not require live organism to be positive.   | T | F |
| 6. It is recommended that the swab provided in the kit be used.  | T | F |
| 7. Testing should ideally be performed 10 minutes after the specimens have been collected.   | T | F |
| 8. Test Results should be read at 5 minutes.   | T | F |
| 9. The appearance of a red Control line and a red Test line is a positive result.  | T | F |
| 10. External Positive and Negative Controls should be assayed once per kit or as deemed necessary by your laboratory.                        | T | F |